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
Environment: Temperature 25.5 °C Humidity 51 %RH

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
Environment: Temperature 25.5 °C Humidity: 51 %RH

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Device data (model and serial number)

ICP OES	PQ9000
Mobile cooling unit	PolyScience
Autosampler	-

Analytik Jena representative information

Name	Mr.Puktija somton
Company	Analytik-Jena Fareast Thailand

Maintenance date

Date	11/10/2022
------	------------

Table of contents

1	Customer and service data	3
2	Maintenance with subsequent Operational Qualification (OQ)	4
3	Maintenance work on the base unit	5
4	Maintenance work on accessories	8
5	Checking basic functions and device parameters	10
5.1	Check of correction functions	10
5.2	Check of safety circuits	10
5.3	Check of neon energy	11
5.4	Adjustment of transfer optics (Mn 257.610 nm)	11
5.5	Verification of the wavelength accuracy, generator robustness and analytical performance	12
6	Comments and objections	16
7	Completing maintenance	17

3 Maintenance work on the base unit

Update firmware and software

	Complies	Does not comply
Save device parameters and neon table	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the firmware version and FPGA	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the version of the ASpect PQ software	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Update software, FPGA and firmware	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Approval for firmware and software updates was given by the customer.		
The customer's was informed that the Operational Qualification (OQ) and software validation may have to be repeated for validated systems after a firmware/software update.		
Update ASpect PQ	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Current version	1.2.3.0	<input type="checkbox"/>
▪ Version after update	—	<input type="checkbox"/>
Update FPGA	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Current version	h09; h00	<input type="checkbox"/>
▪ Version after update	—	<input type="checkbox"/>
Updating the firmware	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Current version	1.36	<input type="checkbox"/>
▪ Version after update	—	<input type="checkbox"/>
Record running time (in hours)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Base unit	3031	<input type="checkbox"/>
▪ RF generator	1770	<input type="checkbox"/>

Service the rear of the device

	Complies	Does not comply
Replace the air filter mat	<input checked="" type="checkbox"/>	<input type="checkbox"/>

2 Maintenance with subsequent Operational Qualification (OQ)

- If an Operational Qualification (OQ) is carried out after maintenance, there is no need to check analytical parameters in the maintenance report.
- The analytical parameters are then checked in the separate OQ report.

Maintenance with subsequent Operational Qualification (OQ)	Yes	No
	<input type="checkbox"/>	<input checked="" type="checkbox"/>

Comments:

Service the spectrometer

		Complies	Does not comply
Clean and lubricate the spindle of the prism and grating with a little oil		<input checked="" type="checkbox"/>	<input type="checkbox"/>

Check gas flows at the gas box

Gas flow	Rated value	Actual value	Complies	Does not comply
Nebulizer gas flow	0.5 L/min ± 15 %	0.512	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Auxiliary gas flow	0.5 L/min ± 15 %	0.506	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Cone gas flow	1,7 ... 2,8 L/min	2.7	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Oxygen is used (optional)			<input type="checkbox"/>	<input type="checkbox"/>
Oxygen flow	00.5 L/min ± 15 %	-	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Service the sample introduction system

	Complies	Does not comply
Remove, disassemble and clean the torch	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Check the surface of the glass parts (injector, inner tube, outer tube and glass holder) for damage, heavily coated and broken parts. Replace glass parts, if necessary		
▪ Clean glass parts in aqua regia		
Check the sealing rings, replace if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
▪ Sealing rings in torch holder (3 pieces) and torch shuttle (2 pieces)		
▪ Sealing rings in injector and mixing chamber mount		
(Check sealing rings for porous and hardened areas, replace if necessary)		
Check pump tubing (sample, waste), replace if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check peristaltic pump, disassemble and lubricate if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the nebulizer for damage and clean it	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Remove and clean the mixing chamber	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the sealing ring of the mixing chamber	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check and clean torch height adjustment drive	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Lubricate the linear drive with a little oil	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the igniting pin for signs of wear	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Comments:

Service the plasma compartment

	Complies	Does not comply
Disassemble the cone, check geometry and clean it if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Take out the windows (axial/radial) and clean or replace them		
Replacing the sealing rings		
Check the surface of the first mirror, clean or replace it if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check the air outlet for deposits and blockages, clean it if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Remove the cover plate below the plasma compartment and clean it	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Renew the filter mat		
Check the coil for its material condition, position, geometry and for the tightness of connections	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Clean or replace bonnet	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Check RF deflector spring, replace it if necessary	<input checked="" type="checkbox"/>	<input type="checkbox"/>