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
Agilent CrossLab Start Up Services

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Agilent Preventive Maintenance provides factory recommended service for your analytical instruments to assure reliable operation and the accuracy of your results.

Delivered by highly trained and certified service engineers using genuine Agilent parts and supplies, Agilent Preventive Maintenance provides everything you need to reduce unplanned downtime and keep your systems operating at their peak. This checklist will be completed at the end of the service and provided to you as a record of the preventive maintenance activities.

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Important Customer Web Links

- For more information about Agilent Technologies services, please visit our website using the following URL: <http://www.agilent.com/en-us/products/crosslab-instrument-services/service-repair>
- The Agilent Community is an excellent place to get answers, collaborate with others about applications and Agilent products, and find in-depth documents and videos relevant to Agilent technologies. Visit <https://community.agilent.com/welcome>.
- To access Agilent University, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- A useful Agilent Resource Center web page is available, which includes short videos on maintenance, quick lists of consumables for new instruments, and other valuable information. Check out the Resource Page here: <https://www.agilent.com/en-us/agilentresources>.
- Need technical support, FAQs, supplies? – visit our [Support Home page](http://www.agilent.com/search/support) <http://www.agilent.com/search/support>.
- Videos about specific preparation requirements for your instrument can be found by searching the Agilent YouTube channel at <https://www.youtube.com/user/agilent>.

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.759722222
© Agilent Technologies, Inc. 2023

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- Contact the customer and ensure that all necessary supplies are available before the preventive maintenance visit.
- Only select those pages that relate to the system or module being serviced.
- Complete empty fields with the relevant information.
- Complete the relevant checkboxes in the checklist using either a "X" or tick mark "✓".
- Check "Section not applicable" checkboxes to indicate services/tasks not delivered, as appropriate.
- Complete the Preventive Maintenance service in the order of the tasks listed.
- Complete the Service Review section together with the customer.
- Complete the fields for page numbers at the foot of each selected page.
- Complete the total number of pages field in the Service Completion section.
- Ask the customer to sign the Service Completion section including the customer's and your signature.

Additional Instruction Notes

- Check for any active service notes for this unit. If there are any applicable "Safety" or "Modification Recommended" Service notes, plan to implement the changes on this unit before doing any qualification service.
- Do not implement firmware updates, unless you get approval from the customer and are sure that they are compatible with the instrument control software.

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.759722222
© Agilent Technologies, Inc. 2023

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System Information

- ☐ Check this box if an instrument configuration report is attached instead of completing the table below.

Instrument System Name and ID	CN17100005
Instrument System Site and Location	UAE

List System Component Product Numbers	List the Serial Numbers of each Component
1. G3950A	CN17100005
2. G4513A	CN17120171
3.	
4.	
5.	
6.	
7.	
8.	
9.	
10.	

Preparation

- ☒ Discuss any specific issues with the customer before starting.
- ☒ Review the instrument logbook for recorded problems and comments.
- ☒ Save instrument control settings before starting the procedure.
- ☒ Perform a general inspection of the system for cleanliness.
- ☒ Check for proper installation of parts, assemblies, sensors etc.
- ☒ Check system for required installation of components, settings as defined by current Service Notes.
- ☐ Check for required firmware updates and verify with customers if they would like them installed.
- ☒ Before starting the following procedures, record the Detector Signal Output(s) in the results table. If the GC is turned OFF or in a service mode, comparing the detector outputs before and after the service is not possible.

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.759722222
© Agilent Technologies, Inc. 2023

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Preventive Maintenance Procedure

Clean and inspect GC

- ☐ Unplug power cord from the power source.
- ☐ Open GC covers and vacuum/remove any dust/debris. Pay particular attention to cooling fans.
- ☐ Inspect internal connectors for proper contact and placement.
- ☐ Reconnect Power to the GC. Power the GC on and verify the power on self-test passed.
- ☐ Verify operation of all instrument fans.

Inlet and detector consumable replacement

- ☒ For the inlet installed, perform inlet maintenance using the built-in procedures accessed from Agilent 9000 touch screen display or web interface.
- ☒ Replace column Compression Bolts.
- ☒ Replace the split vent trap for the Split/Splitless Capillary (SSL) or Multi-Mode Inlet (MMI) using the built-in procedure accessed from Agilent 9000 touch screen display or web interface.
- ☒ If the GC includes a Flame Ionization Detector (FID), replace the Jet. If the Ignitor shows any buildup of sample or corrosion, replace the Ignitor. Examine the FID collector and castle assemblies for contamination – clean as necessary. Use the built-in procedures accessed from Agilent 9000 touch screen display or web interface.
- ☒ Replace the Guard Chip or Jumper Chip for the Split/Splitless Capillary (SSL) or Multi-Mode Inlet (MMI) using the built-in procedure accessed from Agilent 9000 touch screen display or web interface.

Inlet and Detector Tests

- ☒ Zero all pressure sensors.
- ☒ Perform the Inlet pressure decay test.
- ☒ Perform the Inlet leak/restriction test.
- ☒ Perform the FID Jet restriction test if FID installed.
- ☒ Perform the FID leakage test if FID installed.
- ☒ Record if test passed or failed in the results table.

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.7597222222
© Agilent Technologies, Inc. 2023

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ALS Maintenance

- ☒ Check all cabling and configuration settings between GC, tray, and injectors.
- ☒ Vacuum or remove any dust, especially around fans.
- ☒ Check operation of all fans.
- ☒ Check syringe for smooth plunger operation.
- ☒ Check for smooth operation of the needle support rod – clean if necessary
- ☐ Check for correct operation of syringe volume stops

Restore Instrument

- ☒ Restore the normal operating conditions using the Keyboard or Data System.
- ☒ Check and record the post PM detector signal output values. Results should be similar or lower than the detector outputs recorded prior to PM.
- ☒ Perform a chemical checkout. If this is a routine PM, inject the customer's sample using the ALS if applicable. This will act as a final checkout of both the ALS and the GC.

Disclaimer:
If the PM Service is performed prior to a qualification service, then use the qualification procedure as a guide for final instrument set up and checkout.

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.7597222222
© Agilent Technologies, Inc. 2023

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Signature Page

Service Review

- ☒ Attach available reports/printouts of all tests to this documentation.
- ☒ Record the Preventive Maintenance service activity in the customer's records/logbook.
- ☒ Update/reset Instrument maintenance counters as appropriate.
- ☒ Affix the PM sticker to the system or Instrument logbook based on the customer's request.
- ☐ Complete the Service Engineer Comments section if there are additional comments.
- ☒ Review with the customer this service, parts replaced, and test results obtained.
- ☐ If the Instrument firmware was updated, record the details of the change in the Service Engineer's Comments box or if necessary, in the customer's IQ records.
- ☐ Supply the customer with a copy of the Smart Alerts flyer.
- ☐ Describe Smart Alerts to the customer.
- ☐ Install Smart Alerts if requested.

PM Test Results Table

Detector Signal Outputs	Before PM Service	After PM Service
Detector output [D1] FID	30	18
Detector output [D2]		

Tests	Expected Result	Actual Result or N/A
Inlet Leak Test	Pass	Pass
Inlet Restriction Test	Pass	Pass
FID Jet restriction test if FID installed	Pass	Pass
FID leakage test if FID installed	Pass	Pass

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.7597222222
© Agilent Technologies, Inc. 2023

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Intuvo Parts List Table

Note: The following kits are recommended for capillary and MultiMode inlets. If this is a general PM and the customer has a preferred set of consumables, you may use the customer's consumables.

Part Description	Part Number	Product/Model # where used	Quantity Consumed
FID Jet 0.11inch ID	5200-0176	G3950A	1
Inlet PM Kit	5189-6497	G3950A	1
FID Ignitor Glow Plug	19231-60660	G3950A	1
Bus Bolt with Washer	G4581-60260	G3950A	1
Guard Chip for SS Inlet	G4587-60565	G3950A	1
Guard Chip for MMI	G4587-60665	G3950A	1
Jumper Chip for SS Inlet	G4587-60575	G3950A	1
Jumper Chip for MMI	G4587-60675	G3950A	1
Column Compression Bolts	G4581-60260	G3950A	1
Split Vent Trap Filter (2pk)	5188-6497	G3950A	1

Service Engineer Comments


If there are any specific points you wish to note as part of performing the service or other items of interest for the customer, please write include them in this box.

Service Completion

Service request number 6007559675 Date service completed 03 Mar 2025
Agilent signature [Signature] Customer signature _____
Total number of pages in this document _____

Revision: 2.04, Issued: May 1st, 2024
Agile Document Number: D0013614
DE number: 44166.7597222222
© Agilent Technologies, Inc. 2023

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Revision: 2.04, Issued: May 1st, 2024
 Agile Document Number: D0013614
 DE number: 44166.759722222
 © Agilent Technologies, Inc. 2023

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No.	Instrument/Equipment	Manufacturer	Model/Serial No.	Calibration	Calibration No.	Date of Calibration	Due Date of Calibration
1	Agilent 7000 Series GC System	Agilent	7000 Series GC System	Agilent Calibration Services	2503-0437OC-3	1 Mar 25	7 Mar 25

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
 534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL.0-2717-3000-29 FAX.0-2719-9484

Certificate of Calibration

Cert. No.: 25TM579
 Page : 1 of 3

Equipment : Hot Air Oven
 Manufacturer : Memmert
 Model : UF 55
 Serial No. : B212.0411
 ID No. : UAE.WAO.005/2556
 Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
 3 Soi Udornsuk 41, Sukhumvit Road,
 Bangchak, Phrakhanong,
 Bangkok 10260
 Location : Lab Floor 2
 Received Order : 19 March 2025
 Calibration Date : 19 March 2025
 Ambient Temperature : (26 ± 10) °C
 Relative Humidity : (50 ± 30) %
 AC Line Voltage : (220 ± 22) V
 Calibrated by : Man Pattanapongpaiboon
 Approved by : Kunchit
 Approved Signatory
 () Chakrit Waewwanjua
 () Suwit Imjai
 (✓) Kunchit Promprat
 Issue Date : 27 March 2025

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This certificate may not be reproduced other than in full, except with the prior written
 Approval of the head of Corporate Services 3, Equipment Calibration and Testing Services.

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Equipment : Hot Air Oven
 Condition As-Received : Used Item
 Reference : 2503-0437OC-3
 Procedure Used :-

Cert. No.: 25TM579
 Page : 2 of 3

Calibration were conducted using calibration procedure CP-OT02 based on TLAS G-20 according to direct measurement method with Data Acquisition which connected with Resistance Temperature Detector (RTD) and Thermocouple Type T.

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- Reference standard instrument-

Instrument	Serial No.	Cert. No.	Traceable	Due Date
1) Data Acquisition	KY44073381	24LM73	TPA	18 May 2025
- This certificate is valid only to the item calibrated on date and place of calibration.
- This certification is traceable to the International System of Unit.

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
Function of UUC* : Temperature Source

Fresh air setting : Close

Environment during calibration		
	Beginning	Finished
Temp. (°C)	27	28
REL Humid. (%)	49	55
AC Supply (Volt)	221	224

Probe Installation Details : Dimension of Chamber :
 a = 5.0 cm D = 0.50 m
 b = 5.0 cm W = 0.80 m
 c = 5.0 cm H = 0.75 m
 Capacity = 0.30 m³

Ref. Std. ID No.: @ Calibration Point		
Position :	(120,180) °C	(104) °C
1	23-01TC-01	1RTD-2/1
2	23-01TC-02	1RTD-2/2
3	23-01TC-03	22-01RTD-03
4	23-01TC-04	1RTD-2/4
5	23-01TC-05	1RTD-2/5
6	23-01TC-06	1RTD-2/6
7	23-01TC-07	23-01RTD-07
8	23-01TC-08	1RTD-2/8
9 (ref.)	23-01TC-09	23-01RTD-09

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Equipment : Hot Air Oven
Condition As-Received : Used Item
Reference : 2503-04370C-3
Result of Calibration : (*) Without Adjustment
Function of UUC* : Temperature Source
Fresh air setting : Close

Cert. No.: 25TM579
Page : 3 of 3

Calibration Point (°C)	UUC* Setting (°C)	UUC* Reading (°C)	Temperature stability (± °C)	Temperature uniformity (°C)	Overall Variation (°C)	Coverage Factor k
104.0	104.0	104.0	0.040	0.43	0.78	2
120.0	120.0	120.0	0.84	1.3	1.6	2
180.0	180.0	180.0	0.49	1.5	1.8	2

Calibration Point (°C)	Measured Temperature (°C)									Uncertainty (± °C)
	1	2	3	4	5	6	7	8	9 (ref.)	
104.0	104.335	104.135	104.363	104.317	103.849	103.738	104.179	104.229	104.025	0.42
120.0	119.575	119.366	119.807	119.905	118.994	119.194	119.888	119.894	120.064	1.1
180.0	180.286	179.510	180.401	180.551	179.281	179.463	180.196	180.451	180.374	1.2

Average* : The average of 30 values in each position.
Temperature stability : One-half of the greatest maximum difference of measured temperature at any one sensor.
Temperature uniformity : The maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time or at as close an observation time as possible to determine the temperature pattern or homogeneity within the chamber under steady-state conditions.
Overall Variation : The Difference of the maximum and minimum measured temperatures throughout observation.
UUC* : Unit Under Calibration
Note : The reported uncertainty of measurement was included stability and excluded uniformity .

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Certificate of Testing

Cert.No.: 25TW29
Page: 1 of 2

Equipment : DO Meter
Manufacturer : YSI
Model : 5100
Serial No. : 11B 101863
ID No. : UAE.WAO.004/2554
Received Date : 14 February 2025
Test Date : 17 February 2025
Reference : 2502-0473DSC-1
Submitted by : United Analyst and Engineering Consultant Co., Ltd.
3 Sol Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260
Laboratory Condition : Temperature (25 ± 5) °C
Humidity (50 ± 20) %
Test Procedure : In - house method : CP-CHS
by Comparison Technique with Azide Modification Method
Tested by : Walailak Sirithan
Approved by :
Approved Signatory
() Chakrit Waewwanjua
() Ponpan Palpin
(✓) Sathip Meangmal
Issue Date : 18 February 2025

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1. Reference Standard Instruments :
This certification is traceable to the International System of Unit through the reference standards laboratory of Industrial Calibration Center, Technology Promotion Association (Thailand-Japan).

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Result : Dissolved Oxygen Meter Adjustment With Air 100 %
Dissolved Oxygen Probe No.: 24F100202

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Code No.: BL002-25

Page: 1 of 3

Customer Name: United Analyst and Engineering Consultant Co., Ltd.
Address: 3 Soi Udomsuk 41, Sukhumvit Rd., Bang Chak, Phrakhanong, Bangkok 10260

Equipment: Electronic Balance
Manufacturer: Mettler Toledo
Model: AB204-S/FACT
Serial No.: 1129361010
Asset No.: UAE.WAS.002/2552

Building: N/A Floor: 1 Room: 107

Received Date: April 22, 2025

Date of Calibration: April 23, 2025

Calibration Conditions: Temperature 22.8 °C to 23.4 °C
Humidity 54.8 % to 68.9 %
Pressure 756.6 mmHg to 758.2 mmHg

Calibrated by: Sakkarin Srirahang

Approved by: Suwit Chotnok

Signature: *[Signature]*

Issued Date: April 25, 2025

Note: 1) The Uncertainties are for a confidence probability of approximately 95%

2) This Certificate is valid only to the item calibrated on date and place of calibration.

3) This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation Scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the units of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the United Analyst and Engineering Consultant Co., Ltd. (UAE)

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Code No.: BL002-25

Page: 2 of 3

Equipment: Electronic Balance
Model: AB204-S/FACT
Serial No.: 1129361010
Max. Capacity: 220 g
Calibration Date: April 23, 2025
Condition As-Received: In Condition

Manufacturer: Mettler Toledo
Readability: 0.0001 g
ID No.: UAE.WAS.002/2552

Condition of Equipment:

Condition of This Result of Calibration:

1. Calibration Method: This instrument was calibrated by metro UAE/CP/CAL 006 In-House Method based on UKAS Lab 14 : 2022

2. Reference Standards:

Reference Standard:	Model	Serial No.	Calibrated By	Certificate No.	Traceability	Due Date
Standard Weight Class E2 (IDM/L)	1 mg to 1 kg	8799109122	AMARC	25 009355	Mettler Toledo	21-Jan-27
Standard Weight Class F1 (IDM/L)	1 mg to 200 g	11119512	AMARC	24 013840	Mettler Toledo	04-Feb-26
Instrument	Model	Serial No.	Calibrated By	Certificate No.	Traceability	Due Date
Thermo-Hydro-Baro Meter	MHB-3825D	AK-46457	SUCLLS	50-H-00997107	Sakkarin Srirahang	21-Nov-25
Thermo-Hydro-Baro Meter	MHB-3825D	AK-46457	TPA	25P755	TPA	25-Feb-26

3. This certification is traceable to SI Unit

4. This certification was certified only for the instrument we calibrated

5. This result of calibration was found accurate as shown on date and place of calibration only.

6. Through the reference standard laboratory of AMARC 25-009355 Calibration 0152

Calibration Result:




1. Repeatability of Reading:

Nominal Value (g)	Standard Deviation of Reading (g)
200*	0.000045

2. Eccentric or off-center loading

A mass of 100 g was placed and moved to various position on pan

The Balance reading obtained is given in the table.

1 (g)	2 (g)	3 (g)	4 (g)	5 (g)	Maximum Difference (g)
100.0000	99.9996	99.9997	100.0003	100.0005	0.0005

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Code No.: BL002-25

Page: 3 of 3

Equipment: Electronic Balance
Model: AB204-S/FACT
Serial No.: 1129361010
Max. Capacity: 220 g
Calibration Date: April 23, 2025

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3. Error of Indication from nominal or conventional mass value:

Nominal Value (g)	Reference Value (g)	Indication (g)	Correction (g)	Uncertainty ($\pm$ mg)	Coverage Factor k
Unloaded	0.0000000	0.0000	0.0000	0.10	2.05
0.01	0.0100025	0.0099	0.0001	0.10	2.05
0.05	0.0500056	0.0500	0.0000	0.10	2.05
0.1	0.1000012	0.0999	0.0001	0.10	2.05
0.5	0.5000133	0.5000	0.0000	0.10	2.05
1	1.0000105	1.0000	0.0000	0.10	2.05
10	10.000010	10.0000	0.0000	0.11	2.04
40	40.000076	40.0000	0.0000	0.14	2.00
50	50.000056	50.0000	0.0001	0.15	2.00
80	80.000107	80.0000	0.0001	0.18	2.00
100	100.000109	99.9999	0.0002	0.17	2.00
120	120.00015	119.9999	0.0003	0.21	2.00
150	150.000165	149.9998	0.0003	0.24	2.00
160	160.000175	159.9997	0.0005	0.26	2.00
200	200.000129	199.9998	0.0004	0.30	2.00

4. Effect of Tare test:

Tare Load (g)	Test Load (g)	Indication (g)	Correction (g)
100	20.000041	19.9999	0.0001
	40.000076	39.9998	0.0002
	60.000056	59.9997	0.0003
	80.000107	79.9999	0.0002
	100.000168	100.0004	-0.0003

Remark:

The report uncertainty of measurement was based on standard uncertainty multiplied by coverage factor k, providing a level of confidence of approximately 95%.

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Calibration Certificate

Certificate No.: 2502226-001-01
Client name: UNITED ANALYST AND ENGINEERING CONSULTANT CO., LTD.
Address: 3 Soi Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260

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Approved by *[Signature]*
(Mr.Pheraphat Tuanjit)

Manager, Division of Calibration Laboratory
Responsible for the Technical Management Team

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ICS:09 Revision: 01 Date: 20-04-63

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR20SDU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

Date of Calibration: 20 March 2025 Page 2 of 4

Environment Condition: Ambient Temperature: 21.2 ± 0.6 °C Relative Humidity: 48 ± 3.5 %

Place of Calibration: 208 Balance Room, UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.

Condition of Equipment: Good Condition

Condition of This Results of Calibration:

1. Calibration Method: NFI Method W-MA-001 In-House Method based on UKAS Lab 14 : 2019

2. Reference Standards:

Reference Standard	Model	Serial No.	Calibrated By	Certificate No.	Due Date
Standard Weight Class E2	1mg to 200g	8505567572	TCS	M24341005	19 April 2025
Instrument	Model	Serial No.	Calibrated By	Certificate No.	Due Date
Thermo-Hygro Meter	608-H1	N-1.81H 01/7/23	Quakey Reborn	QR25-0y42	10 February 2026

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EC-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR20SDU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

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3. Departure from Nominal Value: (Range: 0 - 82 g ; Resolution: 0.00001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (±g)	Coverage Factor k
Unloaded	0.00000	0.00000	0.00000	0.0000089	2.00
0.001	0.001003	0.00100	0.00000	0.0000092	2.00
0.005	0.005002	0.00500	0.00000	0.0000094	2.00
0.01	0.010003	0.01000	0.00000	0.0000091	2.00
0.05	0.049996	0.05000	0.00000	0.0000098	2.00
0.1	0.100011	0.10000	0.00001	0.000011	2.00
0.5	0.500016	0.50000	0.00002	0.000014	2.00
1	1.000003	1.00001	-0.00001	0.000016	2.00
2	2.000023	2.00005	-0.00003	0.000017	2.00
5	5.000016	5.00005	-0.00003	0.000021	2.00
10	10.00009	10.00005	-0.00004	0.000026	2.00
20	20.000030	20.00012	-0.00009	0.000037	2.00
30	30.000039	30.00012	-0.00008	0.000050	2.00
50	50.000078	50.00014	-0.00011	0.000068	2.00
80	80.000067	80.00020	-0.00013	0.00011	2.00

EC-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR20SDU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

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3. Departure from Nominal Value: (Range: >80 - 200 g ; Resolution: 0.0001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (±g)	Coverage Factor k
90	90.00015	90.0002	-0.00007	0.000015	2.00
100	100.00006	100.0001	0.00004	0.000016	2.00
110	110.00007	110.0001	0.00003	0.000017	2.00
120	120.00009	120.0002	-0.00001	0.000018	2.00
130	130.00010	130.0002	-0.00001	0.000019	2.00
140	140.00013	140.0002	-0.00001	0.000019	2.00
150	150.00009	150.0002	-0.00001	0.000021	2.00
160	160.00010	160.0002	-0.00001	0.000022	2.00
170	170.00012	170.0002	-0.00001	0.000023	2.00
200	200.00012	200.0002	-0.00001	0.000028	2.00

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EC-012 Revision: 01 Date: 20-04-65

Calibration Certificate

Certificate No.: 2502226-002-01
Client name: UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Address: 3 Soi Udomsuk 41, Sukhumvit Road, Bangchack, Prakhnanong, Bangkok 10260

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Approved by *N. Nigrobat*
(Mr.Pheraphat Tuanjit)
Manager, Division of Calibration Laboratory
Responsible for the Technical Management Team

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Page : 1 of 2

CERTIFICATE OF CALIBRATION

Equipment : COD Test Tube Heater
Meter Model : HI839800-02 **Serial No. :** H0185001
Tube Heater : 25 Vial Capacity **Resolution :** 0.1°C
Temperature Range : (-10 to 160)°C **Temperature of Reaction :** 150°C
Manufacturer : Hanna Instruments **Made in :** Romania
Condition As-Received : Used Product **Reference :** RE250401
Ambient Temperature : (25 ± 2)°C **Relative Humidity :** (50 ± 15) % RH
Customer name : United Analyst and Engineering Consultant Co., Ltd.
 3 Soi Udomsuk 41, Sukhumvit Rd., Bangchak,
 Phrakhanong, Bangkok 10260
Received date : 5 March 2025
Calibrate date : 7 March 2025
Issue date : 7 March 2025
Calibrated Location : Hanna Instruments (Thailand) Ltd.
Calibration Procedure : This calibrator was conducted by using in-house calibration procedure
 CP-04 by using certified reference standard instruments.

Calibrated by : ☒ Mr. Pichit Petthong
☐ Mr. Channarong Soinak

Approved by :
 Mr. Anan Suwanchaisakul
 Authorized Signatory

This certificate was certified only for the instrument we calibrated.

This result of calibration was found accurate on date and place of calibration only.

** This certificate may not be reproduced other than in full, except with the prior written **

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Reference Standard Instruments : This certification is traceable to the international unit of unit maintained through:

Instruments	Model	Serial No.	Certificate No.	Traceable
Data Acquisition Switch Unit	34970A	MY44065265	WK2407-141-1	WK Electric Co., Ltd.
Digital Thermo-Hygrometer	HT-771SD	A1.07155	25H171	Technology Promotion Association (Thailand-Japan).

Calibration Result:

Measurement Temperature Source Accuracy for COD Reactor.

Capacity (Vial)	Nominal Value (°C)	Average Value (°C)	Uncertainty of Measurement (±°C)
25 Vial	150.0	150.4	0.47

Unit : °C

(1A)	(2A)	(3A)	(4A)	(5A)
150.407	150.377	150.269	150.402	150.422
(1B)	(2B)	(3B)	(4B)	(5B)
150.426	150.394	150.644	150.690	150.542
(1C)	(2C)	(3C)	(4C)	(5C)
150.477	150.303	150.627	150.257	150.176
(1D)	(2D)	(3D)	(4D)	(5D)
150.462	150.456	150.199	150.406	150.102
(1E)	(2E)	(3E)	(4E)	(5E)
150.185	150.513	150.235	150.460	150.442

Figure: Shows the location of the temperature source.

The report uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor $k = 2$, providing a level of confidence of approximately 95%

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Note: While non-current production AA instrument and or accessory models are not covered specifically in this document it can be used as a basic reference.

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- Contact the customer and ensure that all necessary supplies are available before the preventive maintenance visit.
- Confirm the ability of the instrument to deliver continued safe operation as established via the Agilent AA safe operation flow chart. (Refer directly to the AA 55/240/280 Preventive Maintenance Scope of Work to make this decision.)
- Only select those pages that relate to the system or module being serviced.
- Complete empty fields with the relevant information.
- Complete the relevant checkboxes in the checklist using either a "X" or tick mark "✓".
- Check "Section not applicable" check boxes to indicate services/tasks not delivered, as appropriate.
- Complete the Preventive Maintenance service in the order of the tasks listed.
- Complete the Service Review section together with the customer.
- Complete the fields for page numbers at the foot of each selected page
- Complete the total number of pages field in the Service Completion section
- Ask the customer to sign the Service Completion section including the customer's and your signature.

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Instrument System Name and ID	240 FS RRS
Instrument System Site and Location	United Analyst and Engineering Consultant

List System Component Product Numbers	List the Serial Numbers of each Component
1. G 9432 R	M1 13160001
2.	
3.	
4.	
5.	
6.	
7.	
8.	
9.	

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- ☒ Discuss any specific issues with the customer before starting.
- ☐ For HF application systems, if standard sample introduction system was not installed, ask the customer to install it. หมายเหตุ
- ☒ Review the instrument logbook for recorded problems and comments.
- ☒ Save instrument control settings before starting the procedure.
- ☒ Perform a general inspection of the system for cleanliness.
- ☒ Check for proper installation of parts, assemblies, sensors etc.
- ☒ Check system for required installation of components, settings as defined by current Service Notes
- ☒ Check for required firmware updates and verify with customers if they would like them installed.
- ☒ Use SVD to perform a Full Wavelength Scan for Cu HCL - "As found test_1"
- ☒ Perform a Basic Cu ABS test - "As found test_2"
- ☒ Print the Details page or screen captures of the test results and attach to the end of this checklist.

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- ☐ Section not applicable

Electronic components

- ☒ Review and confirm instrument configuration data in SVD
- ☒ Confirm power supply voltages using the SVD Power Supply diagnostic.
- ☒ For Dual Beam Instruments - Confirm RBC frequency using the SVD RBC frequency diagnostic.

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Optics components

- ☒ Check that external optical surfaces are clean – Clean or replace as required
- ☒ Use SVD and perform Mono Wavelength Correction.
- ☒ Use SVD and perform Slit Calibration.
- ☒ Use SVD and perform Grating Squareness Diagnostic.
- ☒ Use SVD and perform Zero Order Offset/Mono Correction.
- ☒ Use SVD and perform Wavelength Repeatability.
- ☒ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.
- ☒ Check that the signal energy of the D2 and HC lamps track properly. Advise customer if their D2 lamp is showing emission degradation due to age.

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- ☒ Inspect the burner interlock plate to ensure that the interlock pin is secure and correct for the burner type.
- ☒ Clean the burner slot with a clean white card.
- ☒ Check the uniformity of the slot width.
- ☒ Clean the burner if required.
- ☒ Change the burner o-ring.
- ☒ Clean the nebulizer, spray chamber and liquid trap.
- ☒ Change all o-rings and seals in the nebulizer, nebulizer block and spray chamber.
- ☒ Check that the pressure relief bung releases readily.
- ☒ Change o-rings on the fuel and oxidant delivery bars.
- ☒ Leave the liquid trap EMPTY and verify the flame will not ignite in this state.
- ☒ Refill liquid trap and check that overflow drains freely into the drain/waste tube.
- ☒ Check the drain/waste tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel.
- ☒ Check and clean the igniter electrode.

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- ☒ Pressure test for leaks.
- ☒ Leak test gasbox internal components and connections.
- ☒ Check safety interlock status and operation using the *SVD Interlock monitoring diagnostic*.

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- ☐ Inspect the GTA workhead gas hoses and connections for leaks.
- ☐ Pressure test for gas leaks
- ☐ If the cooler system is accessible (stand-alone) check for correct operation and coolant/water level – this includes any temperature and pressure settings plus filter cleaning (air flow and water).
- ☐ Inspect the GTA workhead water hoses and connections for leaks.
- ☐ Check all graphite components and replace if necessary.

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- ☒ Section NOT Applicable
- ☐ Check condition of the PSD capillary – replace if necessary
- ☐ Check condition and operation of PSD syringe – ensure it does not have air locks and bubbles.
- ☐ Change PSD rinse bottle o-ring
- ☐ Check and clean the rinse vessel
- ☐ Check the drain tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel.
- ☐ Ensure that the waste vessel is suitable for use with the furnace system.

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- ☒ Section NOT Applicable
- ☐ Re-torque screws securing the hubs, presser arms and pump rotors.
- ☐ Adjust each roller so that it rotates freely.
- ☐ Wipe clean the pump rotor rollers and pump bands with a dry clean cloth.
- ☐ Ensure that the presser arms and the surfaces near the pump are free from dirt and spills.
- ☐ Remove the pump module rear cover and check for the incursion of liquids and any signs of corrosion.
- ☐ Re-torque the nuts that fasten the motor mounting plates to the chassis
- ☐ Check clips securing the diluents holder and replace if necessary.
- ☐ Disconnect, clean T piece, and reassemble the tubing using the following steps.

Revision: 13.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

- ☐ Remove the T-piece by disconnecting the pump tubes, the pump bands and all other tubing.
- ☐ Place the T-piece in an ultrasonic bath containing strong detergent 1-5% Decon 30 or similar, for approximately 5-10 minutes.
- ☐ Wash the T-piece under a tap with a strong flow of water
- ☐ Rinse with distilled water through all of the inlets in the reverse direction to normal sample flow.
- ☐ Reassemble.

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- ☐ Cleaning the spill tray, rack location mat, end frames and chassis's accessories with a damp soft cloth and diluted mild detergent.
- ☐ Cleaning the autosampler cover panels with domestic window cleaner.
- ☐ Checking the X- axis and Z- axis drive belts for cracks, splits, damaged teeth, excessive fraying, color changes or degradation from fumes..
- ☐ Check the X- axis, Theta- axis and Z- axis FFC cables for cracks, incorrect positioning, damaged edge or damaged connectors.

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- ☒ Section NOT Applicable
- ☐ Check the x-axis and z-axis timing belts – Replace if there is any cracks, splits or color deterioration and belt tension.
- ☐ Check belt tensions - adjust if required
- ☐ Check the lubrication pad for single x-axis shaft. If pad is dry or customer has observed any vibration or erratic movements of the x-axis carriage, add 1 mL of Dow Corning 200 © Fluid, 200 CS into the well.
- ☐ Check the auto-sampler ability to find tube positions - Calibrate if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Revision: 13.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

Vapor generation accessory VGA (hydride generator)

- ☐ Section NOT Applicable
- ☐ Inspect VGA gas supply hose.
- ☐ Inspect/replace VGA pump tubing
- ☐ Check low gas pressure interlock setting – adjust if required.
- ☐ Check precision orifice gas flow setting – adjust if required.
- ☐ Check gas regulator pressure to 46 psi (325 kPa) – adjust if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

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Test Results

Test Description	Expected Test Result	Actual Test Result
Flame optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	49 /
Flame performance test with 5 ppm copper sample		
Air/acetylene, mixing paddle removed	Abs value > 0.5	0.5599
Air/acetylene, mixing paddle installed, 10 replicates	%RSD < 1.0	0.2 /
Deuterium furnace optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	-
Deuterium furnace performance test with 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0%	-
Abs value	≥ 0.15	-
Zeeman furnace analytical performance: 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0%	-
Abs value	≥ 0.10	-
MSR%	≥ 70 %	-

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

AA consumable and parts list table

Part Description	Part Number	Product/Model # where used	PM supplied or Consumable	Instrument Type
Test Solution - Cu 5ppm solution	6610030100	50 55 140 243 280	PM supplied	Common
Test Solution - Blank solution	5190-0001	50 55 140 243 280	PM supplied	Common
Copper, 1000 ug/ml, 100ml	5190-8279	50 55 140 243 280	*	Common
Kit, Nk 70 rings, aqueous, complete set	9910093400	50 55 140 243 280	PM supplied	Flame
Organic Kit	9910093500	50 55 140 243 280	PM supplied	Flame
Wire Nebulizer Cleaning	9910024700	50 55 140 243 280	consumable	Flame
Tubing Capillary Std Nebs	9910024800	50 55 140 243 280	consumable	Flame
Capillary Tube Hvac Neb (3) (organics only)	9910044000	50 55 140 243 280	consumable	Flame
Glass impact beads (5/pkg)	9910025700	50 55 140 243 280	consumable	Flame
Teflon impact beads (5/pkg) (organics only)	9910053300	50 55 140 243 280	consumable	Flame
Burner cleaning strip (100/pkg)	9910053900	50 55 140 243 280	consumable	Flame
Window UV silica - round (right side)	2010082600	50 55 140 243 280	PM supplied	Common
Window UV silica - rectangular (left side)	2010082500	50 55 140 243 280	PM supplied	Common
Pad adhesive window (round)	2010012700	50 55 140 243 280	PM supplied	Common
Pad adhesive window (rectangular)	2010012800	50 55 140 243 280	PM supplied	Common
Electrode kit (1 pr) (D2)	6310063400	G1A120	PM supplied	Furnace
Shroud (D2)	6310063100	G1A120	PM supplied	Furnace
Zeeman electrode kit (1 pr)	6310063500	G1A120	PM supplied	Furnace
Zeeman shroud	6310063600	G1A120	PM supplied	Furnace
O-ring, PSD rinse bottle	6910025900	PSD120	PM supplied	Furnace

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Items classified as PM supplied in the above table are included in the standard PM

Those classified as consumable should be provided by the customer or charged to the customer if supplied by the Agilent service engineer.

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Placeholder for Service Engineer Comments (optional)

Service Completion

Service request number 6003949143 Date service completed 30 Jan 2025

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Total number of pages in this document 13

Revision: 12.00, Issued November 2021

© Agilent Technologies, Inc. 2021

SVD Results Report

Report ID: 4	Diagnostic Start Time: 1/30/2025 9:47:25 AM	Diagnostic End Time: 1/30/2025 9:47:25 AM
Customer: UAE	Service Engineer: Kanyakorn S.	
Address: Soi Udomsuk 41, Sukhumvit Rd. Bangkok	Contact Details: 025376363#1	

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Averaging Period: 30.0

Datapoint Count: 20

Upper Limit:
91.00Highest Measured Frequency:
50.00Average Frequency:
50.00Lower Limit:
49.00Lowest Measured Frequency:
50.00Result: **Passed**

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

Grating: S

Lamp Element(s): Copper
Lamp Socket Function: 3
Lamp Current(mA): 4.00
Slit Width(nm): 0.5
1st Order Wavelength(nm): 324.80
2nd Order Wavelength(nm): 649.60

	Lower Limit (nm)	Actual (nm)	Upper Limit (nm)	Result:
Zero Order	649.60	649.60	649.60	Passed
First Order	324.80	324.80	324.80	Passed
Second Order	162.40	162.40	162.40	Passed

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Wavelength Report: 12V

Lamp Used: Copper
Peak Used(nm): 324.80
Connected to Socket: 3

Lamp Current(mA): 4
Slit Width(nm): 0.2
Slit Height: Normal

Lamp Alignment: **Performed**

Lower Limit(nm)	324.80	Upper Limit(nm)
Sample 1: 324.823		
Sample 2: 324.823		
Sample 3: 324.823		
Sample 4: 324.823		
Sample 5: 324.819		
Sample 6: 324.819		
Sample 7: 324.819		
Sample 8: 324.819		
Sample 9: 324.823		
Sample 10: 324.819		

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Interpret:

Sample Flow: ☒ Working
 H2O Sample Flow: ☒ Working
 Flame Control: ☒ Working
 Gas Control: ☒ Working
 Pressure Control: ☒ Working
 Liquid Control: ☒ Working
 Flame Detect: ☒ Working
 GCU Active: ☒ Working
 Oxidant Pressure: ☒ Working
 Gas Vent Changeover: ☒ Working
 Ignition: ☒ Working

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Auto Lab

Lamp 1: ☒ Working
 Lamp 2: ☒ Working
 Lamp 3: ☒ Working
 Lamp 4: ☒ Working
 Lamp 5: ☒ Working
 Lamp 6: ☒ Working
 Lamp 7: ☒ Working
 Lamp 8: ☒ Working

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
1/30/2025 10:53 AM
Page 1 of 1

SpectrAA

Analyst:
 Date Started: 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
 Worksheet: Sensitivity Test 01
 Comment:
 Methods: Cu
 Computer name: DESKTOP-R9UFRS
 Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL 21910	0.000	38.8	0.0002
Readings	0.0002	0.0003	0.0001
STANDARD 1	5.000	0.1	0.5571
Readings	0.5571	0.5563	0.5575

Curve Fit: Linear Origin
 Characteristic Conc: 0.035 mg/L
 r: 1.0000
 Calculated Conc: 0.002 ± 0.009
 Residuals: -0.007 ± 0.009

Abs = 0.11141 × C

5 ppm Cu	5.000	0.3	0.5568
Readings	0.5512	0.5596	0.5615
			1/30/2025 10:52:54 AM

Sequential by time report


1/30/2025 10:48 AM
Page 1 of 1

SpectrAA

Analyst:
 Date Started: 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
 Worksheet: Precision Test
 Comment:
 Methods: Cu
 Computer name: DESKTOP-R9UFRS
 Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL 21910	0.000	64.1	-0.0002
Readings	-0.0003	-0.0003	-0.0001
STANDARD 1	5.000	0.3	0.6052
Readings	0.6035	0.6073	0.6047


Curve Fit: Linear Origin
 Characteristic Conc: 0.12105 mg/L
 r: 1.0000
 Calculated Conc: -0.002 ± 0.009
 Residuals: 0.002 ± 0.009

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PinAAcle 900F Preventive Maintenance Report

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PinAAcle 900F Preventive Maintenance (PM)

Company Name:	United Analyst and Engineering Consultant Co., LTD.		
Address (Instrument Location):	41 Sukumvit Rd., Phra Khanong, Bangkok 10260		
Serial Number:	PFBS20031902	PM Number:	1 of 2
Customer Name (if applicable):	K. Yinda	Telephone Number:	095-5580049
Customer Support Engineer Name:	K. Chayanan	Service Order Number:	WO-03126047
Date PM Performed: (DD-MMM-YYYY)	29-Apr-2025	Next PM Due Date: (DD-MMM-YYYY)	29-Oct-2025
Standard Labor Hours to Complete PM :		5 hours	

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Copyright Information

This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc.

Copyright © 2013 PerkinElmer, Inc.

Trademarks

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law. PerkinElmer is a registered trademark of PerkinElmer, Inc. All other trademarks and registered trademarks not owned by PerkinElmer, Inc. or its subsidiaries that are depicted herein are the property of their respective owners. Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no Warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

PerkinElmer shall not be liable for incidental or consequential damages in connection with the furnishing or use of this document.

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Parts Lists

Parts Included with the PM		
Part Number (if applicable)	Description	Quantity
B0501696	Fan Filters	2
N3160156	O-Ring Kits for Sampling Introduction (Stainless Steel Nebulizer)	N/A
N3160157	O-Ring Kits for Sampling Introduction (Plastic Nebulizer)	1
N9301714	Replacement Acetylene Filter Cartridge	1
TH001022	Replacement Air Filter Cartridge	1

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Procedure Checklist

Use [✓] to check off those steps in the checklist that have been completed.

- General:**
 - Review the Instrument performance with the customer and document any recent problems.
 - Inspect the customer log book and make any appropriate PM entries.
 - Perform general inspection of system for cleanliness.
- PC Instrument Software:**
 - Instrument Software user files/databases archived, packed, and/or deleted as needed.
- Mechanical:**
 - Inspect and clean all fans and filters. Replace filters if necessary.
 - Inspect all gas lines for leaks and/or wear. Replace if needed.
 - Clean exterior of the Instrument.
 - Inspect the burner head, burner chamber, and nebulizer. Clean if needed as stated in the Hardware Guide.
 - Check burner head dimensions with the feeler gauge as stated in the Hardware Guide in the Maintenance chapter section on cleaning the burner head and checking sloth width. Replace if out of specification.
 - Check the condition of the end cap, burner head, and nebulizer O-rings. Replace if necessary.
 - Check the drain system for signs of wear. Replace worn or damaged parts.
 - Visually check for proper flame conditions when lighting the Air-C₂H₂ and N₂O-C₂H₂ flames (if applicable).
- Electrical:**
 - Inspect PC boards. Clean if necessary.
 - Carefully check all internal and external cable connections.
 - Check instrument firmware revisions upgrade to current levels (if necessary)
 - Run Diagnostics Test within the Advanced function of the Spectrometer page. Check the results in the service log folder in the Spectrometer BM Log Viewer.
- Optics:**
 - Inspect and clean the sample compartment windows, if needed.
 - Inspect optics. Clean or replace if necessary.
- Gasses:**
 - Verify that the Gasses supplied to the instrument are within the pressure and purity specifications found in the PinAAcile 900 Series Pre-Installation Checklist SDB.
 - Verify that the acetylene filter and air filter element is dry. Replace if necessary.

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Parameter	Specification	Test Results	Pass/Fail
Flame Sensor	Air/C ₂ H ₂ Flame correctly shuts down	Active	Passed
Drain Sensor	Air/C ₂ H ₂ Flame correctly shuts down	Active	Passed
Nebulizer Sensor	Air/C ₂ H ₂ Flame correctly shuts down	Active	Passed
C ₂ H ₂ Pressure Sensor	Air/C ₂ H ₂ Flame correctly shuts down	Active	Passed
Air Pressure Sensor	Air/C ₂ H ₂ Flame correctly shuts down	Active	Passed
Burner Head Sensor	Choosing Nitrous Oxide as the oxidant should trigger an interlock shuts down	Active	Passed

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Parameter	Specification	Certificate Value at 533.6 nm (Abs.)	Test Results	Pass/Fail
1.0 A ND Filter	± 5% from Cert.	0.9668	0.9878	Passed
0.2 A ND Filter	± 5% from Cert.	0.1953	0.1876	Passed

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Description: Verifies the instrument's ability to compensate for Background absorption.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.010	0.009	Passed

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10. Review:

- Review with the customer PM work performed.
- Review with the customer routine maintenance procedures.
- Discuss recommended customer supplied materials to have on hand.
- Attach PM sticker.

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Review

The preventive maintenance checks and (if applicable) performance tests for PinAAcle 900F have been completed.

This PinAAcle 900F ☒ Passes ☐ Fails the preventive maintenance.

Review of Preventive Maintenance:

Authorized PerkinElmer Representative:	Date:
Chayman K.	29 Apr 2025
Authorized Customer Representative:	Date:
Chayman K.	29 Apr 2025

PinAAcle 900F Preventive Maintenance Report (PM)

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Work Description	
• PM 20 (Annual) • Cleaning Cell, Chamber, Filter • Wavelength Calibration : Pass • Wavelength Scan As Co./Ka/KM : Pass • Align cell with Hg : OK	
Start Date	End Date
29/04/2025	05/05/2025

Tools Used					
Quantity	Calibrated Tool	Description	Serial Number	Last Calibration Date	Next Calibration Date
*** No Calibrated Tools Used ***					

Materials Used				
Part Number	Part Description	Note	Lot/Serial Number	Quantity
*** No Parts Used ***				

Labor Details				
Part Number	Part Description	Start Date	Quantity	
SV000013	Preventive Maintenance	29/04/2025	4	

Work Complete	Customer Signature	Technician Signature
Yes ☒ No ☐		
PM/ODM/PLS with Customer		
Yes ☐ No ☒	Chayman K.	Chayman K.
	852568 K. Chayman (chayman@chayman.com)	852568 Kanin, Chayman

Terms & Conditions	

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Customer Acknowledgment of receipt of the above report / replacement.
Special Terms and Conditions: This is not an invoice.
Terms will be applied to your invoice if applicable.

Atomic Absorption/FIAS 100/400 Preventive Maintenance (PM)

Company Name:	United Analyst and Engineering Consultant Co., LTD.		
Address (Instrument Location):	41 Sukumvit Rd., Phra Khanong, Bangkok 10260		
Room Number:	Lab		
Asset Number (if applicable):	2 of 2W	Customer System ID:	K. Yanida
Service Engineer Name:	K. Chayman	Service Order Number:	WO-03051971
Date PM Performed: (DD-MMM-YYYY)	29-Apr-2025	Next PM Due Date: (DD-MMM-YYYY)	29-Oct-2025

Part Number	Release	Publication Date	PerkinElmer
09370005	C	January 2013	

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General Instructions:

Always check with the customer before making any changes that may affect the customer's analysis or calibration.
The completed document should be signed by an authorized PerkinElmer and customer representative and left with the customer.
Update the PM sticker and instrument logbook as required.

Copyright Information

This document contains proprietary information that is protected by copyright. All rights are reserved.
No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc. Copyright © 2013 PerkinElmer, Inc.

Trademarks

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law. PerkinElmer is a registered trademark of PerkinElmer, Inc. All other trademarks and registered trademarks not owned by PerkinElmer, Inc. or its subsidiaries that are depicted herein are the property of their respective owners.
Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.
PerkinElmer shall not be liable for incidental or consequential damages in connection with the furnishing or use of this document.

Atomic Absorption/FIAS 100/400 Preventive Maintenance

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Parts Lists

Parts Included with the PM				
Part Number (if applicable)	Description	Quantity	Batch/Lot #	Expiration Date (MM/YY)
B050 2706	Fan Filter	1	N/A	N/A

Additional Tools Required for PM				
Part Number (if applicable)	Description	Quantity	Serial #	Calibration Due Date (MM/YY)
N/A	Digital Volt Meter	1	N/A	N/A
Additional Reagents and Standards Required for PM				
Part Number (if applicable)	Description	Quantity	Batch/Lot #	Expiration Date (MM/YY)
N/A	N/A	N/A	N/A	N/A

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1. General:

- ☒ Review the instrument performance with the customer and document any recent problems.

- ☒ Is the Working Environment Acceptable? If not, document.

- ☐ Visual Damage (if yes, describe)

- ☒ Check incoming AC line voltage for proper levels and grounding.
☒ Verify Voltage switch on back of instrument is correct.
☒ Perform general inspection of system for cleanliness. Clean if needed.
☒ Gas supply cylinders secured, lines leak checked and argon or nitrogen supply pressure verified (45 – 58 psi).
☒ Inspect the customer log book and make any appropriate PM entries.
☒ Fan checked and filter cleaned
☒ Heating mantle or Universal Cell Holder checked

2. Instrument components

- ☒ Non-return valve checked/repaired/replaced if needed (B019 8111). Clean the valve if there is any liquid in it. Replace the rubber sleeve (B013 5123) if it is worn. Check the flow meter for any signs of fluid in it. Clean the flow meter if needed.
☒ Verify condition of pump pressure adjustment levers (B050 7794 - look for cracks or problems with the springs), pump rollers (B300 0251 check for wear), and thumb screws (B050 7796).
☒ Check the Multiport valve for proper switching, flow, and insure there are no leaks. Clean valve parts and replace o-rings if needed (large o-ring: B050 1250, small o-ring: B004 5095). Use a squirt bottle & fishing line to try to dislodge clogs.
☒ Firmware Version checked. Latest is 2.20.

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- ☒ Mixing separator assembly checked
☒ Filter/membrane checked (B050 8306)
☒ Condition of the pump tubing (replace if necessary), correct pump tubing for the solutions being run. Make sure the correct magazines are being used. B050 7791 for 0.13 – 1.80 mm tubing; B050 7792 for 1.60 – 3.18 mm tubing.

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4. Cell, Cell Windows, Transfer Line:

- ☒ Cell checked
☒ Cell windows checked
☒ Transfer line checked for moisture (if moisture is a problem, the Nafion dryer might be needed)

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6. Review:

- ☒ Review with the customer PM work performed.
☒ Review with the customer routine maintenance procedures.
☒ Discuss recommended customer-supplied materials to have on hand.
☒ Attach PM sticker.
☒ Update Logbook.

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Atomic Absorption FIAS 100/400 Preventive Maintenance Page 6 of 6

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Work Description	
- PM 52 Warranty - Cleaning Port Valve, Mandrel, Turret - Run fig test; Pass	
Start Date	End Date
28/04/2558	28/04/2558
Work Description	

Tools Used				
Quantity	Calibrated Tool	Description	Serial Number	Last Calibration Date
*** No Calibrated Tools Used ***				

Material Used				
Part Number	Part Description	Note	Lot/Serial Number	Quantity
*** No Parts Used ***				

Labour Details			
Part Number	Part Description	Start Date	Quantity
SV000013	Preventive maintenance	28/04/2558	1
SV000002	Service Travel	28/04/2558	1

Work Complete	Customer Signature	Technician Signature
Yes ☒ No ☐ PM/CM/PPV Left with Customer Yes ☐ No ☒	 85/2558 K. Nithakorn นนทิ์น	 85/2558 Kanan, Chayanan

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Customer Acknowledgment of receipt of the above repair / replacement.
Special Terms and Conditions: This is not an invoice.
Taxes will be applied to your invoice if applicable.

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Terms & Conditions
Customer Acknowledgment of receipt of the above repair / replacement. Special Terms and Conditions: This is not an invoice. Taxes will be applied to your invoice if applicable.

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CERTIFICATE OF CALIBRATION

Certificate No. : SP25-019 Page 1 of 5

Customer : United Analyst and Engineering Consultant Co., Ltd. (Head Office)

Address : 3 Soi Udomsuk 41, Sukhumvit Road, Hangehak, Phrakhanong, Bangkok 10260

Location of calibration : Instrument room (207)

Equipment : UV-Vis Spectrophotometer

Manufacturer : Agilent Technologies

Model : Cary 60

Serial No. : MY15410009

ID No. : UAE.WAT.020/2558

Received Date : 26 May 2025

Calibration Date : 26 May 2025

Issue Date : 29 May 2025

Condition Instrument : Good

Calibrated by : Mr. Tanawat Rittidach Approved by : Ms. Chonticha Sangsorn

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The measurement capability of the laboratory and its responsibility to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the DQE Services Co., Ltd.

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Traceability : This certification is traceable to the International System of Unit maintained at National -

Institute of Standards and Technology (NIST) through Starna Scientific Limited

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REPORT OF CALIBRATION

Certificate No. : SP25-019 Page 3 of 5

Calibration Results : Without adjustment

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Remark : - UUC = Unit Under Calibration
- N/A = Not Available
- The result expanded uncertainty of measurement U is stated as the standard uncertainty of measurement multiplied by the coverage factor k, which for a normal distribution corresponds to a coverage probability of approximately 95%

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CERTIFICATE OF CALIBRATION

Certificate No. : SP25-021 Page 1 of 5

Customer : United Analyst and Engineering Consultant Co., Ltd. (Head Office)

Address : 3 Soi Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260

Location of calibration : Instrument room (207)

Equipment : UV-Vis Spectrophotometer

Manufacturer : HITACHI

Model : U-2900

Serial No. : 21E22-009

ID No. : UAE.WAT.051/2564

Received Date : 26 May 2025

Calibration Date : 26 May 2025

Issue Date : 29 May 2025

Condition Instrument : Good

Calibrated by : ปณิฏ Approved by : จุฬจิรา
(Mr.Tanawat Ritidach) (Ms.Chonicha Sangnern)
Technical Manager Quality Manager

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REPORT OF CALIBRATION

Certificate No. : SP25-021 Page 2 of 5

Environment Condition : Ambient Temperature 25 ± 5 °C
Relative humidity 55 ± 20 %RH

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Traceability : This certification is traceable to the International System of Unit maintained at National Institute of Standards and Technology (NIST) through Starna Scientific Limited

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Remark : - UUC = Unit Under Calibration

- N/A = Not Available

- The result expanded uncertainty of measurement U is stated as the standard uncertainty of measurement multiplied by the coverage factor λ.

which for a normal distribution corresponds to a coverage probability of approximately 95%

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Certificate No. : SP25-024

Page 1 of 5

Customer : United Analyst and Engineering Consultant Co., Ltd. (Head Office)

Address : 3 Soi Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260

Location of calibration : Instrument room (207)

Equipment : UV-Vis Spectrophotometer

Manufacturer : HITACHI

Model : U-5100

Serial No. : 23A4-008

ID No. : UAE.WAS.010/2567

Received Date : 17 June 2025

Calibration Date : 17 June 2025

Issue Date : 20 June 2025

Condition Instrument : Good

Calibrated by : ธันว์
(Mr.Tanawat Ritidach)
Technical Manager

Approved by : ชลจิษฐ์
(Ms.Choniticha Sangngern)
Quality Manager

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Traceability : This certification is traceable to the International System of Unit maintained at National -

Institute of Standards and Technology (NIST) through Sarna Scientific Limited

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Remark : *UUC = Unit Under Calibration

*N/A = Not Available

*The result expanded uncertainty of measurement U is stated as the standard uncertainty of measurement multiplied by the coverage factor k.

which for a normal distribution corresponds to a coverage probability of approximately 95%

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WORK ORDER INFORMATION			
Top-Level		Order Type	Preventive Maintenance
Installed Product ID	IB-00105024	Billing Type	Chargeable
Product	SKALAR 2SAN59000 SAN++ Classic CFA 230V 2SAN59000	PO No.	HPO-250400209
Serial No.	182688	Warranty No.	
		Contract No.	

PRODUCTS SERVICED		
Installed Product Id	Serial Number	Product
IB-00105024	182688	SKALAR 2SAN59000 SAN++ Classic CFA 230V 2SAN59000

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Line Number	Engineer	Start Date And Time	End Date And Time	Billable Labor Hour	Billable Travel Hour	Travel KM
WL-00342192	Yongyuth Chanphong	05/23/2025 9:30 AM	05/23/2025 6:00 PM	8.5		
Total				8.5	0	0

I Reach us at DKSH Service-Hotline : +66 2 639 7000
2533 Sukhumvit Road, Bangkok, 10260, Phrakhanong, Bangkok, Thailand
Phone +66 2 639 7000 Fax +66 2 333 1026

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EXPENSES

Part No	Expense Type	Description	Line Quantity
---------	--------------	-------------	---------------

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Results

Item	Characteristic	Before	After	Remark
1	Visual inspect	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
2	Power supply (220-240 VAC)	☒ 220 VAC	☒ 220 VAC	
3	Computer	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
4	Program	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
5	Auto sampler	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
6	Module holder			
	- Motor pump	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Pump tube	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Air-injection	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Chemistry manifolds, Switching valve, Coil, Membrane	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
7	Detector			
	- Filter	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Flow cell	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Lamp	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
8	Interface	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
9	Rinsing valves	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
10	Temperature / Reactor	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
11	Flame photometer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
12	UPS / Stabilizer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	

Warning and Error Checked

Item	Event	Before	After
13	Error list	☒ None ☐ Approv	☒ None ☐ Approv

DKSH Technology Limited (Head office)
2533 Sukhumvit Road, Bangkok, Phrakhanong, Bangkok, 10260
Phone +66 2 639 7000, Mobile +66 91 813 8681, yongyuth.yc@dksh.com, www.dksh.com

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Page 1/2

Check with Standard

Item	Characteristic	Before	After	Remark
14	Base Line Test	☐ Pass ☐ Fail ☐ N/A	☐ Pass ☐ Fail ☐ N/A	
15	Detector Signal Test	☐ Pass ☐ Fail ☐ N/A	☐ Pass ☐ Fail ☐ N/A	

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Page 6/2

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
 534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL 0-2717-3000-29 FAX 0-2719-9484

Certificate of Calibration

Cert.No.: 25CH408
 Page: 1 of 3

Equipment :	pH Meter
Manufacturer :	EcoSense
Model :	pH100A
Serial No. :	JC04744
ID No. :	UAE.EFM.058/2566 (EFM.pH.01/66)
Condition As-Received:	Used Item
Received Date :	01 April 2025
Calibration Date :	03 April 2025
Reference :	2504-0031WSC-1
Submitted by :	United Analyst and Engineering Consultant Co., Ltd. 3 Soi Udumuk 41, Sukhumvit Road, Bangkok, Phrakhanong, Bangkok 10260
Ambient Temperature :	(25 ± 2.5) °C
Relative Humidity :	(50 ± 15) %
Calibration Procedure :	In - house method : - CP-CH5 by direct measurement with DC voltage standard and direct measurement with certified reference material (CRM) - CP-CH8 by comparison with temperature standard
Calibrated by :	Walalak Sirithan
Approved by :	 Approved Signatory
() Chakrit Waewwanjua () Ponpan Paipim (✓) Sathip Meangmai	
Issue Date :	4 April 2025

The Uncertainties are for a confidence probability of approximately 95%

This certificate may not be reproduced other than in full, except with the prior written
 Approval of the head of Corporate Services 3: Equipment Calibration and Testing Services.

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Cert.No.: 25CH408
Page: 2 of 3

Condition of this calibration result

1. Reference Standard Instrument

Instrument	Serial No.	ID No.	Cert. No.	Due Date
1) Document Process Calibrator	54030049	130RC116	24E2759	25 Aug 2025
2) Ref. Standard Thermometer	4982054	110RC044	24I757	14 July 2025

- This Certification is traceable to SI Through Technology Promotion Association (Thailand - Japan)

2. Certified Reference Materials :The measurement results are traceable to SI through Hach Lenge GmbH Ltd.,
Deutsche Akkreditierungsstelle, Accredited No.D-RM-15184-01-00
:The measurement results are traceable to SI through CPA chem Ltd.,
ANSI-ASQ National Accreditation Board, Accredited No. AR-1835

Buffer Solution	Manufacturer	Lot No.	Exp. date
pH 4.007	CPA chem	1066665	18 Jan 2027
pH 6.999	Hach Lenge GmbH	C03220	29 Oct 2026
pH 10.010	CPA chem	1066669	18 Jan 2026

3. This certificate is valid only to the item calibrated on date and place of calibration.

Calibration Results

Function : mV Measurement

Performing standard curve by Document Process Calibrator at pH (4,7)(7,10)

Unit Under Calibration	Nominal Value	Standard Voltage Input	Actual Reading		Uncertainty of Measurement ($\pm$ mV)	Coverage factor k
	pH	mV	mV	pH		
pH Meter S/N: JC04744	4.00	177.48	177	4.01	0.58	2.00
	7.00	0.00	0	7.00	0.58	2.00
	7.00	0.00	0	7.00	0.58	2.00
	10.00	-177.48	-177	10.01	0.58	2.00

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
534/4 PATTANAKARN ROAD SOI 18, SUKHUANG, SUKHUANG BANGKOK 10250
TEL 0-2717-3000-29 FAX 0-2719-9484

Certificate of Calibration

Cert. No.: 25TM205
Page : 1 of 3

Equipment : BOD Incubator

Manufacturer : Arco

Model : UC4-1320

Serial No. : 13URC4S013201

ID No. : UAE.WAO.015/2561

Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
3 Soi Udomsuk 41, Sukhumvit Road
Bangchak, Phraekhanong
Bangkok 10260

Location : Lab. Floor 2

Received Order : 08 February 2025

Calibration Date : 08 February 2025

Ambient Temperature : (26 ± 10) $^{\circ}$ C

Relative Humidity : (50 ± 30) %

AC Line Voltage : (220 ± 22) V

Calibrated by : Krisda Malee

Approved by :
() Chakrit Waewwanjua
() Suwit Imjai
(✓) Kunchit Promprat

Issue Date : 21 February 2025

The Uncertainties are for a confidence probability of approximately 95%

This certificate may not be reproduced other than in full, except with the prior written
Approval of the head of Corporate Services 3: Equipment Calibration and Testing Services.

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Equipment : BOD Incubator
Condition As-Received : Used Item
Reference : 2502-01660C-1
Procedure Used :-

Calibration were conducted using calibration procedure CP-OT02 based on TLAS G-20 according to direct measurement method with Data Acquisition which connected with Resistance Temperature Detector (RTD).
The temperature scale used was based on ITS-90.

Condition of this result of calibration

1. Reference standard Instrument:-

Instrument	Serial No.	Cert. No.	Traceable	Due Date
1) Data Acquisition	MY57013823	24LM71	TPA	12 May 2025

2. This certificate is valid only to the item calibrated on date and place of calibration.

3. This certification is traceable to the International System of Unit.


Remark : TPA : Technology Promotion Association (Thailand - Japan)

Result of Calibration : (*) Without Adjustment

Function of UUC* : Temperature Source

Fresh air setting : Not Available

Environment during calibration		
	Beginning	Finished
Temp. ($^{\circ}$ C)	26	25
REL.Humid. (%)	49	52
AC Supply (Volt)	221	220

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Equipment : BOD Incubator
 Condition As-Received : Used Item
 Reference : 2502-01660C-1
 Result of Calibration : (*) Without Adjustment
 Function of UUC* : Temperature Source
 Fresh air setting : Not Available

Cert.No.: 25TM205
 Page : 3 of 3

Calibration Point (°C)	UUC* Setting (°C)	UUC* Reading (°C)	Temperature stability (± °C)	Temperature uniformity (°C)	Overall Variation (°C)	Coverage Factor k
20.0	20.0	19.9	0.36	0.56	0.99	2

Calibration Point (°C)	Measured Temperature (°C)									Uncertainty (±°C)
	Position									
	1	2	3	4	5	6	7	8	9 (ref.)	
20.0	19.841	19.714	20.110	19.862	19.747	19.710	19.676	19.789	19.695	0.54

Average* : The average of 30 values in each position.
 Temperature stability : One-half of the greatest maximum difference of measured temperature at any one sensor.
 Temperature uniformity : The maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time or at as close an observation time as possible to determine the temperature pattern or homogeneity within the chamber under steady-state conditions.
 Overall Variation : The Difference of the maximum and minimum measured temperatures throughout observation.
 UUC* : Unit Under Calibration

Note : The reported uncertainty of measurement was included stability and excluded uniformity .

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
 53/4 PATANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL.0-2717-3000-29 FAX.0-2719-9484

ILAC-MRA

Certificate of Calibration

Cert.No.: 25CH862
 Page.: 1 of 3

Equipment : pH Meter
 Manufacturer : Horiba
 Model : LAQUA-PH210
 Serial No. : HA1F0002
 ID No. : UAE.EFM.200/2564(EFM.pH.08/64)
 Condition As-Received : Used Item
 Received Date : 22 July 2025
 Calibration Date : 23 July 2025
 Reference : 2507-0753WSC-2
 Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
 3 Sol Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260

Ambient Temperature : (25 ± 2.5) °C
 Relative Humidity : (50 ± 15) %
 Calibration Procedure : In - house method :
 - CP-CH5 by direct measurement with DC voltage standard and direct measurement with certified reference material (CRM)
 - CP-CH8 by comparison with temperature standard

Calibrated by : Walalak Sirithean

Approved by :

Saithip
 Approved Signatory

() Chakrit Waowwanjua
 () Porpan Palpim
 (✓) Saithip Moangmai

Issue Date : 24 July 2025

The Uncertainties are for a confidence probability of approximately 95 %

This certificate may not be reproduced either in full, except with the prior written Approval of the head of Corporate Services 3: Equipment Calibration and Testing Services.

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Cert.No.: 25CH862
 Page.: 2 of 3

Condition of this calibration result

1. Reference Standard Instrument

Instrument	Serial No.	ID No.	Cert.No.	Due Date
1) Document Process Calibrator	54030049	130RC116	24E2759	25 Aug 2025
2) Ref. Standard Thermometer	4982054	110RC044	251708	03 July 2026

- This measurement result is traceable to SI through Technology Promotion Association (Thailand - Japan)

2. Certified Reference Materials

:The measurement results are traceable to SI through Hach Lange GmbH Ltd.
 Deutsche Akkreditierungsstelle, Accredited No.D-RM-15184-01-00
 :The measurement results are traceable to SI through CPA chem Ltd.,
 ANSI-ASQ National Accreditation Board, Accredited No. AR-1835

Buffer Solution	Manufacturer	Lot No.	Exp. date
pH 4.007	CPA chem	1088665	18 Jan 2027
pH 7.000	Hach Lange GmbH	C03232	02 Dec 2026
pH 10.010	CPA chem	1114385	08 June 2028

3. This certificate is valid only to the item calibrated on date and place of calibration.

Calibration Results

Function : mV Measurement

Performing standard curve by Document Process Calibrator at pH (4.7)(7.10)

Unit Under Calibration	Nominal Value	Standard Voltage Input	Actual Reading		Uncertainty of Measurement (±mV)	Coverage factor k
	pH	mV	mV	pH		
pH Meter S/N: HA1F0002	4.00	177.48	177.5	4.01	0.056	2.00
	7.00	0.00	0.0	7.02	0.056	2.00
	7.00	0.00	0.0	7.02	0.056	2.00
	10.00	-177.48	-177.5	10.01	0.058	2.00

Cert.No.: 25CH862
 Page.: 3 of 3

Calibration Results

Function : pH Measurement

Performing three buffers standard curve by using buffer nominal pH (4.7)(7.10)

Unit Under Calibration	Standard pH Buffer Solution	Actual pH Reading	Actual mV Reading (mV)	Uncertainty of pH Measurement (±)	Coverage factor k
pH Electrode S/N: Q9AG0067	4.007	4.01	153.1	0.0085	2.05
	7.000	6.99	-18.1	0.0092	2.00
	7.000	7.01	-18.3	0.0085	2.00
	10.010	10.01	-188.7	0.0062	2.00

Function : Temperature Measurement

(*) Without adjustment

This equipment was connected with Temperature Probe:

- Model : 9652-10D
 - Serial No. : Q9AG0067
 Dimension of probe
 - Length : 103 mm.
 - Diameter : 15 mm.
 - Immersion Depth : 90 mm.

Calibration Point (°C)	Standard Temperature (°C)	UUC* Reading (°C)	Error (°C)	Uncertainty of measurement (± °C)	Coverage factor k
15.0	15.002	15.0	-0.002	0.13	2.00
30.0	30.002	30.0	-0.002	0.13	2.00
45.0	45.003	44.9	-0.103	0.13	2.00

Remark : - UUC* = Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

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Calibration Certificate

Certificate No.: 2502226-002-01
Client name: UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Address: 3 Soi Udomsuk 41, Sukhumvit Road,
Bangchack, Prakhonong, Bangkok 10260

Page 3 of 4

Equipment: Electronic Balance
Manufacturer: METTLER TOLEDO
Model: XSR205DU
Serial No.: C210685394
ID No.: UAE.WAO.010/2565
Order No.: 2502226
Operation No.: 2502226-002
Date of Receipt: 19 March 2025
Date of Calibration: 20 March 2025

Calibrated by Mr.Yothin Charoensuk
Scientist

Approved by *for N. Ningsubatt*
(Mr.Pheraphat Tuanjil)
Manager, Division of Calibration Laboratory
Responsible for the Technical Management Team

Date of Issue: 25 March 2025

The uncertainties are for a confidence probability of approximately 95%

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation Scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the units of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full, except with the prior written approval of the National Food Institute.

FCS-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-002-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C210685394
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.010/2565

Page 2 of 4

Date of Calibration: 20 March 2025
Environment Condition: Ambient Temperature: 21.2 ± 0.6 °C Relative Humidity: 48 ± 3.5 %
Place of Calibration: 208 Balance Room, UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Condition of Equipment: Good Const on
Condition of This Results of Calibration:

1. Calibration Method: NFI Method W-MA-001 In-House Method based on UKAS Lab 14: 2019

2. Reference Standards:

Reference Standard	Model	Serial No.	Calibrated By	Certificate No.	Due Date
Standard Weight Class E2	1mg to 200g	8505567572	TCS	M24041005	19 April 2025
Instrument	Model	Serial No.	Calibrated By	Certificate No.	Due Date
Thermo-Hygro Meter	608-H1	NFI.BTH D1723	Quality Reborn	QR25-0547	10 February 2026

3. This certification is traceable to SI UNIT

4. This certificate was certified only for the instrument we calibrated.

5. This result of calibration was found accurate as shown on date and place of calibration only

Calibration Results:

1. Repeatability of Reading:

Nominal Value (g)	Standard Deviation of Reading (g)
40	0.0000012
80	0.0000012
100	0.0000000
200	0.0000000

2. Off-Center Error:

A mass of 100 g was placed and moved to various position on pan
The balance reading obtained is given in the table.

1	2	3	4	5	6	(Maximum Difference)
(g)	(g)	(g)	(g)	(g)	(g)	(g)
100.0001	100.0001	100.0001	100.0001	100.0001	100.0001	0.0000

FCS-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-002-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C210685394
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.010/2565

Page 3 of 4

Date of Calibration: 20 March 2025

Calibration Results: (Continued)

Calibration Range: 0-80 g

Calibration Adjustment: Internal Calibration

3. Departure from Nominal Value: (Range: 0 - 82 g; Resolution: 0.00001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (± g)	Coverage Factor k
Unload	0.000000	0.000000	0.000000	0.00000037	2.00
0.001	0.001000	0.001000	0.000000	0.00000030	2.00
0.005	0.005000	0.005001	-0.000001	0.00000032	2.00
0.01	0.010000	0.010002	-0.000002	0.00000039	2.00
0.05	0.050000	0.050001	-0.000001	0.00000056	2.00
0.1	0.100011	0.100002	0.000009	0.0000011	2.00
0.5	0.500016	0.500004	-0.000002	0.0000014	2.00
1	1.000033	1.000005	-0.000005	0.0000016	2.00
2	2.000023	2.000006	-0.000004	0.0000017	2.00
5	5.000015	5.000006	-0.000005	0.0000020	2.00
10	10.000039	10.000005	-0.000004	0.0000026	2.00
20	20.000030	20.000007	-0.000004	0.0000037	2.00
30	30.000029	30.000009	-0.000005	0.0000030	2.00
50	50.000028	50.000008	-0.000005	0.0000038	2.00
80	80.000027	80.000013	-0.000005	0.0000041	2.00

FCS-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-002-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C210685394
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.010/2565

Page 4 of 4

Date of Calibration: 20 March 2025

Calibration Results: (Continued)

Calibration Range: >80-200 g

Calibration Adjustment: Internal Calibration

3. Departure from Nominal Value: (Range: >80 - 200 g; Resolution: 0.00001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (± g)	Coverage Factor k
90	90.00010	90.00005	-0.00005	0.000015	2.00
100	100.00006	100.00001	-0.00005	0.000016	2.00
110	110.00007	110.00002	-0.00005	0.000017	2.00
120	120.00009	120.00002	-0.00007	0.000018	2.00
130	130.00010	130.00003	-0.00007	0.000019	2.00
140	140.00013	140.00002	-0.00011	0.000019	2.00
150	150.00009	150.00002	-0.00007	0.000021	2.00
160	160.00010	160.00002	-0.00008	0.000022	2.00
170	170.00012	170.00002	-0.00010	0.000023	2.00
200	200.00013	200.00002	-0.00011	0.000028	2.00

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

FCS-012 Revision: 01 Date: 20-04-65

No.	Instrument/Equipment	Parameter	Manufacturer	Model/Serial No.	Calibrator	Certification No.	Date of Calibration	Due date of Calibration*
1	Atomic Absorption Spectrometer	CHROMIUM CHROMIUM HEXAVALENT COPPER LEAD MANGANESE NICKEL ZINC	Agilent Technologies	AA240PS / M113160001	Agilent Technologies (Thailand) Co., Ltd.	Preventive Maintenance Checklist	30/1/2025	29/1/2026
2	Atomic Absorption Spectrometer	ARSENIC	Perkin Elmer	ProAAle 900P / P820001902	Perkin Elmer Co., Ltd.	Preventive Maintenance Report	24/02/2025	20/02/2026
3	Analysis Balance	TOTAL DISOLVED SOLIDS	Mettler Toledo	XS205DU / C2150683394	National Food Institute, Ministry of Industry, Thailand	2502226-002-01	20/09/2025	19/09/2026
4	Continuous Flow Analyzer(CFA)	CYANIDE	Solar Analytical B.V., the Netherlands	Sara-3300-02 / 180688	DOH (Thailand) Ltd.	WO-00076079	23/05/2025	22/05/2026
5	Digestion Units	TOTAL ALUMINUM, NITROGEN	Foss Tecator	2520 Auto / 91794669	National Food Institute, Ministry of Industry, Thailand	2501410-001-01	27/11/2025	26/11/2026
6	Mercury Analyzer	MERCURY	NVC Japan	R4-500 / 17180278	Coast Group Corporation Ltd.	Preventive Maintenance Report	27/07/2025	1/7/2026
7	pH Meter	pH	HANNA	LAQUA-34-E10 / A11-0002	Technology Promotion Association (Thailand-Japan)	25C-062	23/07/2025	22/07/2026

Due Date of Calibration*: Based on the annual calibration plan. At least 1 time per year.

United Analyst and Engineering Consultant Co., Ltd. (UAE)
Certified Laboratory ISO/IEC 17025

Certificate Page 1 of 1

Agilent 55 240 280 Series Atomic Absorption Spectroscopy Systems

Preventive Maintenance Checklist

Agilent Preventive Maintenance provides factory recommended service for your analytical systems to assure reliable operation and the accuracy of your results.

Delivered by highly trained and certified service engineers using genuine Agilent parts and supplies, Agilent Preventive Maintenance provides everything you need to reduce unplanned downtime and keep your systems operating at their peak. This checklist will be completed at the end of the service and provided to you as a record of the installation.

Note: While non-current production AA instrument and/or accessory models are not covered specifically in this document it can be used as a basic reference.For more information about Agilent Technologies services please visit our web site using the following URL: <http://www.agilent.com/en-us/services>

Introduction

Customer Information

- Customers should provide all necessary operating supplies upon request of the engineer.
- A customer representative should be available to the engineer while performing the preventive maintenance procedures.
- Any parts, not included in the Parts Lists section of this document, are not part of the recommended Preventive Maintenance service, nor are they included in the price of this service.
- If a system requires the use of extra or special procedures and/or parts for the maintenance service, then these must be ordered separately and charged as a repair, which may incur additional costs.

Revision: 10-00, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Important Customer Web Links

- For more information about Agilent Technologies services, please visit our website using the following URL: <http://www.agilent.com/en-us/products/crosslab-instrument-services/service-repair>
- To access Agilent University, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- A useful Agilent Resource Center web page is available, which includes short videos on maintenance, quick lists of consumables for new instruments, and other valuable information. Check out the Resource Page here: <https://www.agilent.com/en-us/agilentresources>
- Need technical support, FAQs, supplies? – visit our Support Home page at <http://www.agilent.com/search/support>
- Get answers. Share insights. Build connections. Join the Agilent Community at <https://community.agilent.com/welcome>

Service Engineer's Responsibilities

- Contact the customer and ensure that all necessary supplies are available before the preventive maintenance visit.
- Confirm the ability of the instrument to deliver continued safe operation as established via the Agilent AA safe operation flow chart. (Refer directly to the AA 55/240/280 Preventive Maintenance Scope of Work to make this decision.)
- Only select those pages that relate to the system or module being serviced.
- Complete empty fields with the relevant information.
- Complete the relevant checkboxes in the checklist using either a "X" or tick mark "✓".
- Check "Section not applicable" check boxes to indicate services/tasks not delivered, as appropriate.
- Complete the Preventive Maintenance service in the order of the tasks listed.
- Complete the Service Review section together with the customer.
- Complete the fields for page numbers at the foot of each selected page.
- Complete the total number of pages field in the Service Completion section.
- Ask the customer to sign the Service Completion section including the customer's and your signature.

This information is subject to change without notice.

Revision: 10-00, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Instrument Maintenance

System Information

- ☐ Check this box if an instrument configuration report is attached instead of completing the table.

Instrument System Name and ID	240 FS AAS
Instrument System Site and Location	United Analyst and Engineering Consultant

List System Component Product Numbers	List the Serial Numbers of each Component
1 G 9432 A	M113160001
2	
3	
4	
5	
6	
7	
8	
9	

Preparation, Safe operation and Initial performance checks

Revision: 10-00, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Agilent AA safe operation flow chart inspections (to determine if the PM can be performed).

NOTE: If by following the flow chart the instrument is deemed to be unsafe for continued use you **MUST NOT** continue PM work. Inform the customer immediately of the Agilent recommendation that use of the instrument be discontinued.

- ☒ Discuss any specific issues with the customer before starting.
- ☐ For HF application systems, if standard sample introduction system was not installed, ask the customer to install it. 8119
- ☒ Review the instrument logbook for recorded problems and comments.
- ☒ Save instrument control settings before starting the procedure.
- ☒ Perform a general inspection of the system for cleanliness.
- ☒ Check for proper installation of parts, assemblies, sensors etc.
- ☒ Check system for required installation of components, settings as defined by current Service Notes
- ☒ Check for required firmware updates and verify with customers if they would like them installed.
- ☒ Use SVD to perform a Full Wavelength Scan for Cu HCL - "As found test.1"
- ☒ Perform a Basic Cu ABS test - "As found test.2"
- ☒ Print the Details page or screen captures of the test results and attach to the end of this checklist.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Preventive Maintenance Procedures

FLAME SYSTEM section

- ☐ Section not applicable

Electronic components

- ☒ Review and confirm instrument configuration data in SVD
- ☒ Confirm power supply voltages using the *SVD Power Supply diagnostic*.
- ☒ For Dual Beam instruments - Confirm RBC frequency using the *SVD RBC frequency diagnostic*.

Mechanical components

- ☒ Check the burner adjuster controls for complete and free movement. If the burner adjuster needs lubrication, use Molykote 321 or mineral based molybdenum disulphide grease
- ☒ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☒ Monochromator drive
 - ☒ Slit drive
 - ☒ Lamp selector
 - ☐ ABA

Optics components

- ☒ Check that external optical surfaces are clean – Clean or replace as required
- ☒ Use SVD and perform *Mono Wavelength Correction*.
- ☒ Use SVD and perform *Slit Calibration*.
- ☒ Use SVD and perform *Grating Squareness Diagnostic*.
- ☒ Use SVD and perform *Zero Order Offset/Mono Correction*.
- ☒ Use SVD and perform *Wavelength Repeatability*.
- ☒ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.
- ☒ Check that the signal energy of the D2 and HC lamps track properly. Advise customer if their D2 lamp is showing emission degradation due to age.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Sample Introduction and Atomization

- ☒ Inspect the burner interlock plate to ensure that the interlock pin is secure and correct for the burner type.
- ☒ Clean the burner slot with a clean white cord.
- ☒ Check the uniformity of the slot width.
- ☒ Clean the burner if required.
- ☒ Change the burner o-ring.
- ☒ Clean the nebulizer, spray chamber and liquid trap.
- ☒ Change all o-rings and seals in the nebulizer, nebulizer block and spray chamber.
- ☒ Check that the pressure relief bung releases readily.
- ☒ Change o-rings on the fuel and oxidant delivery bars.
- ☒ Leave the liquid trap EMPTY and verify the flame will not ignite in this state.
- ☒ Refill liquid trap and check that overflow drains freely into the drain/waste tube.
- ☒ Check the drain/waste tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel
- ☒ Check and clean the igniter electrode

Gas handling components and safety interlocks

- ☒ Pressure test for leaks
- ☒ Leak test gasbox internal components and connections
- ☒ Check safety interlock status and operation using the *SVD interlock monitoring diagnostic*.

Analytical performance for Flame systems

- ☒ Ignite a flame.
- ☒ Check that you can adjust the nebulizer uptake rate from 4 to 6.5 mL per minute.
- ☒ Optimize the instrument ready to perform Cu sensitivity test.
- ☒ Create a manual method to perform a Basic Cu ABS test - "Final Performance Testing"
- ☒ Run a PM completed sensitivity test for a 5 ppm copper sample and record the results in the AA PM Performance test results and measurements table.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

FURNACE SYSTEM section

- ☒ Section not applicable

Electronic components

- ☐ Review and confirm instrument configuration data in SVD
- ☐ Confirm power supply voltages using the *SVD Power Supply diagnostic*.

Mechanical components

- ☐ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☐ Monochromator drive
 - ☐ Slit drive
 - ☐ Lamp selector

Optics components

- ☐ Check that external optical surfaces are clean – Clean or replace as required
- ☐ Use SVD and perform *Mono Wavelength Correction*.
- ☐ Use SVD and perform *Slit Calibration*.
- ☐ Use SVD and perform *Grating Squareness Diagnostic*.
- ☐ Use SVD and perform *Zero Order Offset/Mono Correction*.
- ☐ Use SVD and perform *Wavelength Repeatability*.
- ☐ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.

Gas handling, water system and workhead component checks

- ☐ Inspect the GTA workhead gas hoses and connections for leaks.
- ☐ Pressure test for gas leaks
- ☐ If the cooler system is accessible (stand-alone) check for correct operation and coolant/water level – this includes any temperature and pressure settings plus filter cleaning (air flow and water)
- ☐ Inspect the GTA workhead water hoses and connections for leaks.
- ☐ Check all graphite components and replace if necessary.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Tube
- ☐ Electrodes
- ☐ Shroud

- ☐ Check and clean the end windows on the workhead.
- ☐ Check safety interlock operation.

Analytical performance for Furnace systems

- ☐ Optimize the instrument ready to perform Cu sensitivity test.
- ☐ Run the sensitivity test for a 25 ppb copper sample and record the results in the results table.

PSD autosampler accessory for Furnace systems

- ☒ Section NOT Applicable
- ☐ Check condition of the PSD capillary – replace if necessary
- ☐ Check condition and operation of PSD syringe – ensure it does not have air locks and bubbles.
- ☐ Change PSD rinse bottle o-ring.
- ☐ Check and clean the rinse vessel
- ☐ Check the drain tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel.
- ☐ Ensure that the waste vessel is suitable for use with the furnace system.

Sample introduction pump system (SIPS) accessory

- ☒ Section NOT Applicable
- ☐ Re-torque screws securing the hubs, presser arms and pump rotors.
- ☐ Adjust each roller so that it rotates freely
- ☐ Wipe clean the pump rotor rollers and pump bands with a dry clean cloth
- ☐ Ensure that the presser arms and the surfaces near the pump are free from dirt and spills
- ☐ Remove the pump module rear cover and check for the incursion of liquids and any signs of corrosion.
- ☐ Re-torque the nuts that fasten the motor mounting plates to the chassis.
- ☐ Check clips securing the diluents holder and replace if necessary.
- ☐ Disconnect, clean T-piece, and reassemble the tubing using the following steps.

Revision 10/09, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Remove the T-piece by disconnecting the pump tubes, the pump bands and all other tubing.
- ☐ Place the T-piece in an ultrasonic bath containing strong detergent 1-5% Decon 30 or similar, for approximately 5-10 minutes
- ☐ Wash the T-piece under a tap with a strong flow of water
- ☐ Rinse with distilled water through all of the inlets in the reverse direction to normal sample flow.
- ☐ Reassemble.

Sample preparation system (SPS 4) accessory

- ☒ Section NOT Applicable

The Agilent SPS 4 autosampler is designed to need minimal maintenance.

The following maintenance requirements are suggested to maintain the performance of the autosampler.

- ☐ Cleaning the spill tray, rack location mat, end frames and chassis accessories with a damp soft cloth and diluted mild detergent.
- ☐ Cleaning the autosampler cover panels with domestic window cleaner.
- ☐ Checking the X- axis and Z- axis drive belts for cracks, splits, damaged teeth, excessive fraying, color changes or degradation from fumes.
- ☐ Check the X axis, Theta axis and Z axis FFC cables for cracks, incorrect positioning, damaged edge or damaged connectors.

NOTE: The autosampler requires no extra lubrication throughout its lifetime.
For further details refer to the SPS 4 service manual G8410-90050.

Sample preparation system (SPS 3) accessory

- ☒ Section NOT Applicable

- ☐ Check the x-axis and z-axis timing belts – Replace if there is any cracks, splits or color deterioration and belt tension.
- ☐ Check belt tensions - adjust if required
- ☐ Check the lubrication pad for single x-axis shaft. If pad is dry or customer has observed any vibration or erratic movements of the x-axis carriage, add 1 mL of Dow Corning 200 0 Fluid, 200 CS into the well.
- ☐ Check the auto-sampler ability to find tube positions - Calibrate if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Revision 10/09, Issued November 2021

© Agilent Technologies, Inc. 2021

Vapor generation accessory VGA (hydride generator)

- ☒ Section NOT Applicable
- ☐ Inspect VGA gas supply hose.
- ☐ Inspect/replace VGA pump tubing
- ☐ Check low gas pressure interlock setting – adjust if required.
- ☐ Check precision orifice gas flow setting – adjust if required.
- ☐ Check gas regulator pressure to 46 psi (325 kPa) – adjust if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

UltraAA lamp accessory (external)

- ☒ Section NOT Applicable
- ☐ Check the condition of the power cable.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Restore System

- ☐ If you have altered the customer's instrumentation during the course of PM, restore to the original status to allow the customer to conduct their normal activities (e.g., reload the customer's method.)

Guidance

If the PM service is performed prior to a qualification service, then use the qualification procedure as a guide for final instrument set up and checkout.

Revision 10/09, Issued November 2021

© Agilent Technologies, Inc. 2021

Signature Page

Service Review

- ☒ Attach available reports/printouts of all tests to this documentation.
- ☒ Record the Preventive Maintenance service activity in the customer's records/logbook.
- ☒ Update/reset instrument maintenance counters as appropriate.
- ☒ Affix the PM sticker to the system or instrument logbook based on the customer's request.
- ☒ Complete the Service Engineer Comments section if there are additional comments.
- ☒ Review this service, parts replaced, and test results obtained with the customer.
- ☒ If the instrument firmware was updated, record the details of the change in the Service Engineer's Comments box or if necessary, in the customer's IQ records

Test Results

Test Description:	Expected Test Result	Actual Test Result
Flame optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	49 %
Flame performance test with 5 ppm copper sample		
Air /acetylene, mixing paddle removed	Abs value > 0.5	0.559%
Air /acetylene, mixing paddle installed, 10 replicates	%RSD < 1.0	0.2 %
Deuterium furnace optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	—
Deuterium furnace performance test with 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0%	—
Abs value	≥ 0.15	—
Zeeman furnace analytical performance: 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0%	—
Abs value	≥ 0.10	—
MSR%	≥ 70 %	—

Revision 10/09, Issued November 2021

© Agilent Technologies, Inc. 2021

AA consumable and parts list table

Part Description	Part Number	Product/Model	PM supplied or Consumable	Instrument Type
Test Solution - Cu 5ppm solution	6610050100	50 55 140 243 280	PM supplied	Common
Test Solution - Blank solution	5190-0001	50 55 140 243 280	PM supplied	Common
Copper, 100C ug/ml, 100ml	5190-8279	50 55 140 243 280	*	Common
Kit, Ask 7 O-rings, aqueous, complete set	9910053400	50 55 140 243 280	PM supplied	Flame
Organic Kit	9910053500	50 55 140 243 280	PM supplied	Flame
Wire Nebulizer Cleaning	9910024700	50 55 140 243 280	consumable	Flame
Tubing Capillary Std Nebs	9910024800	50 55 140 243 280	consumable	Flame
Capillary Tube Hvac Neb (3) (organics only)	9910044000	50 55 140 243 280	consumable	Flame
Glass impact beads (5/pk)	9910025700	50 55 140 243 280	consumable	Flame
Teflon impact beads (5/pk) (organics only)	9910053300	50 55 140 243 280	consumable	Flame
Burner cleaning strip (100/pk)	9910053900	50 55 140 243 280	consumable	Flame
Window UV silica - round (right side)	2010082500	50 55 140 243 280	PM supplied	Common
Window UV silica - rectangular (left side)	2010082500	50 55 140 243 280	PM supplied	Common
Pail adhesive window (round)	4910012700	50 55 140 243 280	PM supplied	Common
Pad adhesive window (rectangular)	4910012800	50 55 140 243 280	PM supplied	Common
Electrode kit (1 pr) (D2)	63100C3400	GTA120	PM supplied	Furnace
Shroud (D2)	63100C3100	GTA120	PM supplied	Furnace
Zeeman electrode kit (1 pr)	63100C3500	GTA120	PM supplied	Furnace
Zeeman shroud	63100C3600	GTA120	PM supplied	Furnace
O-ring PDS rinse bottle	6910025900	PDS120	PM supplied	Furnace

* For engineers who only service AA instruments 5190-8279 can be used as a cheaper alternative for 6610030100.

Items classified as PM supplied in the above table are included in the standard PM. Those classified as consumable should be provided by the customer or charged to the customer if supplied by the Agilent service engineer.

Revision: 10/00, Issued November 2021

© Agilent Technologies, Inc. 2021

Service Engineer Comments (optional)

If there are any comments or issues with the instrument or parts, please enter them in the comments section of the report.

Service Completion

Service request number: 6007549143 Date service completed: 30 Jan 2023
 Agilent signature: Kanyakorn S. Customer signature: Sanda Y
 Total number of pages in this document: 13

Revision: 10/00, Issued November 2021

© Agilent Technologies, Inc. 2021

SVD Results Report

Report ID: <u>1</u>	Diagnostic Start Time: <u>1/30/2023 9:14:28 AM</u>	Diagnostic End Time: <u>1/30/2023 9:16:00 AM</u>
Customer: UAE	Service Engineer: Kanyakorn S.	
Address: Soi Udomsuk 41, Sukhumvit Rd, Bangkok	Contact Details: 026376363#1	

Configuration:

Serial Number: MY13160001 Turret Type: Automatic
 Instrument Model: Varian AA140/240/280 Number Of Lamps: 4
 Flame Instrument: True Mono Type: Automatic
 Furnace Instrument: True Gasbox Type: Y Gas Box
 Zeeman Present: False Auto Burner Adjuster: False
 Internal Zeeman: False Mains Frequency: 50
 Internal UltraAA: False Firmware Version: 2.11
 Optics Type: Double Beam Photomultiplier Type: Normal(900nm)
 D2 CS Correction Filter: True PWB Version: 45
 Box Block Version: 1.00

EEPROM Data:

Instrument Run Hours: 69918.180 D2 Run Hours: 53396.500
 Zero Wavelength Offset: 0.133 D2 Serial Number: not set
 Mono Correction: 0.770 D2 Install Date: 1/1/1970
 Flame Hours: 32411.634 D2 Original Intensity: 1.000
 D2 Last Intensity: 475.000

Frequency:

Averaging Period: 30.0
 Datapoint Count: 20

Upper Limit: 51.00 Highest Measured Frequency: 50.00
 Average Frequency: 50.00
 Lower Limit: 49.00 Lowest Measured Frequency: 50.00

Result: **Passed**

Power Supply:

Averaging Period: 30.0
 Datapoint Count: 20



	Lower Limit (V)	Actual (V)	Upper Limit (V)	Result:
12.00 V Rail	10.80	12.12	13.20	Passed
-12.00 V Rail	-13.20	-11.90	-10.80	Passed
5.00 V Rail	4.50	5.04	5.50	Passed
310.00 V Rail	270.60	320.00	341.00	Passed

Optics

Beam Balance:

Lamp Type: Copper
Lamp Socket Used: 3

Peak Selected: 324.80
Lamp Alignment: **Performed**

Grating: S

Lamp Element(s): Copper
Lamp Socket Used: 3
Lamp Current(mA): 4.00
Slit Width(nm): 0.5
Wavelength(nm): 324.80
Lamp Alignment: **Performed**

	Lower Limit(nm)	Upper Limit(nm)	Result:
Zero Order	649.97	0.10	Passed
First Order	324.80	325.15	Passed
Second Order	0.10	649.97	Passed

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Interlocks:

Gas Flow Interlock: **Working**
Flame Detect: **Working**
HSC Gas Flow Interlock: **Working**
Flame Out Gas Flow Interlock: **Working**
Gas Flow Interlock: **Working**
Pressure Interlock: **Working**
Liquid Gas Interlock: **Working**
Flame Detect: **Working**
ECU Active: **Working**
Oxidant Pressure: **Working**
Flame Changeover: **Untested**
Ignition: **Working**

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Auto Burn-in

Lamp 1: 324.80 nm (Copper)
Lamp 2: 324.80 nm (Copper)
Lamp 3: 324.80 nm (Copper)
Lamp 4: 324.80 nm (Copper)
Lamp 5: 324.80 nm (Copper)
Lamp 6: 324.80 nm (Copper)
Lamp 7: 324.80 nm (Copper)
Lamp 8: 324.80 nm (Copper)
Lamp 9: 324.80 nm (Copper)
Lamp 10: 324.80 nm (Copper)

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Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Sensitivity Test101
Comment
Methods Cu
Computer name DESKTOP-9SUFPS
Serial Number: MY13160001


Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL 2770	0.000	38.8	0.0002
Readings			
	0.0002	0.0003	0.0001

1/30/2025 10:51:40 AM

STANDARD 1	Conc. mg/L	%RSD	Mean Abs
Readings			
	0.000	0.1	0.5571

1/30/2025 10:52:22 AM

Curve Fit = Linear Origin
Characteristic Conc = 0.036 mg/L
r = 1.0000
Calculated Conc = 0.002 5100
Residuals = -0.002 0.000

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Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Precision Test
Comment
Methods Cu
Computer name DESKTOP-9SUFPS
Serial Number: MY13160001


Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL 2770	0.000	64.1	-0.0002
Readings			
	-0.0003	-0.0003	-0.0001

1/30/2025 10:46:52 AM

STANDARD 1	Conc. mg/L	%RSD	Mean Abs
Readings			
	0.000	0.3	0.6062

1/30/2025 10:47:24 AM

Curve Fit = Linear Origin
Characteristic Conc = 0.105 mg/L
r = 1.0000
Calculated Conc = -0.002 5000
Residuals = 0.002 0.000

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Pump speed Medium

Key to tube colors

- Sample
- Calibration
- Calibration/QC
- Sample/QC
- Not Assigned

Align Probe

Reset

Stop

Peak

Optimization Lamp

HC Lamp

1.30

1.00

0.50

0.00

0.917

Optimize

Optimize

Rescale

Set Zero

Gain 49 %

Sensitivity Check 15 mg/L gives about 0.2 Abs at 324.8 nm, A/A burner

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United Analyst and Engineering Consultant Co., Ltd.

Automatic Mercury Analyzer

Model : RA-4500

Preventive Maintenance Report

SERIAL No. RA-4500 : 11780278

Soft version : Ver 2.0.5

ROM version : Ver 2.0.2

DATE : 02 JULY 2025

DUE DATE : 02 JULY 2026

INSPECTED BY : *Natthaphong P.*
(Natthaphong P.)APPROVED BY : *Kitichai S.*
(Kitichai S.)

Kinetic Solutions Company Limited.
2, Soi Lat Krabang 1, Lat Krabang Subdistrict,
Lat Krabang District, Bangkok 10520
Tel. (+66) 062-769-5221

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Inspection

ITEMS	SPECIFICATIONS		RESULT	JUDGE
1. Quantity	-	Accessories are completed.	GOOD	OK
2. Appearance	2.1 Overall Appearance	No visible damage.	GOOD	OK
	2.2 Parts / Cables	Correctly placed and connected.	GOOD	
3. Indication	Nameplate / Label	Plate and Label are indicated.	GOOD	OK
4. Self check				
4.1 GLP Counter	Mercury Lamp	5000 hours	110 hr	OK
	Membrane Filter	2000 hours or 1 year after replace	< 1 hr	OK
	Main Pump Tube	750 hours or 1 year after replace	< 1 hr	OK
	Absorbed Hg	1500 mg	< 1 hr	OK
	P1 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P2 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P3 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P4 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P5 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P6 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	P7 Tube	2000 hours or 1 year after replace	< 1 hr	OK
	Heater	2000 hours	42 hr	OK
	4.2 Check/Test	Flow rate Adjustment	Flow rate 0.14 - 0.20 L/min	0.18 L/min
Signals Detector		V.SIG is 3.5 - 4.5 V. V.REF is 3.5 - 4.5 V.	4.03 V. 4.03 V.	OK OK
Cooling Fan		Check the operation of cooling fan	PASS	OK
Color Sensor		signals (R,G,B) at least one nonzero	PASS	OK
Radiation Thermometer		a positive valve form thermometer	PASS	OK
Heater		heater temp rises 4 °C within 5 min.	PASS	OK
5. Heater		Temperature	At 85°C ± 2°C with 30 min.	94.7°C
6. Calibration Curve	no pretreatment	0-10ng : Max.Dev. 5.0%	0.90%	OK
7. Repeatability	100 µg/L (n=5)	Average: 100 ± 5 µg/L	101.4 µg/L	OK
		C.V. ≤ 5.0%	0.71%	OK
8. Blank	no pretreatment	Less than 0.5 (AREA)	0.103 AREA	OK

Apparatus

NAME	Date Certified	Expiration
Mercury ICP Standard (1000 µg/mL) AccuStandard, Inc. Lot 223035027	March 10, 2023	March 10, 2028

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Remark

1. Cleaning Mercury Analyzer RA-4500

- Body case, Cell Detector, Table of Sample

2. Replace consumable one year set for RA-4500

3. Checking Mercury Analyzer RA-4500

- Mercury Lamp	- Outlet Activated Carbon Filter
- Membrane Filter	- Flow Sensor
- Main Pump Tube	- Valves
- Motor Pump	- Heater
- Motor Reagent pump	- Cooling Fan
- Motor of Arm	- Stirrer motor
- Motor of Table	- Color sensor
- Motor of Lift	- Instrument Performance

4. Adjust Optical system

5. Inspection Calibration curve and Repeatability of Mercury Analyzer RA-4500

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

Method

Method (Pretreatment: without)
111112504
100 µg 5.0 mL
Measurement Time (sec) : 120sec


STD	No.	STD (µg/L)	SVOL (mL)	CVOL (mL)	DVOL (mL)	STD (ng)	AREA (ON)	MEAS (ng)	Dev (%)	Note
0.000	1	0.000	5.000	5.000	5.000	0.0000	0.1107	0.0000	-	
100.000	2	100.000	0.025	5.000	5.000	2.500	3.2690	2.4250	3.0	
100.000	3	100.000	0.050	5.000	5.000	5.000	6.5434	4.9392	1.2	
100.000	4	100.000	0.075	5.000	5.000	7.500	9.8649	7.4895	0.1	
100.000	5	100.000	0.100	5.000	5.000	10.000	13.2068	10.0570	0.6	

SMP	No.	NAME	SVOL (mL)	CVOL (mL)	DVOL (mL)	AREA (ON)	MEAS (ng)	CONC (µg/L)	Note
1	STD 100 ppb		0.050	5.000	5.000	6.7938	5.0854	101.708	
2	STD 100 ppb		0.050	5.000	5.000	6.7129	5.0633	101.385	
3	STD 100 ppb		0.050	5.000	5.000	6.6370	5.0110	100.220	
4	STD 100 ppb		0.050	5.000	5.000	6.7489	5.0970	101.540	
5	STD 100 ppb		0.050	5.000	5.000	6.7514	5.0989	101.578	
6	Blank		5.000	5.000	5.000	0.1035	-0.0055	0.000	

Statistics	No.	NAME	TRY	AV (µg/L)	SD (µg/L)	Cv (%)
1	STD 100 ppb		5	101.4464	0.724979	0.71

Self Check

Heat check: PASS!! (25.7degC[05:00] -> 29.7degC[02:29])
Sensor check: PASS!! (25.0 - 15= 238)
Leak check: PASS!! (0.18L/min)
Sig/Ref check: PASS!! (Sig: 4.03V, Ref: 4.03V)
Drift check: PASS!! (0.0000037 - -0.0000026 = 0.0000063)

125 Main Street
New Haven, CT 06511
USA

AccuStandard, Inc.
CERTIFICATE OF ANALYSIS

AccuTrace™ Reference Standard

Catalog No: ICP-34N-1
Description: Mercury ICP Standard
Element: Mercury (Hg)
SRM: 3133
Lot: 222035027
Matrix: 10% Nitric acid
Hazards: Refer to SDS for complete safety information

Date Certified: Mar 10, 2023
Expiration: Mar 10, 2028
Density: 1.052 g/mL
Sample Size: 100 mL
Components: 1
Storage Condition: Ambient (>5 °C)

Certified Reference Material

Signal Word: Danger

Certified Concentration: 1000 µg/mL

Trace Elements in µg/mL

Ag nd<0.02	Ce nd<0.2	Od nd<0.02	Lu nd<0.02	Pb nd<0.2	Se nd<0.02	Ti nd<0.02
Al nd<0.02	Co nd<0.02	Ge nd<0.2	Mg nd<0.02	Pd nd<0.2	Be nd<0.2	Tl nd<0.2
As nd<0.2	Cr nd<0.02	Hf nd<0.02	Mn nd<0.02	Pr nd<0.2	Bi N/A	Tm nd<0.02
Au nd<0.02	Cs N/A	Hg *	Mo nd<0.02	Pt nd<0.2	Sm nd<0.2	U nd<0.2
B nd<0.02	Cu nd<0.02	He nd<0.02	Nb nd<0.02	Rb N/A	Sb nd<0.02	V nd<0.02
Ba nd<0.02	Dy nd<0.02	Ir nd<0.2	Nd nd<0.02	Ra nd<0.2	Sr nd<0.02	W nd<0.2
Be nd<0.02	Er nd<0.02	La nd<0.02	Rh nd<0.2	Ta nd<0.2	Y nd<0.02	
Bi nd<0.2	Eu nd<0.02	K nd<0.2	Ru nd<0.02	Tb nd<0.02	Yb nd<0.02	
Ce nd<0.02	Fe nd<0.02	Li nd<0.02	Sb nd<0.2	Th nd<0.02	Zn nd<0.02	
Cd nd<0.02	Ga nd<0.02	Lu nd<0.02	P N/A	Th nd<0.02	Zr nd<0.02	

This Certified Reference Material was verified in accordance with ISO/IEC 17025 (A1-1308) and ISO 17034 (AR-1403).
The solution was assayed gravimetrically, using a balance calibrated against weight sets, ID #66270, traceable to NIST.
A product with a purity of 1A, 2B, etc. is not a reference material. It is not intended to be used as a reference material without the full traceability and documentation of its production date, storage, and use conditions.
This product contains mercury and MUST be disposed of in accordance with all federal, state and local regulations.
The gravimetric uncertainty for this product is 0.24%.
In order to verify the concentration(s), the final solution was checked by plasma emission spectrometry (ICP) against material traceable to the above listed NIST SRMs.
We use the highest purity raw materials available to minimize impurity levels in the final solution. Typically 99.999% pure starting materials are used as well as high purity acids and ASTM Type 1 18 megohm deionized water.
All trace level elements impurities were determined via plasma emission spectrometry on the concentrate.
All weights are traceable through NIST. Test No. 654291344-18 & 654292805-18.
All glassware used in preparation is Class A.
All bottles are acid washed and triple rinsed with deionized water prior to use.
Shake bottle prior to use and do not pipette directly out of the bottle. Use only cleaned Class A volumetric glassware. Keep bottle tightly capped.

Certified By: *[Signature]*
Vongyuth Chanphong, Manager of Quality

Page 1 of 1

For use in routine laboratory analysis

AccuStandard is accredited to ISO 17034, ISO/IEC 17025 and certified to ISO 9001:2015

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Delivering Growth - in Asia and beyond

1. Quality Standards:

ISO 17034:2016 - General Requirements for the Competence of Reference Material Producers
ISO/IEC 17025:2017 - General Requirements for the Competence of Testing And Calibration Laboratories
ISO 9001:2015 - Quality Management System - Requirements
Eagle Registrations

2. **Intended Use:** The product covered by this certificate is designed for calibration or for use in quality control procedures for the specified chemical compounds listed on the reverse side. This product can be used for quantification and/or identification. This product can also be used as a reference material to validate analytical procedures, subject to the conditions under Section 7.

3. **Manufacturing:** All balances are calibrated daily using an in-house procedure with weights that are compared annually to master weights and traceable to NIST. The balances are also calibrated annually by an ISO/IEC 17025 accredited calibration laboratory. Please refer to the NIST test number listed on the front of this certificate. Class A glassware is used in the manufacture and quality control of all standards. Good Laboratory Practices have been used throughout the preparation of this Standard.

4. **Homogeneity:** This product is sufficiently homogeneous and any sample size would be within the uncertainty budget.

5. **Stability:** The manufacturer guarantees the stability of this solution through the expiration date stated on the label, when handled and stored according to the conditions stated on the label.

6. **Uncertainty:** The uncertainty values as stated on the face of this certificate have been determined using the EURACHEM/CITAC Guide. We report a combined expanded uncertainty equal to the positive square root of the total variance of the uncertainty of the components using the following formula: $u_c = \sqrt{u_1^2 + u_2^2 + u_3^2 + u_4^2 + u_5^2 + u_6^2 + u_7^2 + u_8^2 + u_9^2 + u_{10}^2}$. This formula represents uncertainty components from the mass, volume, short-term stability, long-term stability and homogeneity factors associated with the production of this product. The expanded uncertainty assumes a normal distribution and a coverage factor of k=2 is chosen using approximately a 95% confidence level.

7. **Legal Notice and Limit of Liability:** This product is for routine laboratory analysis and research purposes only. The company's liability will be limited to replacement of product or refund of purchase price. Notice of claims must be made within thirty (30) days from date of delivery.

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Delivering Growth - in Asia and beyond

Test Report

Customers	United Analyst and Engineering Consultant Co., Ltd.		
Equipment	Continuous Flow Analyzer	Manufacturer	SKALAR
Controller Model	SA5000	Auto Sample Model	SA1052
Controller Serial No.	182688	Auto Sample Serial No.	181729
Date of test	23-May-2025	Period	12 Month
Environment temperature	24.3 °C	Humidity	54.4 %RH

Results

Item	Characteristic	Before	After	Remark
1	Visual inspect	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
2	Power supply (210 - 240 VAC)	220 VAC	220 VAC	
3	Computer	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
4	Program	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
5	Auto sampler	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
6	Module holder			
	- Motor pump	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Pump tube	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Air-injection	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Chemistry manifolds, Switching valve, Coli, Membrane	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
7	Detector			
	- Filter	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Flow cell	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Lamp	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
8	Interface	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
9	Rinsing valves	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
10	Temperature / Reactor	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
11	Flame photometer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
12	UPS / Stabilizer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	

Warning and Error Checked

Item	Event	Before	After
13	Error list	☐ None ☐ Appear	☐ None ☐ Appear

Check with Standard

Item	Characteristic	Before	After	Remark
14	Base Line Test	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
15	Detector Signal Test	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	

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
Test By : 
(Mr. Yongyuth Champhong)

Approved by : 
(Mr. Eknasong Wanliang)


Position : Supervisor, Technical Service

Position : Manager, Technical Services


Base Line Test : Reagent_Baseline_CN_Phenol_NH3


Detector Signal Test : Phenol


FAV4.: Result Window

Detector Signal Test : NH3

FAV4.: Result Window


DKSH Technology Limited (Head office)
 2533 Sukhumvit Road, Bangkok, Phrakhanong, Bangkok, 10260
 Phone +66 2 639 7000, Mobile +66 91 813 8581, yongyuth.y@dksh.com, www.dksh.com

Delivering Growth – In Asia and Beyond

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
Page 5/2

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
 534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL 0-2717-3000-29 FAX 0-2719-9484

Certificate of Calibration

Cert.No.: 25CH166
 Page: 1 of 3

Equipment: Conductivity Meter
 Manufacturer: Horiba
 Model: LAQUA-EC210
 Serial No.: HCSL0011
 ID No.: UAE.EFM.009/2583(EFM.SCT.03/63)
 Condition As-Received: Used Item
 Received Date: 04 February 2025
 Calibration Date: 05 February 2025
 Reference: 2502-0107WSC-2
 Submitted by: United Analyst and Engineering Consultant Co., Ltd.
 3 Soi Udomsuk 41, Sukhumvit Road,
 Bangkok, Phrakhanong, Bangkok 10260

Ambient Temperature: (25 ± 2.5) °C
 Relative Humidity: (50 ± 15) %
 Calibration Procedure: In-house method:
 - CP-CH6 by direct measurement
 with certified reference material (CRM)
 - CP-CH8 by comparison with temperature standard

Calibrated by: Warakorn Lemgatrakul

Approved by:

() Chakrit Waewwanjua
 () Ponpan Paipim
 (✓) Sailthip Meangmal

Issue Date:

06 February 2025

The Uncertainties are for a confidence probability of approximately 95%

This certificate may not be reproduced other than in full, except with the prior written
 Approval of the head of Corporate Services 3: Equipment Calibration and Testing Services.

Cert.No.: 25CH166
 Page: 2 of 3

Condition of this result of calibration

1. Reference Standard Instrument :-

Instrument	Serial No.	ID No.	Certificate No.	Due date
1) Thermometer	1963678	130RC096	241995	09 Sep 2025
2) Ref. Std. Thermometer	4982054	110RC044	241757	14 July 2025

- This Certification is traceable to SI Through Technology Promotion Association (Thailand - Japan)

2. Certified Reference Materials :-

- Conductivity calibration solution, CPA chem Ltd., The measurement results are traceable to SI
 through CPA chem Ltd., ANSI-ASQ National Accreditation Board, Accredited No. AR-1835

Conductivity Solution	Manufacturer	Lot No.	Exp. date
1412.9 µS/cm	CPA Chem	1005307	15 June 2025
12.881 mS/cm	CPA Chem	1005308	15 June 2025

- Control Conductivity calibration solution temperature by Water bath (25 ± 0.1) °C

3. This certificate is valid only to the item calibrated on date and place of calibration.

Calibration results

Function: Conductivity Measurement

(*) After Adjustment at 1412.9 µS/cm

Conductivity Electrode Serial No.: 9B9F0045

Standard Conductivity Solution	Before Adjustment UUC* Reading	After Adjustment UUC* Reading	Uncertainty of Measurement (±)	Coverage factor k
1412.9 µS/cm	1455 µS/cm	1412 µS/cm	9.2 µS/cm	2.00
12.881 mS/cm	13.05 mS/cm	12.66 mS/cm	0.086 mS/cm	2.00

Remark: - UUC* = Unit Under Calibration

Cert.No.: 25CH166
Page.: 3 of 3

Calibration Results

Function : Temperature Measurement

This equipment was connected with Temperature Probe;

- Model : 9383
- Serial No. : 9B9F0045
Dimension of probe;
- Length : 110 mm
- Diameter : 16 mm
- Immersion Depth : 100 mm

Calibration Result : Without adjustment

Calibration Point (°C)	Standard Temperature (°C)	UUC* Reading (°C)	Error (°C)	Uncertainty of Measurement (± °C)	Coverage factor k
15.0	15.003	15.0	-0.003	0.13	2.00
30.0	30.004	30.0	-0.004	0.13	2.00
45.0	45.003	45.0	-0.003	0.13	2.00

Remark : - UUC* = Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

-o0o-

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
TEL 0-2717-3000-29 FAX 0-2719-5484

ILAC-MRA

Certificate of Calibration

Cert.No.: 25CH588
Page.: 1 of 3

Equipment : pH Meter
Manufacturer : Horiba
Model : LAQUA-PH210
Serial No. : HA0D0082
ID No. : UAE.EFM.072/2564(EFM.pH.05/54)
Condition As-Received: Used Item
Received Date : 20 May 2025
Calibration Date : 21 May 2025
Reference : 2505-0602WSC-3
Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
3 Soi Udomsuk 41, Sukhumvit Road, Bangkok,
Phra Khanong, Bangkok 10260
Ambient Temperature : (25 ± 2.5) °C
Relative Humidity : (50 ± 15) %
Calibration Procedure : In-house method :
- CP-CH5 by direct measurement with DC voltage
standard and direct measurement with
certified reference material (CRM)
- CP-CH8 by comparison with temperature standard
Calibrated by : Watalek Sinithean
Approved by :
Saitthip Meangmai
Issue Date : 23 May 2025

The Uncertainties are for a confidence probability of approximately 95%

This certificate may not be reproduced other than in full, except with the prior written
Approval of the head of Corporate Services 3 : Equipment Calibration and Testing Services.

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Cert.No.: 25CH588
Page.: 2 of 3

Condition of this calibration result

1. Reference Standard Instrument

Instrument	Serial No.	ID No.	Cert. No.	Due Date
1) Document Process Calibrator	54030049	130RC116	24E2759	25 Aug 2025
2) Ref. Standard Thermometer	4982054	110RC044	24I757	14 July 2025

- This measurement result is traceable to SI through Technology Promotion Association (Thailand - Japan)

2. Certified Reference Materials

The measurement results are traceable to SI through Hoch Lenge GmbH Ltd.,
Deutsche Akkreditierungsstelle, Accredited No.D-RM-15184-01-00
The measurement results are traceable to SI through CPA chem Ltd.,
ANSI-ASQ National Accreditation Board, Accredited No. AR-1635

Buffer Solution	Manufacturer	Lot No.	Exp. date
pH 4.007	CPA chem	1066665	18 Jan 2027
pH 7.000	Hach Lenge GmbH	C03232	02 Dec 2026
pH 10.010	CPA chem	1066669	18 Jan 2026

3. This certificate is valid only to the item calibrated on date and place of calibration.

Calibration Results

Function : mV Measurement

Performing standard curve by Document Process Calibrator at pH (4.7)(7.10)

Unit Under Calibration	Nominal Value	Standard Voltage Input		Actual Reading		Uncertainty of Measurement (± mV)	Coverage factor k
	pH	mV	mV	mV	pH		
pH Meter S/N.: HA0D0082	4.00	177.48	177.3	4.01	0.058	2.00	
	7.00	0.00	0.0	7.00	0.058	2.00	
	7.00	0.00	0.0	7.00	0.058	2.00	
	10.00	-177.48	-177.3	10.01	0.058	2.00	

Cert.No.: 25CH588
Page.: 3 of 3

Calibration Results

Function : pH Measurement

Performing three buffers standard curve by using buffer nominal pH (4.7)(7.10)

Unit Under Calibration	Standard pH Buffer Solution	Actual pH Reading	Actual mV Reading (mV)	Uncertainty of pH Measurement (±)	Coverage factor k
pH Electrode S/N.: Q9AA0036	4.007	4.01	160.4	0.0086	2.05
	7.000	7.00	-13.5	0.0095	2.00
	7.000	7.01	-12.6	0.0096	2.00
	10.010	10.00	-186.0	0.0092	2.00

Function : Temperature Measurement

(*) Without adjustment

This equipment was connected with Temperature Probe;

- Model : 9652-10D
- Serial No. : Q9AA0036
Dimension of probe
- Length : 103 mm.
- Diameter : 16 mm.
- Immersion Depth : 80 mm.

Calibration Point (°C)	Standard Temperature (°C)	UUC* Reading (°C)	Error (°C)	Uncertainty of measurement (± °C)	Coverage factor k
15.0	15.002	15.0	-0.002	0.13	2.00
30.0	29.999	30.0	0.001	0.13	2.00
45.0	45.001	45.0	-0.001	0.13	2.00

Remark : - UUC* = Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

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Equipment : Hot Air Oven
 Condition As-Received : Used Item
 Reference : 2503-04370C-3
 Result of Calibration :- (*) Without Adjustment
 Function of UUC* : Temperature Source
 Fresh air setting : Close

Cert. No.: 25TM579
 Page : 3 of 3

Calibration Point (°C)	UUC* Setting (°C)	UUC* Reading (°C)	Temperature stability (± °C)	Temperature uniformity (°C)	Overall Variation (°C)	Coverage Factor k
104.0	104.0	104.0	0.040	0.43	0.78	2
120.0	120.0	120.0	0.64	1.3	1.6	2
180.0	180.0	180.0	0.49	1.5	1.8	2

Calibration Point (°C)	Measured Temperature (°C)									Uncertainty (± °C)
	1	2	3	4	5	6	7	8	9 (ref.)	
104.0	104.335	104.135	104.383	104.317	103.849	103.738	104.179	104.229	104.025	0.42
120.0	119.575	119.366	119.807	119.905	118.994	119.194	119.888	119.994	120.084	1.1
180.0	180.288	179.510	180.401	180.551	179.281	179.403	180.196	180.451	180.374	1.2

Average* : The average of 30 values in each position.
 Temperature stability : One-half of the greatest maximum difference of measured temperature at any one sensor.
 Temperature uniformity : The maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time or at as close an observation time as possible to determine the temperature pattern or homogeneity within the chamber under steady-state conditions.
 Overall Variation : The Difference of the maximum and minimum measured temperatures throughout observation.
 UUC* : Unit Under Calibration
 Note : The reported uncertainty of measurement was included stability and excluded uniformity .

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3 : EQUIPMENT CALIBRATION AND TESTING SERVICES
 534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL. 0-2717-3000 FAX. 0-2719-9484

Certificate of Testing

Cert.No.: 25TW20
 Page.: 1 of 2

Equipment : DO Meter
 Manufacturer : Horiba
 Model : LAQUA-DO210
 Serial No. : HEBM0003
 ID No. : UAE.EFM.015/2563 (EFM.DO.04/63)
 Received Date : 22 January 2025
 Test Date : 23 January 2025
 Reference : 2501-0774WSC-3
 Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
 3 Soi Udomsuk 41, Sukhumvit Road, Bangkok, Phrakhanong, Bangkok 10260
 Laboratory Condition : Temperature (25 ± 5) °C
 Humidity (50 ± 20) %
 Test Procedure : In-house method : CP-CH9
 by Comparison Technique with Azide Modification Method
 Tested by : Watalek Sirithean
 Approved by :
 Approved Signatory
 () Pornhippa Tameyakul
 () Ponpan Palpim
 (✓) Salilip Meangmai
 Issue Date : 24 January 2025

Cert.No.: 25TW20
 Page.: 2 of 2

Condition of this result of calibration

1. Reference Standard Instruments :
 This certification is traceable to the International System of Unit through the reference standards laboratory of Industrial Calibration Center, Technology Promotion Association (Thailand-Japan).

Instruments	Serial No.	ID No.	Certificate No.	Due Date
1. Burette	-	130BU10	23CG1172	22 Mar 2025
2. Balance	14233821	110RC001	24MM131	04 July 2025

2. Standard Material :-

Material	Manufacturer	Lot.No.	Assay
Sodium Thiosulfate 5-Hydrate AR	KEMAUS	2203162447	99.6%

Result : Dissolved Oxygen Meter Adjustment With Air 100 %
 Dissolved Oxygen Probe No.: 9K2B0020

Titration Method (Azide Modification Method) (mg/L)	DO Meter Reading (mg/L)	Standard Deviation (mg/L)
8.20	8.20	0.0071

This report was certified only for the instrument we tested. It is allowable to use for study
 Intend to use for advertising and referral purpose is prohibited. This report may not be reproduced
 other in full, without written approval of the laboratory

-o0o-

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
 CORPORATE SERVICES 3 : EQUIPMENT CALIBRATION AND TESTING SERVICES
 534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
 TEL. 0-2717-3000-29 FAX. 0-2719-9484

Certificate of Calibration

Cert. No.: 25LM12
 Page.: 1 of 2

Equipment : DO Meter with Sensor
 Manufacturer : Horiba
 Model : LAQUA-DO210
 Serial No. : HEBM0003
 ID No. : UAE.EFM.015/2563 (EFM.DO.04/63)
 Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
 3 Soi Udomsuk 41, Sukhumvit Road, Bangkok, Phrakhanong, Bangkok 10260
 Location : TPA On Site Calibration Laboratory
 Received Order : 22 January 2025
 Calibrated Date : 23 January 2025
 Ambient Temperature : (26 ± 10) °C
 Relative Humidity : (50 ± 30) %
 AC Line Voltage : (220 ± 22) V
 Calibrated by : Warakorn Lemgagtrakul

 Approved Signatory
 Approved by :
 () Chakrit Waewwanjua
 () Suwit Imjai
 (✓) Kunchit Promprat
 Issue Date : 29 January 2025

The Uncertainties are for a confidence probability of approximately 95%

This certificate may not be reproduced other than in full, except with the prior written
 Approval of the head of Corporate Services 3 : Equipment Calibration and Testing Services.

Equipment : DO Meter with Sensor
Condition As-Received : Used Item
Reference : 2501-0774WSC-4
Procedure Used :-

Cert. No.: 25LM12
Page.: 2 of 2

Calibration were conducted using in-house calibration procedure CP-OT01 according to comparison with Industrial Platinum Resistance Thermometer (IPRT) into Temperature Bath.

The temperature scale used was based on ITS-90.

Condition of this result of calibration

1. Reference standard instrument:-

Instrument	Serial No.	Cert. No.	Traceable	Due Date
1) Digital Thermometer	2188080	2411022	TPA	17 Sep 2025

2. This certificate is valid only to the item calibrated on date and place of calibration.

3. This certification is traceable to the International System of Unit.

Remark : TPA : Technology Promotion Association (Thailand - Japan)

Result of Calibration :- (*) Without Adjustment

Function : Temperature measurement.

This instrument was connected with temperature sensor, S/N.: 9K2B0020

Calibration Point (°C)	Immersion Depth (mm)	Standard Temperature (°C)	UUC* Reading (°C)	Error (°C)	Uncertainty (± °C)	Coverage Factor k
15.0	90	15.003	15.0	-0.003	0.16	2.00
30.0	90	30.004	30.0	-0.004	0.16	2.00
45.0	90	45.002	45.0	-0.002	0.16	2.00

UUC* : Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k , providing a level of confidence of approximately 95 %.

-000-

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3 : EQUIPMENT CALIBRATION AND TESTING SERVICES
534/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
TEL. 0-2717-3000 FAX. 0-2719-9484

Certificate of Testing

Cert.No.: 25TW29
Page.: 1 of 2

Equipment : DO Meter
Manufacturer : YSI
Model : 5100
Serial No. : 11B 101863
ID No. : UAE.WAO.004/2554
Received Date : 14 February 2025
Test Date : 17 February 2025
Reference : 2502-0473DSC-1
Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
3 Soi Udomsuk 41, Sukhumvit Road, Bangchak,
Phrakhanong, Bangkok 10260
Laboratory Condition : Temperature (25 ± 5) °C
Humidity (50 ± 20) %
Test Procedure : In - house method : CP-CH9
by Comparison Technique with Azide Modification Method
Tested by : Walalak Sirihean
Approved by :
Approved Signatory
() Chakrit Waewwanjua
() Ponpan Palpim
(✓) Sathip Meangmai
Issue Date : 18 February 2025

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Cert.No.: 25TW29
Page.: 2 of 2

Condition of this result of calibration

1. Reference Standard Instruments :

This certification is traceable to the International System of Unit through the reference standards laboratory of Industrial Calibration Center, Technology Promotion Association (Thailand-Japan).

Instruments	Serial No.	ID No.	Certificate No.	Due Date
1. Burette	-	130BU10	23CG1172	22 Mar 2025
2. Balance	14233821	110RC001	24MM131	04 July 2025

2. Standard Material :-

Material	Manufacturer	Lot.No.	Assay
Sodium Thiosulfate 5-Hydrate AR	KEMAUS	2203162447	99.6%

Result : Dissolved Oxygen Meter Adjustment With Air 100 %
Dissolved Oxygen Probe No.: 24F10C202

Titration Method (Azide Modification Method) (mg/L)	DO Meter Reading (mg/L)	Standard Deviation (mg/L)
8.22	8.22	0.0055

This report was certified only for the instrument we tested. It is allowable to use for study
Intend to use for advertising and referral purpose is prohibited. This report may not be reproduced
other in full, without written approval of the laboratory

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Calibration Certificate

Certificate No.: 2502226-001-01
Client name: UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Address: 3 Sol Udomsuk 41, Sukhumvit Road,
Bangchack, Prakhonong, Bangkok 10260

Page 1 of 4

Equipment: Electronic Balance
Manufacturer: METTLER TOLEDO
Model: XSR205DU
Serial No.: C009071872
ID No.: UAE.WAO.012/2563
Order No.: 2502226
Operation No.: 2502226-001
Date of Receipt: 19 March 2025
Date of Calibration: 20 March 2025

Calibrated by Mr.Yothin Charoensuk
Scientist
Approved by *N. Nijphatt*
(Mr.Pheraphat Tuanjit)
Manager, Division of Calibration Laboratory
Responsible for the Technical Management Team
Date of Issue: 25 March 2025

The uncertainties are for a confidence probability of approximately 95%

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation Scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the units of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the National Food Institute.

FCS-009 Revision 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

Page 2 of 4

Date of Calibration: 20 March 2025
Environment Condition: Ambient Temperature: 21.2 °C Relative Humidity: 48 ± 3.5 %
Place of Calibration: 208 Balance Room, UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Condition of Equipment: Good Condition
Condition of This Results of Calibration:
1. Calibration Method: NFI Method W-MA-001 In-House Method based on UKAS Lab 14: 2019
2. Reference Standards:
Reference Standard **Model** **Serial No.** **Calibrated By** **Certificate No.** **Due Date**
Standard Weight Class E7 1mg to 200g B155557572 TCS M24041006 19 April 2025
Instrument **Model** **Serial No.** **Calibrated By** **Certificate No.** **Due Date**
Thermo-Hygro Meter 608 H1 NFI-B11H 012/23 Quality Reborn QR25-03-2 10 February 2026
3. This certificate is traceable to SI UNIT
4. This certificate was certified only for the instrument we calibrated.
5. This result of calibration was found accurate as shown on date and place of calibration only

Calibration Results:


1. Repeatability of Reading:

Nominal Value (g)	Standard Deviation of Reading (g)
40	0.0000052
80	0.0000042
160	0.0000000
200	0.0000000

2. Off-Center Error:

A mass of 100 g was placed and moved to various position on pan.

The balance reading obtained is given in the table.

1	2	3	4	5	6	(Maximum Difference)
(g)	(g)	(g)	(g)	(g)	(g)	(g)
100.0001	100.0001	100.0001	100.0001	100.0001	100.0002	0.0001

FCS-012 Revision 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

Date of Calibration: 20 March 2025 Page 3 of 4

Calibration Results: (Continued)

Calibration Range: 0-80 g

Calibration Adjustment: Internal Calibration

3. Departure from Nominal Value: (Range: 0 - 82 g; Resolution: 0.00001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (±g)	Coverage Factor
Unload	0.00000	0.00000	0.00000	0.0000089	2.00
0.001	0.001003	0.00100	0.00000	0.0000092	2.00
0.005	0.005002	0.00500	0.00000	0.0000091	2.00
0.01	0.010003	0.01000	0.00000	0.0000091	2.00
0.05	0.049995	0.05000	0.00000	0.0000098	2.00
0.1	0.100011	0.10000	0.00001	0.000011	2.00
0.5	0.500016	0.50000	0.00002	0.000014	2.00
1	1.000003	1.00001	-0.00001	0.000016	2.00
2	2.000023	2.00005	-0.00003	0.000017	2.00
5	5.000015	5.00005	-0.00003	0.000021	2.00
10	10.000009	10.00005	-0.00004	0.000026	2.00
20	20.000030	20.00012	-0.00009	0.000037	2.00
30	30.000039	30.00012	-0.00008	0.000050	2.00
50	50.000078	50.00014	-0.00011	0.000068	2.00
80	80.000067	80.00020	-0.00013	0.00011	2.00

FCS-012 Revision 01 Date: 20-04-65

Calibration Report

Certificate No.: 2502226-001-01
Equipment: Electronic Balance
Model: XSR205DU
Serial No.: C009071872
Capacity: 82 g / 220 g
Manufacturer: METTLER TOLEDO
Resolution: 0.00001 g / 0.0001 g
ID No.: UAE.WAO.012/2563

Page 4 of 4

Calibration Results: (Continued)

Calibration Range: >80-200 g

Calibration Adjustment: Internal Calibration

3. Departure from Nominal Value: (Range: >80 - 200 g; Resolution: 0.00001 g)

Nominal Value (g)	Standard Value (g)	Average Reading (g)	Correction (g)	Uncertainty (±g)	Coverage Factor
90	90.00010	90.00002	-0.00008	0.000015	2.00
100	100.00006	100.00001	-0.00005	0.000016	2.00
110	110.00007	110.00001	-0.00006	0.000017	2.00
120	120.00009	120.00002	-0.00007	0.000018	2.00
130	130.00010	130.00002	-0.00008	0.000019	2.00
140	140.00013	140.00002	-0.00011	0.000020	2.00
150	150.00009	150.00002	-0.00007	0.000021	2.00
160	160.00010	160.00002	-0.00008	0.000022	2.00
170	170.00012	170.00002	-0.00010	0.000023	2.00
200	200.00013	200.00002	-0.00011	0.000028	2.00

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor k, providing a level of confidence of approximately 95 %.

FCS-012 Revision 01 Date: 20-04-65

Calibration Certificate

Certificate No.: 2503287-002-01
Client name: UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Address: 3 Soi Udomsuk 41, Sukhumvit Road,
Bangchack, Prakhonong, Bangkok 10260

Page 1 of 3

Equipment: CHAMBER (Incubator)
Manufacturer: BINDER
Model: KB 400
Serial No.: 20220000000391
ID No.: UAE.MIC.029/2565
Order No.: 2503287
Operation No.: 2503287-002
Date of Receipt: 5 June 2025
Date of Calibration: 5 June 2025

Calibrated by Mr. Phraphat Tuanjit
Scientist
Approved by P. Jengphat
(Miss Phraphat Jaengkarnkit)
Vice President, Department of Laboratory Services
Responsible for the Technical Management Team
Date of Issue: 11 June 2025

The uncertainties are for a confidence probability of approximately 95 %.
This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the units of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full extent with the prior written approval of the National Food Institute.

FCS-009 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2503287-002-01
Equipment: CHAMBER (Incubator)
Model: KB 400 Serial No.: 20220000000391
Resolution: 0.1 °C ID No.: UAE.MIC.029/2565
Manufacturer: BINDER

Date of Calibration: 5 June 2025 Page 2 of 3

Location: Room 302 Microbiology Laboratory, UNITED ANALYST AND ENGINEERING CONSULTANT CO.,LTD.
Environment Condition:
Ambient Temperature (18 ± 1) °C
Relative Humidity (53 ± 6) %
Line Voltage (230 ± 5) Volt

Condition of this results of Calibration:

- This instrument was calibrated by insert 13 standard thermometer into its chamber and calibration according to WTE-014 Based on TIAS G-20-1/02-06 (E); Guidelines for Calibration and Checks of Temperature Controlled Enclosures.
- The temperature scale used was based on ITS - 90.
- All data show below were final values and the initial data may be obtained upon request.

2. Reference Standard Instrument:

Instrument	Model	Serial No./ID No.	Certificate No.	Due Date	Through
Digital Thermometer with sensor	34972A	MY5903377	2501166-001-01	13 January 2026	NATIONAL FOOD INSTITUTE
RTD		(Our lot: 203 / RTD+101 203)			

- This certificate is traceable to International System of Units (SI Units).
- This certificate was certified only for the instrument we calibrated.
- This result of calibration was found accurate as shown on date and place of calibration only.
- Condition of Calibrated item: Good

UUC Description:

Time of Record: 1 Hour 9 Minute At: 35.0 °C
Fresh Air Damper: - Open Position -
X Close Fan 100%
- Not Available

- Result of Calibration: X Without adjustment After adjustment

FCS-012 Revision: 01 Date: 20-04-65

Calibration Report

Certificate No.: 2503287-002-01
Equipment: CHAMBER (Incubator)
Model: KB 400 Serial No.: 20220000000391
Resolution: 0.1 °C ID No.: UAE.MIC.029/2565
Manufacturer: BINDER

Date of Calibration: 5 June 2025
Calibration point: 35.0 °C Page 3 of 3

Calibration Condition	Temperature (°C)	Relative Humidity (%)	Line Voltage (Volt)
MIN	17.4	48	225.0
MAX	18.5	59	235.0

Table 1: Reporting of Temperature

Calibration point (°C)	Measured Temperature (°C) @ Sensor No. (Sensor No.13 is REF)													Uncertainty ± (°C)
	#1	#2	#3	#4	#5	#6	#7	#8	#9	#10	#11	#12	#13	
35.0	35.08	35.11	35.01	35.13	35.17	35.09	34.98	34.89	35.15	35.05	35.06	34.89	35.06	0.27

Table 2: Reporting of Characterization Result

UUC Setting (°C)	UUC Reading (°C)			Temperature Stability ± (°C)	Temperature Uniformity (°C)	Overall Variation (°C)
	MIN	MAX	Average			
35.0	35.0	35.1	35.0	0.038	0.17	0.34

Note: The quoted uncertainty include " Stability " and " Loading effect (20% of Temp Uniformity) " .
UUC = Unit Under Calibration
Stability = One-half of the greatest maximum difference of measured temperatures at any one sensors, for at least half an hour after reaching steady state.
Uniformity = The maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time.
Overall Variation = The difference of the maximum and minimum measured temperatures throughout observation time.
The report uncertainty of measurement was based on standard uncertainty multiplied by coverage factor k=2, providing a level of confidence of approximately 95 %.

FCS-012 Revision: 01 Date: 20-04-65

List of Instrument Certificates for Environmental Quality Analysis

No.	Instrument/Equipment	Parameter	Manufacturer	Model/Serial No.	Calibrator	Certification No.	Date of Calibration	Due date of Calibration
1	Atomic Absorption Spectrometer	ARSENIC	Agilent Technologies	AA200PS / MY13100001	Agilent Technologies (Thailand) Co.,Ltd	Preventive Maintenance	30/1/2025	29/1/2026
2	Analytical Balance	SUPERFIND SOLIDS	Mettler Toledo	XS2050N / CO0971812	National Food Institute Ministry of Industry, Thailand	2502235-001-01	20/2/2025	19/2/2026
3	Continuous Flow Analyser(CFA)	TOTAL AMMONIA	Solar Analytical	Sen+5000-02 / 182688	DKSH (Thailand) Ltd	WD-00074079	23/5/2025	22/5/2026
4	Cold Vapor Atomic Fluorescence Spectrometer (Mercury Analyser)	MERCURY	Analyst's Jena AG	Mercur duo pur / K 7040153	Analyst's Jena GmbH Thailand Ltd.	91400	3/2/2025	2/2/2026
5	DO Meter	DO	Horba	LAQUA-DO210 / HEM0003	Technology Promotion Association (Thailand-Japan)	25TV20	23/1/2025	21/1/2026
6	DO Meter	BIOCHEMICAL OXYGEN DEMAND	YSI	5100 / 118 101863	Technology Promotion Association (Thailand-Japan)	25TV29	17/2/2025	16/2/2026
7	SCF Meter	SALINITY	Horba	LAQUA-EC10 / HCS1011	Technology Promotion Association (Thailand-Japan)	25CH166	5/2/2025	3/2/2026
8	For Air Oven	SUPERFIND SOLIDS	Mettler	UF55 / 82124011	Technology Promotion Association (Thailand-Japan)	25TM679	19/2/2025	18/2/2026
9	Incubator	FECAL COLIFORM BACTERIA	Binder	43400 / 20220000000391	National Food Institute Ministry of Industry, Thailand	2503287-002-01	5/6/2025	4/6/2026
10	pH Meter	pH	Horba	LAQUA-PH10 / HAD00682	technology promotion association (Thailand-Japan)	25CH188	21/2/2025	20/2/2026
11	UV-Vis Spectrophotometer	NITRATE NITROGEN PHOSPHATE PHOSPHORUS	Haashi	U-2000 / 21222-009	DOE Services Co.,Ltd	5945-021	26/5/2025	25/5/2026
12	Turbidity Meter (Portable)	TURBIDITY (NTU)	Dicon Instrument(China)	T100F / 1120501017	Technology Promotion Association (Thailand-Japan)	25CH1053	9/1/2025	7/1/2026

Due Date of Calibration*: Based on the annual calibration plan. At least 1 time per year.

Certificate Page 1 of 1

United Analyst and Engineering Consultant Co., Ltd. (UAE)
Certified Laboratory ISO/IEC 17025

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Serial-No.: K170A0153 Customer-No.: C04-006
Date: 3 February 2025 Carried out by: Mr. Srichai Fak-On

Maintenance with following Operational Qualification (OQ)
(requires a separate OQ protocol)

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Document is not controlled
Document is not controlled

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Maintenance works basic unit

tightness visual check inside the Mercur ☒
visual check if gold-traps are broken ☒
visual check if spectrometer is contaminated ☒
visual check of the fluorescence cell ☒
visual check of the absorption cell, incl. window ☒
reactor cleaning ☒
check pump-hose, if necessary change it ☒
check swivel drive (SEV) ☒
check drying-hose, output gas-liquid-separator ☒
test Bubble-Sensor ☒
check gas flows ☒
check volume flows, reagents ☒
recording stray light values ☒
measurement with 30 ng/l ☒

Maintenance works Autosampler

Serial No.: N/A

lubricate the dosing-winding (Teflon-grease-spray) ☒
clean the dosing cylinder, if necessary exchange it ☒
lubricate the winding system of the height drive with some drops of oil ☒
check the toothed belt ☒
check the position of the mechanical stopper (height: 13mm) ☒
check the pump rate of mixing pump (<14s AS52, typ.7s/<20s AS52S, typ.10s) ☒
check the pump rate of washing cup ☒
check the electrical hose connections for good contact ☒
check the connectors of the magnetic valves ☒
check the dosing hose for buckling, if necessary exchange it ☒

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Services are subject to the General Terms and Conditions of Analytik Jena AG, which will be sent on request.

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3/02/2025 13:44 Page 1/3

Mercur

Report file: C:\WINAAS\TMP\2025result\WO\Pro_039
Program version: 4.7.9.0 Printed on: 3/02/2025 13:44
Recording started on 3/02/2025 13:32 GMT+7.0
Operator:
Laboratory:
Code:
Remarks:

Method parameters

Method Without enrichment / FBR 100 ng/L PM_12-02-2024
Created on 12/02/2024 Time 11:54
Program ---

Parameters Mercur Technique: Hg absorption

Line	253.7 nm		
Lamp type	Hg-LP		
Integr. mode	Peak height	Integr. time	40 s
PMT	236 V	Peak smoothing	12/5
AZ time	5 s		
Delay	0 s		
Working mode	w/o enrich.	System cleaning	Acid
FBR technique	off	Wash time acid	15 s
Pump speed	4	Soaking time	20 s
Sample load time	8 s	Gas load time	10 NL/h
Reaction time	12 s		
Waiting time AZ	15 s		
Purge time 1	40 s		
QC parameters			
QC type	Conc. check	QC check samp. 2	---
QC check samp. 1	---	Conc.	---
Conc.	---	Error limit	---
Error limit	---	Reaction	flag + continue
Rep. measurement	off	QC std.2 no.	1(100.00 ng/L)
QC std.1 no.	1(100.00 ng/L)	QC std.2 limit	± 0.00%
QC std.1 limit	± 50.00%	Reaction	flag + continue
QC std. act.	flag + continue	Reaction	off
Expect. blank abs.	0.0100± 0.0100	QC Recal.factor	Off
QC precision	off		

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3/02/2025 13:44 Page 2/3

Calibration settings

Calib. meth	Standard calib.	Calibr. unit	ng/L
No. standards	1	Conversion fac.	1000000
Type of standards	---	Standard prep.	Premixed
		Blank correct.	---
		Recalibr. std. no.	---
Output unit	µg/L	Conversion fac.	1000
Calib. stat.	Mean	Meas. cycles	3
		Blind cycles	1
Stock sol. 1	---	Stock sol. 2	---
Stock sol. 3	---	Stock sol. 4	---
Type of cal. curve	linear	Intercept	calculated
Weighted cal.	off	Grubbs stat.	off
Check of cal. curve	no outlier test		

Sample statistics

Stat. mode	Mean	Meas. cycles	2
Confid. level	95.4 %	Blind cycles	1
Grubbs stat.	---		

Calibration standards

No	Name	State	Pos	Conc./ ng/L	Abs	SD	RSD/%
1	Cal-Zero	(-)	##	0.00	H: 0.000544 A: 0.005800	0.000115 0.004748	21.26 81.87
2	Cal-Std1	(-)	##	100.00	H: 0.003251 A: 0.042341	0.000094 0.003312	2.921 7.824

Calibration function 1 3/02/2025 13:43 Calibration (Peak height)

Abs=k1+k2*conc


k1=0.000544 k2=0.000027

Recal. factor: ---

Slope	0.00003 Abs/(ng/L)	R2-adjusted	1.0000
sc0	1.00000 ng/L	Charact. conc.	161.087 (ng/L)/1%
Lower limit	0 ng/L	Upper limit	110. ng/L
Detection limit	---	Deter. limit	---

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

Measurements and events (sorted by time)

Hg	Without enrichment / FBR 100 ng/L PM_12-02-2024					3/02/2025	13:32
ID	Conc.	Abs	BG	SD	RSD/%	Int. type	Time
Cal-Zero		0.000564				PkH	13:35
		0.000420					13:36
		0.000649					13:37
	0ng/L	0.000544		0.00011577	21.26		13:37
Cal-Std1		0.003268				PkH	13:40
		0.003336					13:41
		0.003148					13:43
	100.ng/L	0.003251		0.000094975	2.921		13:43
Calibration	Calibration function: 01						13:43
Peak plots						Hg	

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Mercur

Report file: C:\WinAAS\TMP\2025\result\WO1Pro_040
 Program version: 4.7.9.0 Printed on: 3/02/2025 17:37
 Recording started on 3/02/2025 17:26 GMT+7.0
 Operator:
 Laboratory:
 Code:
 Remarks:

Method parameters

Method Without Enrichment / FBR / 30 µg/L_PM_3-02-2025
 Created on 3/02/2025 Time 10:33
 Program ---

Parameters Mercur Technique: Hg fluorescence

Line	253.7 nm		
Lamp type	Hg-LP		
Integr. mode	Peak height	Integr. time	35 s
PMT	451 V		
AZ time	5 s	Peak smoothing	12/5
Delay	0 s		
Working mode	w/o enrich.	System cleaning	Off
FBR technique	on	Wash time acid	10 s
Pump speed	3	Soaking time	20 s
Sample load time	12 s	Gas load time	10 NL/h
Reaction time	12 s		
Waiting time AZ	5 s		
Delay	0 s		
Purge time1	30 s		
Purge time2	15 s	Gas wash time2	10 NL/h

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QC parameters

QC type	Conc. check	QC check samp. 2	---
QC check samp. 1	---	Conc.	---
Error limit	---	Error limit	---
Rep. measurement	off	Reaction	flag + continue
QC std. 1 no.	1(30.000 ng/L)	QC std. 2 no.	3(0.100 ng/L)
QC std. 1 limit	± 20.00%	QC std. 2 limit	± 20.00%
QC std. act.	flag + continue		
Expect. blank abs.	0.0100± 0.0100	Reaction	flag + continue
QC precision	off	Reaction	off
		QC Recal.factor	Off

Calibration settings

Calib. meth	Standard calib.	Calibr. unit	ng/L
No. standards	1	Conversion fac.	1000000
Type of standards	---	Standard prep.	Premixed
		Blank correct.	---
		Recalib. std. no.	---
Output unit	µg/L	Conversion fac.	1000
Calib. stat.	Mean	Meas. cycles	3
		Blind cycles	1
Stock sol. 1	---	Stock sol. 2	---
Stock sol. 3	---	Stock sol. 4	---
Type of cal. curve	linear	Intercept	Zero
Weighted cal.	off	Grubbs stat.	off
Check of cal. curve	no outlier test		

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

Mercur

Report file: C:\WinAAS\TMP\2025\resu\W01Pro_041
 Program version: 4.7.9.0 Printed on: 3/02/2025 17:59
 Operator: Recording started on 3/02/2025 17:46 GMT+7.0
 Laboratory:
 Code:
 Remarks:

Method parameters

Hg

Method With Enrichment / FBR / 30 µg/L_PM_3-02-2025
 Created on 3/02/2025 Time 12:45
 Program ---

Parameters Mercur Technique: Hg fluorescence

Line	253.7 nm		
Lamp type	Hg-LP		
Integr. mode	Peak height	Integr. time	20 s
PMT	444 V		
AZ time	5 s	Peak smoothing	12/5
Delay	0 s		
Working mode	Enr. w/o reload.	System cleaning	Off
FBR technique	on	Wash time acid	10 s
Pump speed	3	Soaking time	20 s
Sample load time	10 s	Gas load time	5 NL/h
Reaction time	10 s		
Waiting time AZ	5 s		
Delay	0 s		
Purge time1	20 s		
Purge time2	15 s	Gas wash time2	5 NL/h
Purge time3	10 s	Gas wash time3	10 NL/h
Heat.time coll.1	20 s	Cool. time coll.1	25 s

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QC parameters

QC type	Conc. check	QC check samp. 2	---
QC check samp. 1	---	Conc.	---
Conc.	---	Error limit	---
Error limit	---	Reaction	flag + continue
Rep. measurement	off	QC std.2 no.	1(30.000 µg/L)
QC std.1 no.	1(30.000 µg/L)	QC std.2 limit	± 50.00%
QC std.1 limit	± 50.00%	Reaction	flag + continue
QC std. act.	flag + continue	Reaction	off
Expect. blank abs.	0.0100± 0.0100	QC Recal.factor	Off
QC precision	off		

Calibration settings

Calib. meth	Standard calib.	Calibr. unit	µg/L
No. standards	1	Conversion fac.	1000
Type of standards	---	Standard prep.	Premixed
		Blank correct.	---
		Recalib. std. no.	---
Output unit	µg/L	Conversion fac.	1000
Calib. stat.	Mean	Meas. cycles	3
		Blind cycles	1
Stock sol. 1	---	Stock sol. 2	---
Stock sol. 3	---	Stock sol. 4	---
Type of cal. curve	linear	Intercept	Zero
Weighted cal.	off	Grubbs stat.	off
Check of cal. curve	no outlier test		

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

Agilent 55 240 280 Series Atomic Absorption Spectroscopy Systems

Preventive Maintenance Checklist

Agilent Preventive Maintenance provides factory recommended service for your analytical systems to assure reliable operation and the accuracy of your results.

Delivered by highly trained and certified service engineers using genuine Agilent parts and supplies, Agilent Preventive Maintenance provides everything you need to reduce unplanned downtime and keep your systems operating at their peak. This checklist will be completed at the end of the service and provided to you as a record of the installation.

Note: While non-current production AA instrument and/or accessory models are not covered specifically in this document it can be used as a basic reference.

For more information about Agilent Technologies services please visit our web site using the following URL: <http://www.agilent.com/en-us/services>

Introduction

Customer Information

- 1 Customers should provide all necessary operating supplies upon request of the engineer.
- 2 A customer representative should be available to the engineer while performing the preventive maintenance procedures.
- 3 Any parts, not included in the Parts Lists section of this document, are not part of the recommended Preventive Maintenance service, nor are they included in the price of this service.
- 4 If a system requires the use of extra or special procedures and/or parts for the maintenance service, then these must be ordered separately and charged as a repair, which may incur additional costs.

Revision: 10/03, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Important Customer Web Links

- For more information about Agilent Technologies services, please visit our website using the following URL: <http://www.agilent.com/en-us/products/crosslab-instrument-services/service-repair>
- To access Agilent University, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- A useful Agilent Resource Center web page is available, which includes short videos on maintenance, quick lists of consumables for new instruments, and other valuable information. Check out the Resource Page here: <https://www.agilent.com/en-us/agilentresources>
- Need technical support, FAQs, supplies? – visit our Support Home page at <http://www.agilent.com/search/support>
- Get answers. Share insights. Build connections. Join the Agilent Community at <https://community.agilent.com/welcome>

Service Engineer's Responsibilities

- Contact the customer and ensure that all necessary supplies are available before the preventive maintenance visit.
- Confirm the ability of the instrument to deliver continued safe operation as established via the Agilent AA safe operation flow chart. (Refer directly to the AA 55/240/280 Preventive Maintenance Scope of Work to make this decision.)
- Only select those pages that relate to the system or module being serviced.
- Complete empty fields with the relevant information.
- Complete the relevant checkboxes in the checklist using either a "X" or tick mark "✓".
- Check "Section not applicable" check boxes to indicate services/tasks not delivered, as appropriate.
- Complete the Preventive Maintenance service in the order of the tasks listed.
- Complete the Service Review section together with the customer.
- Complete the fields for page numbers at the foot of each selected page.
- Complete the total number of pages field in the Service Completion section.
- Ask the customer to sign the Service Completion section including the customer's and your signature.

This information is subject to change without notice.

Revision 10/03, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Instrument Maintenance

System Information

- ☐ Check this box if an instrument configuration report is attached instead of completing the table.

Instrument System Name and ID	240 FS AAS
Instrument System Site and Location	United Analyst and Engineering Consultant

List System Component Product Numbers	List the Serial Numbers of each Component
1 G 8432 A	M1 13160001
2	
3	
4	
5	
6	
7	
8	
9	

Preparation, Safe operation and Initial performance checks

Revision 10/03, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Agilent AA safe operation flow chart inspections (to determine if the PM can be performed).

NOTE: If by following the flow chart the instrument is deemed to be unsafe for continued use you MUST NOT continue PM work. Inform the customer immediately of the Agilent recommendation that use of the instrument be discontinued.

- ☒ Discuss any specific issues with the customer before starting.
- ☐ For HF application systems, if standard sample introduction system was not installed, ask the customer to install it. **NA**
- ☒ Review the instrument logbook for recorded problems and comments.
- ☒ Save instrument control settings before starting the procedure.
- ☒ Perform a general inspection of the system for cleanliness.
- ☒ Check for proper installation of parts, assemblies, sensors etc.
- ☒ Check system for required installation of components, settings as defined by current Service Notes
- ☒ Check for required firmware updates and verify with customers if they would like them installed.
- ☒ Use SVD to perform a Full Wavelength Scan for Cu HCL - "As found test_1"
- ☒ Perform a Basic Cu ABS test - "As found test_2"
- ☒ Print the Details page or screen captures of the test results and attach to the end of this checklist.

Revision 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Preventive Maintenance Procedures

FLAME SYSTEM section

- ☐ Section not applicable

Electronic components

- ☒ Review and confirm instrument configuration data in SVD
- ☒ Confirm power supply voltages using the **SVD Power Supply diagnostic**.
- ☒ For Dual Beam Instruments - Confirm RBC frequency using the **SVD RBC frequency diagnostic**.

Mechanical components

- ☒ Check the burner adjuster controls for complete and free movement. If the burner adjuster needs lubrication, use Molykote 321 or mineral based molybdenum disulphide grease.
- ☒ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☒ Monochromator drive
 - ☒ Slit drive
 - ☒ Lamp selector
 - ☐ ABA

Optics components

- ☒ Check that external optical surfaces are clean – Clean or replace as required
- ☒ Use SVD and perform **Mono Wavelength Correction**.
- ☒ Use SVD and perform **Slit Calibration**.
- ☒ Use SVD and perform **Grating Squareness Diagnostic**.
- ☒ Use SVD and perform **Zero Order Offset/Mono Correction**.
- ☒ Use SVD and perform **Wavelength Repeatability**.
- ☒ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.
- ☒ Check that the signal energy of the D2 and HC lamps track properly. Advise customer if their D2 lamp is showing emission degradation due to age.

Revision 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Sample Introduction and Atomization

- ☒ Inspect the burner interlock plate to ensure that the interlock pin is secure and correct for the burner type.
- ☒ Clean the burner slot with a clean white card.
- ☒ Check the uniformity of the slot width.
- ☒ Clean the burner if required.
- ☒ Change the burner o-ring
- ☒ Clean the nebulizer, spray chamber and liquid trap.
- ☒ Change all o-rings and seals in the nebulizer, nebulizer block and spray chamber.
- ☒ Check that the pressure relief bung releases readily.
- ☒ Change o-rings on the fuel and oxidant delivery barbs
- ☒ Leave the liquid trap EMPTY and verify the flame will not ignite in this state.
- ☒ Refill liquid trap and check that overflow drains freely into the drain/waste tube.
- ☒ Check the drain/waste tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel
- ☒ Check and clean the igniter electrode

Gas handling components and safety interlocks

- ☒ Pressure test for leaks
- ☒ Leak test gasbox internal components and connections
- ☒ Check safety interlock status and operation using the **SVD interlock monitoring diagnostic**.

Analytical performance for Flame systems

- ☒ Ignite a flame.
- ☒ Check that you can adjust the nebulizer uptake rate from 4 to 6.5 mL per minute.
- ☒ Optimize the instrument ready to perform Cu sensitivity test.
- ☒ Create a manual method to perform a Basic Cu ABS test - "Final Performance Testing"
- ☒ Run a PM completed sensitivity test for a 5 ppm copper sample and record the results in the AA PM Performance test results and measurements table.

Revision 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

FURNACE SYSTEM section

- ☒ Section not applicable

Electronic components

- ☐ Review and confirm instrument configuration data in SVD
- ☐ Confirm power supply voltages using the **SVD Power Supply diagnostic**.

Mechanical components

- ☐ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☐ Monochromator drive
 - ☐ Slit drive
 - ☐ Lamp selector

Optics components

- ☐ Check that external optical surfaces are clean – Clean or replace as required
- ☐ Use SVD and perform **Mono Wavelength Correction**.
- ☐ Use SVD and perform **Slit Calibration**.
- ☐ Use SVD and perform **Grating Squareness Diagnostic**.
- ☐ Use SVD and perform **Zero Order Offset/Mono Correction**.
- ☐ Use SVD and perform **Wavelength Repeatability**.
- ☐ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.

Gas handling, water system and workhead component checks

- ☐ Inspect the GTA workhead gas hoses and connections for leaks.
- ☐ Pressure test for gas leaks
- ☐ If the cooler system is accessible (stand-alone) check for correct operation and coolant/water level – (this includes any temperature and pressure settings plus filter cleaning (air flow and water).
- ☐ Inspect the GTA workhead water hoses and connections for leaks.
- ☐ Check all graphite components and replace if necessary.

Revision 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Tube
- ☐ Electrodes
- ☐ Shroud

- ☐ Check and clean the end windows on the workhead.
- ☐ Check safety interlock operation.

Analytical performance for Furnace systems

- ☐ Optimize the instrument ready to perform Cu sensitivity test.
- ☐ Run the sensitivity test for a 25 ppb copper sample and record the results in the results table.

PSD autosampler accessory for Furnace systems

☒ Section NOT Applicable

- ☐ Check condition of the PSD capillary – replace if necessary
- ☐ Check condition and operation of PSD syringe – ensure it does not have air locks and bubbles
- ☐ Change PSD rinse bottle o-ring
- ☐ Check and clean the rinse vessel
- ☐ Check the drain tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel
- ☐ Ensure that the waste vessel is suitable for use with the furnace system.

Sample introduction pump system (SIPS) accessory

☒ Section NOT Applicable

- ☐ Re-torque screws securing the hubs, presser arms and pump rotors.
- ☐ Adjust each roller so that it rotates freely.
- ☐ Wipe clean the pump rotor rollers and pump bands with a dry clean cloth
- ☐ Ensure that the presser arms and the surfaces near the pump are free from dirt and spills
- ☐ Remove the pump module rear cover and check for the incursion of liquids and any signs of corrosion.
- ☐ Re-torque the nuts that fasten the motor mounting plates to the chassis.
- ☐ Check clips securing the diluents holder and replace if necessary.
- ☐ Disconnect, clean T piece, and reassemble the tubing using the following steps.

Revision: 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

- ☐ Remove the T-piece by disconnecting the pump tubes, the pump bands and all other tubing.
- ☐ Place the T-piece in an ultrasonic bath containing strong detergent 1-5% Decon 30 or similar, for approximately 5-10 minutes
- ☐ Wash the T-piece under a tap with a strong flow of water.
- ☐ Rinse with distilled water through all of the inlets in the reverse direction to normal sample flow.
- ☐ Reassemble.

Sample preparation system (SPS 4) accessory

☒ Section NOT Applicable

The Agilent SPS 4 autosampler is designed to need minimal maintenance.

The following maintenance requirements are suggested to maintain the performance of the autosampler.

- ☐ Cleaning the spill tray, rack location mat, end frames and chassis accessories with a damp soft cloth and diluted mild detergent.
- ☐ Cleaning the autosampler cover panels with domestic window cleaner.
- ☐ Checking the X- axis and Z- axis drive belts for cracks, splits, damaged teeth, excessive fraying, color changes or degradation from fumes.
- ☐ Check the X- axis, Theta- axis and Z- axis FFC cables for cracks, incorrect positioning, damaged edge or damaged connectors.

NOTE: The autosampler requires no extra lubrication throughout its lifetime.
For further details refer to the SPS 4 service manual 68410-90050.

Sample preparation system (SPS 3) accessory

☒ Section NOT Applicable

- ☐ Check the x-axis and z-axis timing belts – Replace if there is any cracks, splits or color deterioration and belt tension.
- ☐ Check belt tensions - adjust if required
- ☐ Check the lubrication pad for single x-axis shaft. If pad is dry or customer has observed any vibration or erratic movements of the x-axis carriage, add 1 mL of Dow Corning 200 @ Fluid, 200 CS into the well.
- ☐ Check the auto-sampler ability to find tube positions - Calibrate if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Revision: 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

Vapor generation accessory VGA (hydride generator)

☒ Section NOT Applicable

- ☐ Inspect VGA gas supply hose.
- ☐ Inspect/replace VGA pump tubing
- ☐ Check low gas pressure interlock setting – adjust if required.
- ☐ Check precision orifice gas flow setting – adjust if required.
- ☐ Check gas regulator pressure to 45 psi (325 kPa) – adjust if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

UltraAA lamp accessory (external)

☒ Section NOT Applicable

- ☐ Check the condition of the power cable.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Restore System

- ☐ If you have altered the customer's instrumentation during the course of PM, restore to the original status to allow the customer to conduct their normal activities (e.g., reload the customer's method.)

Guidance

If the PM service is performed prior to a qualification service, then use the qualification procedure as a guide for final instrument set up and checkout.

Revision: 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

Signature Page

Service Review

- ☒ Attach available reports/printouts of all tests to this documentation.
- ☒ Record the Preventive Maintenance service activity in the customer's records/logbook.
- ☒ Update/reset instrument maintenance counters as appropriate.
- ☒ Affix the PM sticker to the system or instrument logbook based on the customer's request.
- ☒ Complete the Service Engineer Comments section if there are additional comments.
- ☒ Review this service, parts replaced, and test results obtained with the customer.
- ☒ If the instrument firmware was updated, record the details of the change in the Service Engineer's Comments box or if necessary, in the customer's IQ records.

Test Results

Test Description	Expected Test Result	Actual Test Result
Flame optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	49 /
Flame performance test with 5 ppm copper sample		
Air/acetylene, mixing paddle removed	Abs value > 0.5	0.559%
Air/acetylene, mixing paddle installed, 10 replicates	%RSD < 1.0	0.2 /
Deuterium furnace optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	—
Deuterium furnace performance test with 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0 %	—
Abs value	> 0.15	—
Zeeman furnace analytical performance: 25 ppb copper sample (327.4 nm)		
Precision %RSD	≤ 4.0 %	—
Abs value	≥ 0.10	—
MSR%	≥ 70 %	—

Revision: 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

AA consumable and parts list table

Part Description	Part Number	Product/Model #s where used	PM supplied or Consumable	Instrument Type
Test Solution - Cu ppm solution	6610030100	50 55 140 240 280	PM supplied	Common
Test Solution - Blank solution	5190-7001	50 55 140 240 280	PM supplied	Common
Copper, 1000 ug/ml, 100ml	5190-8279	50 55 140 240 280	*	Common
Kit, Mk 7 O-rings, aqueous, complete set	9910093400	50 55 140 240 280	PM supplied	Flame
Organic Kit	9910093500	50 55 140 240 280	PM supplied	Flame
Wire Nebulizer Cleaning	9910024700	50 55 140 240 280	consumable	Flame
Tubing-Capillary Std Neba	9910024900	50 55 140 240 280	consumable	Flame
Capillary Tube Hvac Ncb (3) (organics only)	9910044000	50 55 140 240 280	consumable	Flame
Glass impact beads (5/pk)	9910025700	50 55 140 240 280	consumable	Flame
Teflon impact beads (5/pk): (organics only)	9910053300	50 55 140 240 280	consumable	Flame
Burner cleaning strip (100/pk)	9910053900	50 55 140 240 280	consumable	Flame
Window UV silica - round (right side)	2010082600	50 55 140 240 280	PM supplied	Common
Window UV silica - rectangular (left side)	2010082500	50 55 140 240 280	PM supplied	Common
Pad adhesive window (round)	4910012700	50 55 140 240 280	PM supplied	Common
Pad adhesive window (rectangular)	4910012800	50 55 140 240 280	PM supplied	Common
Electrode kit (1 pr) (D2)	6310003400	GTA120	PM supplied	Furnace
Shroud (D2)	6310003100	GTA120	PM supplied	Furnace
Zeeman electrode kit (1 pr)	6310003500	GTA120	PM supplied	Furnace
Zeeman shroud	6310003600	GTA120	PM supplied	Furnace
O-ring PSD rinse bottle	6910025900	PSD120	PM supplied	Furnace

* For engineers who only service AA instruments 5190-8279 can be used as a cheaper alternative for 6610030100.

Items classified as PM supplied in the above table are included in the standard PM. Those classified as consumable should be provided by the customer or charged to the customer if supplied by the Agilent service engineer.

Revision: 10/00, Issued: November 2021

© Agilent Technologies Inc. 2021

Service Engineer Comments (optional)

If there are any comments, please provide them in this field, along with the instrument or user's name, to help us better understand the problem and solve it.

Service Completion

Service request number: 6007049143

Date service completed: 30 Jan 2025

Agilent signature: Kanyakorn S.

Customer signature: Sunda Y.

Total number of pages in this document: 13

Revision: 10/00, Issued: November 2021

© Agilent Technologies Inc. 2021

SVD Results Report

Report ID: 1	Diagnostic Start Time: 1/30/2025 9:47:25 AM	Diagnostic End Time: 1/30/2025 9:47:25 AM
Customer: UAE	Service Engineer: Kanyakorn S.	
Address: Soi Udomsuk 41, Sukhumvit Rd. Bangkok	Contact Details: 026375363#1	

Configuration:

Serial Number: MY13160001 Turret Type: Automatic
Instrument Model: Varian AA140/240/280 Number Of Lamps: 4
Flame Instrument: True Mono Type: Automatic
Furnace Instrument: True Gasbox Type: 'Y' Gas Box
Zeeman Present: False Auto Burner Adjuster: False
Internal Zeeman: False Mains Frequency: 50
Internal UltraAA: False Firmware Version: 2.11
Optics Type: Double Beam Photomultiplier Type: Normal(900nm)
D2 ES Correction Filter: True PWB Version: 45
Box: Block Version: 1.09

EEPROM Data:

Instrument Run Hours: 60919.180 D2 Run Hours: 63399.600
Zero Wavelength Offset: 0.133 D2 Serial Number: not set
Mono Correction: 0.770 D2 Install Date: 1/1/1970
Flame Hours: 32441.834 D2 Original Intensity: 1.000
D2 Last Intensity: 475.000

Frequency:

Averaging Period: 30.0

Datapoint Count: 20

Upper Limit:
61.00Average Frequency:
50.00Highest Measured Frequency:
50.00Lower Limit:
49.00Lowest Measured Frequency:
50.00Result: **Passed**

Power Supply:

Averaging Period: 30.0



Datapoint Count: 20

	Lower Limit (V)	Actual (V)	Upper Limit (V)	Result:
12.00 V Rail	10.80	12.12	13.20	Passed
-12.00 V Rail	-13.20	-11.96	-10.80	Passed
5.00 V Rail	4.59	5.04	5.50	Passed
310.00 V Rail	270.00	320.00	341.00	Passed

Beam Balance:

Lamp Type: Copper
Lamp Socket Used: 3

Peak Selected: 324.80
Lamp Alignment: **Performed**

Grating: S

Lamp Element(s): Copper
Lamp Socket Position: 3
Lamp Current(mA): 4.00
Lamp Wavelength(nm): 0.5
Lamp Wavelength(nm): 324.80
Lamp Alignment: **Performed**

	Lower Limit (nm)	Upper Limit (nm)	Result:
Zero Current	0.10	0.10	Passed
First Order	0.10	0.10	Passed
Second Order	0.10	0.10	Passed

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Wavelength Report:

Lamp Used: Copper
Peak Used(nm): 324.80
Connected to Socket: 3
Lamp Current(mA): 4
Slit Width(nm): 0.2
Slit Height: Normal

Lower Limit(nm)	Upper Limit(nm)
324.77	324.88
Sample 1: 324.82	Sample 2: 324.823
Sample 3: 324.823	Sample 4: 324.823
Sample 5: 324.823	Sample 6: 324.819
Sample 7: 324.819	Sample 8: 324.819
Sample 9: 324.819	Sample 10: 324.819

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Interlock:

Flame Detect: **Working**
HCO Driver: **Working**
Flame Driver: **Working**
Oxidant Pressure: **Working**
Pressure: **Working**
Ignition: **Working**

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Auto Setup:

Lamp 1: Copper
Lamp 2: Copper
Lamp 3: Copper
Lamp 4: Copper
Lamp 5: Copper
Lamp 6: Copper
Lamp 7: Copper
Lamp 8: Copper
Lamp 9: Copper
Lamp 10: Copper

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
Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Sensitivity Test 01
Comment
Method Cu
Computer name DESKTOP-R9UFRS
Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL ZERO	0.000	38.8	0.0002
Readings			
	0.0002	0.0003	0.0001

STANDARD 1 5.000 0.1 0.5571 1/30/2025 10:51:46 AM

Readings 0.5574 0.5563 0.5575 1/30/2025 10:52:22 AM

Curve Fit = Linear Origin
Characteristic Conc = 0.035 mg/L
r = 1.0000
Calculated Conc = 0.002 ± 0.03
Residuals = -0.002 0.000

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
Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Precision Test
Comment
Method Cu
Computer name DESKTOP-R9UFRS
Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL ZERO	0.000	64.1	-0.0002
Readings			
	-0.0003	-0.0003	-0.0001

STANDARD 1 5.000 0.3 0.6052 1/30/2025 10:46:52 AM

Readings 0.6038 0.6073 0.6047 1/30/2025 10:47:24 AM



Curve Fit = Linear Origin
Characteristic Conc = 0.174 mg/L
r = 1.0000
Calculated Conc = -0.002 ± 0.00
Residuals = 0.002 0.000

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Flame Optimization

Down height 0 (mm)
Pump speed Medium

Key to tube colors
Sample
Calibration
Calibration/QC
Sample/QC
Not Assigned

Sampler Offline

Gold Tube
Back 1
Type 1
Gold Tube

Align Probe
Burn
Stop
Park

Optimization: Lamp

HC Lamp

1.30
1.00
0.50
0.00
0.917

Optimization
Optimization

Rescale
Auto Zero

Gain 49 %

Sensitivity Check 15 mg/L gives about 0.2 Abs at 324.8 nm, A/A burner

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Certificate No.: SP25-021

Page 1 of 5

Customer: United Analyst and Engineering Consultant Co., Ltd. (Head Office)

Address: 3 Soi Udomsuk 41, Sukhumvit Road, Bangchak, Phrakhanong, Bangkok 10260

Location of calibration: Instrument room (207)

Equipment: UV-Vis Spectrophotometer

Manufacturer: HITACHI

Model: U-2900

Serial No.: 21E22-009

ID No.: UAE.WAT.051/2564

Received Date: 26 May 2025

Calibration Date: 26 May 2025

Issue Date: 29 May 2025

Condition Instrument: Good

Calibrated by: (Mr. Tanawat Rittidach)

Approved by: (Ms. Chonticha Sangsorn)

Technical Manager

Quality Manager

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

The measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the DQE Services Co., Ltd.

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REPORT OF CALIBRATION

Certificate No. : SP25-021
Page 2 of 5

Environment Condition : Ambient Temperature 25 ± 5 °C

Relative humidity 55 ± 20 %RH

Calibration method : In-house method CP-01 Based on ASTM E275-08

Certified Reference Materials :

Material	Serial No.	Certificate No.	Due date
Absorbance Standard set	25760	115663	25 October 2025
Absorbance Standard set	25757	115638	25 October 2025
Wavelength Standard set	25806	115657	25 October 2025
Wavelength Standard set	25758	115665	25 October 2025

Traceability : This certification is traceable to the International System of Unit maintained at National -
Institute of Standards and Technology (NIST) through Starna Scientific Limited

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WORK ORDER INFORMATION

Top-Level	Order Type	Preventive Maintenance
Installed Product ID	IB-00105024	Billing Type
Product	SKALAR 2SAN59000 SAN++ Classic CFA 230V 2SAN59000	PO No.
Serial No.	182688	Warranty No.
		Contract No.

PRODUCTS SERVICED

Installed Product Id	Serial Number	Product
IB-00105024	182688	SKALAR 2SAN59000 SAN++ Classic CFA 230V 2SAN59000

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REMARKS

Travel Time Disclaimer:

Please note that the travel time in this report only includes time taken to reach the installed equipment location. It does not include our engineer's return travel time.

Customer Signature:

K Kamphong

Technician: Yongyuth Chanphong
Job Title: Service Manager
Email: yongyuth.yc@dksh.com

Customer Signature

Date: 06/06/2025

Reach us at DKSH Service-Hotline : +66 2 639 7000
2533 Sukhumvit Road, Bangkok, 10260, Phrakhanong, Bangkok, Thailand
Phone +66 2 639 7000 Fax +66 2 333 1028

Delivering Growth - in Asia and beyond

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Job No. WO-00074079

Test Report

Customers	United Analyst and Engineering Consultant Co., Ltd.
Equipment	Continuous Flow Analyzer
Controller Model	SA5000
Controller Serial No.	182688
Date of test	23-May-2025
Environment temperature	24.3 °C
	Humidity
	54.4 %RH

Results

Item	Characteristic	Before	After	Remark
1	Visual inspect	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
2	Power supply (210 - 240 VAC)	220 VAC	220 VAC	
3	Computer	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
4	Program	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
5	Auto sampler	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
6	Module holder			
	- Motor pump	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Pump tube	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Air-injection	☒ Pass ☐ Fail	☒ Pass ☐ Fail	*
	- Chemistry manifolds, Switching valve, Coil, Membrane	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
7	Detector			
	- Filter	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Flow cell	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
	- Lamp	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
8	Interface	☒ Pass ☐ Fail	☒ Pass ☐ Fail	
9	Rinsing valves	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
10	Temperature / Reactor	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
11	Flame photometer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	
12	UPS / Stabilizer	☒ Pass ☐ Fail ☐ N/A	☒ Pass ☐ Fail ☐ N/A	

Warning and Error Checked

Item	Error list	Before	After
13	Error list	☐ None ☐ Appear	☐ None ☐ Appear

DKSH Technology Limited (Head office)
2533 Sukhumvit Road, Bangkok, Phrakhanong, Bangkok, 10260
Phone +66 2 639 7000, Mobile +66 93 813 8681, yongyuth.yc@dksh.com, www.dksh.com

Delivering Growth - in Asia and Beyond

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Standard Equipment Used

Equipment	Equipment I.D.
Digital multi meter	S/N 57600591 Due date: 19-Jun-2025
Thermo hygrometer	S/N 39520444/904 Due date: 27-Dec-2025

Test By:

(Mr. Yongyuth Chanphong)

Approved by:

(Mr. Ewasong Wanklang)

Position:

Supervisor, Technical Service

Position:

Manager, Technical Services


DKSH Technology Limited (Head office)
2533 Sukhumvit Road, Bangkok, Phrakhanong, Bangkok, 10260
Phone +66 2 639 7000, Mobile +66 93 813 8681, yongyuth.yc@dksh.com, www.dksh.com

Delivering Growth - in Asia and Beyond

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

Page 2/2



Base Line Test : Reagent_Baseline_CN_Phenol_NH3



Detector Signal Test : Phenol

Detector Signal Test : NH3

Detector Signal Test : CN

List of Instrument Certificates for Environmental Quality Analysis

No.	Instrument/Equipment	Parameter	Manufacturer	Model/Serial No.	Calibrator	Certification No.	Date of Calibration	Due date of Calibration
1	Atomic Absorption Spectrophotometer	CADMIUM CHROMIUM TRIVALENT LEAD MERCURY NICKEL ZINC	Agilent Technologies	AAS085 / M13160001	Agilent Technologies (Thailand) Co., Ltd.	Preventive Maintenance Checklist	30/7/2025	29/7/2026
2	Atomic Absorption Spectrophotometer	ARSENIC	Perkin Elmer	PinAAcle 900F / PE20001702	Perkin Elmer Co., Ltd.	Preventive Maintenance Report	29/4/2025	28/4/2026
3	Spectrophotometer	CHROMIUM HEXAVALENT	Agilent	Cary 60 G8860A / M151410009	DQS Services Co., Ltd.	SP25-019	26/5/2025	25/5/2026

Due Date of Calibration: Based on the annual calibration plan, at least 1 time per year.

United Analyst and Engineering Consultant Co., Ltd. (UAE)
Certified Laboratory (COPEC) 1023

Certificate Page 1 of 1

Agilent 55 240 280 Series Atomic Absorption Spectroscopy Systems

Preventive Maintenance Checklist

Agilent Preventive Maintenance provides factory recommended service for your analytical systems to assure reliable operation and the accuracy of your results.

Delivered by highly trained and certified service engineers using genuine Agilent parts and supplies, Agilent Preventive Maintenance provides everything you need to reduce unplanned downtime and keep your systems operating at their peak. This checklist will be completed at the end of the service and provided to you as a record of the installation.

Note: While non-current production AA instrument and/or accessory models are not covered specifically in this document it can be used as a basic reference.

For more information about Agilent Technologies services please visit our web site using the following URL: <http://www.agilent.com/en-us/services>

Introduction

Customer Information

- Customers should provide all necessary operating supplies upon request of the engineer.
- A customer representative should be available to the engineer while performing the preventive maintenance procedures.
- Any parts, not included in the Parts Lists section of this document, are not part of the recommended Preventive Maintenance service, nor are they included in the price of this service.
- If a system requires the use of extra or special procedures and/or parts for the maintenance service, then those must be ordered separately and charged as a repair, which may incur additional costs.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Important Customer Web Links

- For more information about Agilent Technologies services, please visit our website using the following URL: <http://www.agilent.com/en-us/products/crosslab-instrument-services/service-repair>
- To access Agilent University, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- A useful Agilent Resource Center web page is available, which includes short videos on maintenance, quick lists of consumables for new instruments, and other valuable information. Check out the Resource Page here: <https://www.agilent.com/en-us/agilentresources>
- Need technical support, FAQs, supplies? – visit our Support Home page at <http://www.agilent.com/search/support>
- Get answers. Share insights. Build connections: Join the Agilent Community at <https://community.agilent.com/welcome>

Service Engineer's Responsibilities

- Contact the customer and ensure that all necessary supplies are available before the preventive maintenance visit.
- Confirm the ability of the instrument to deliver continued safe operation as established via the Agilent AA safe operation flow chart. (Refer directly to the AA 55/240/280 Preventive Maintenance Scope of Work to make this decision.)
- Only select those pages that relate to the system or module being serviced.
- Complete empty fields with the relevant information.
- Complete the relevant checkboxes in the checklist using either a "X" or tick mark "✓".
- Check "Section not applicable" check boxes to indicate services/tasks not delivered, as appropriate.
- Complete the Preventive Maintenance service in the order of the tasks listed.
- Complete the Service Review section together with the customer.
- Complete the fields for page numbers at the foot of each selected page.
- Complete the total number of pages field in the Service Completion section.
- Ask the customer to sign the Service Completion section including the customer's and your signature.

This information is subject to change without notice.

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

Instrument Preventive Maintenance Checklist

Instrument Maintenance

System Information

☐ Check this box if an instrument configuration report is attached instead of completing the table.

Instrument System Name and ID	240 FS AAS
Instrument System Site and Location	United Analyst and Engineering Consultant

List System Component Product Numbers	List the Serial Numbers of each Component
1. G 9432 A	M13160001
2.	
3.	
4.	
5.	
6.	
7.	
8.	
9.	

Preparation, Safe operation and Initial performance checks

Revision: 10.00, Issued November 2021

© Agilent Technologies, Inc. 2021

- ☐ Agilent AA safe operation flow chart inspections (to determine if the PM can be performed).

NOTE: If by following the flow chart the instrument is deemed to be unsafe for continued use you **MUST NOT** continue PM work. Inform the customer immediately of the Agilent recommendation that use of the instrument be discontinued.

- ☒ Discuss any specific issues with the customer before starting.
- ☐ For HF application systems, if standard sample introduction system was not installed, ask the customer to install it. ~~888~~
- ☒ Review the instrument logbook for recorded problems and comments.
- ☒ Save instrument control settings before starting the procedure.
- ☒ Perform a general inspection of the system for cleanliness.
- ☒ Check for proper installation of parts, assemblies, sensors etc
- ☒ Check system for required installation of components, settings as defined by current Service Notes
- ☒ Check for required firmware updates and verify with customers if they would like them installed.
- ☒ Use SVD to perform a Full Wavelength Scan for Cu HCL - "As found test_1"
- ☒ Perform a Basic Cu ABS test - "As found test_2"
- ☒ Print the Details page or screen captures of the test results and attach to the end of this checklist.

Revision 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

Preventive Maintenance Procedures

FLAME SYSTEM section

- ☐ Section not applicable

Electronic components

- ☒ Review and confirm instrument configuration data in SVD
- ☒ Confirm power supply voltages using the *SVD Power Supply diagnostic*.
- ☒ For Dual Beam Instruments - Confirm RBC frequency using the *SVD RBC frequency diagnostic*.

Mechanical components

- ☒ Check the burner adjuster controls for complete and free movement. If the burner adjuster needs lubrication, use Molykote 321 or mineral based molybdenum disulphide grease.
- ☒ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☒ Monochromator drive
 - ☒ Slit drive
 - ☒ Lamp selector
 - ☐ ABA

Optics components

- ☒ Check that external optical surfaces are clean – Clean or replace as required.
- ☒ Use SVD and perform *Mono Wavelength Correction*.
- ☒ Use SVD and perform *Slit Calibration*.
- ☒ Use SVD and perform *Grating Squareness Diagnostic*.
- ☒ Use SVD and perform *Zero Order Offset/Mono Correction*.
- ☒ Use SVD and perform *Wavelength Repeatability*.
- ☒ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.
- ☒ Check that the signal energy of the D2 and HC lamps track properly. Advise customer if their D2 lamp is showing emission degradation due to age.

Revision 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

Sample Introduction and Atomization

- ☒ Inspect the burner interlock plate to ensure that the interlock pin is secure and correct for the burner type.
- ☒ Clean the burner slot with a clean white cord.
- ☒ Check the uniformity of the slot width.
- ☒ Clean the burner if required.
- ☒ Change the burner o-ring.
- ☒ Clean the nebulizer, spray chamber and liquid trap.
- ☒ Change all o-rings and seals in the nebulizer, nebulizer block and spray chamber.
- ☒ Check that the pressure relief bung releases readily.
- ☒ Change o-rings on the fuel and oxidant delivery bars.
- ☒ Leave the liquid trap EMPTY and verify the flame will not ignite in this state.
- ☒ Refill liquid trap and check that overflow drains freely into the drain/waste tube.
- ☒ Check the drain/waste tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel.
- ☒ Check and clean the igniter electrode

Gas handling components and safety interlocks

- ☒ Pressure test for leaks
- ☒ Leak test gasbox internal components and connections
- ☒ Check safety interlock status and operation using the *SVD interlock monitoring diagnostic*.

Analytical performance for Flame systems

- ☒ Ignite a flame.
- ☒ Check that you can adjust the nebulizer uptake rate from 4 to 6.5 mL per minute.
- ☒ Optimize the instrument ready to perform Cu sensitivity test.
- ☒ Create a manual method to perform a Basic Cu ABS test - "Final Performance Testing"
- ☒ Run a PM completed sensitivity test for a 5 ppm copper sample and record the results in the AA PM Performance test results and measurements table.

Revision 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

FURNACE SYSTEM section

- ☐ Section not applicable

Electronic components

- ☐ Review and confirm instrument configuration data in SVD
- ☐ Confirm power supply voltages using the *SVD Power Supply diagnostic*.

Mechanical components

- ☐ Run SVD tests to exercise all motor drives over the full range of their travel:
 - ☐ Monochromator drive
 - ☐ Slit drive
 - ☐ Lamp selector

Optics components

- ☐ Check that external optical surfaces are clean – Clean or replace as required
- ☐ Use SVD and perform *Mono Wavelength Correction*.
- ☐ Use SVD and perform *Slit Calibration*.
- ☐ Use SVD and perform *Grating Squareness Diagnostic*.
- ☐ Use SVD and perform *Zero Order Offset/Mono Correction*.
- ☐ Use SVD and perform *Wavelength Repeatability*.
- ☐ Physically inspect selected HC lamps (customer to supply per their choice) and measure the % Gain for each lamp. Advise customer if lamps are showing emission degradation due to age.

Gas handling, water system and workhead component checks

- ☐ Inspect the GTA workhead gas hoses and connections for leaks.
- ☐ Pressure test for gas leaks
- ☐ If the cooler system is accessible (stand-alone) check for correct operation and coolant/water level – this includes any temperature and pressure settings plus filter cleaning (air flow and water)
- ☐ Inspect the GTA workhead water hoses and connections for leaks.
- ☐ Check all graphite components and replace if necessary.

Revision 10.00, Issued November 2021

© Agilent Technologies Inc. 2021

- ☐ Tube
- ☐ Electrodes
- ☐ Shroud

- ☐ Check and clean the end windows on the workhead.
- ☐ Check safety interlock operation.

Analytical performance for Furnace systems

- ☐ Optimize the instrument ready to perform Cu sensitivity test.
- ☐ Run the sensitivity test for a 25 ppb copper sample and record the results in the results table.

PSD autosampler accessory for Furnace systems

- ☒ Section NOT Applicable
- ☐ Check condition of the PSD capillary – replace if necessary
- ☐ Check condition and operation of PSD syringe – ensure it does not have air locks and bubbles.
- ☐ Change PSD rinse bottle o-ring.
- ☐ Check and clean the rinse vessel.
- ☐ Check the drain tube for good drainage. It should not have tight bends, kinks or loops and the lower end must be above the liquid level in the waste vessel.
- ☐ Ensure that the waste vessel is suitable for use with the furnace system.

Sample introduction pump system (SIPS) accessory

- ☒ Section NOT Applicable
- ☐ Re-torque screws securing the hubs, presser arms and pump rotors.
- ☐ Adjust each roller so that it rotates freely.
- ☐ Wipe clean the pump rotor rollers and pump bands with a dry clean cloth.
- ☐ Ensure that the presser arms and the surfaces near the pump are free from dirt and spills
- ☐ Remove the pump module rear cover and check for the incursion of liquids and any signs of corrosion.
- ☐ Re-torque the nuts that fasten the motor mounting plates to the chassis.
- ☐ Check clips securing the diluents holder and replace if necessary.
- ☐ Disconnect, clean T-piece, and reassemble the tubing using the following steps.

Revision: 10.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

- ☐ Remove the T-piece by disconnecting the pump tubes, the pump bands and all other tubing.
- ☐ Place the T-piece in an ultrasonic bath containing strong detergent 1-5% Decon 30 or similar, for approximately 5-10 minutes.
- ☐ Wash the T-piece under a tap with a strong flow of water.
- ☐ Rinse with distilled water through all of the inlets in the reverse direction to normal sample flow.
- ☐ Reassemble.

Sample preparation system (SPS 4) accessory

- ☒ Section NOT Applicable

The Agilent SPS 4 autosampler is designed to need minimal maintenance.

The following maintenance requirements are suggested to maintain the performance of the autosampler.

- ☐ Cleaning the spill tray, rack location mat, end frames and chassis accessories with a damp soft cloth and diluted mild detergent.
- ☐ Cleaning the autosampler cover panels with domestic window cleaner.
- ☐ Checking the X- axis and Z- axis drive belts for cracks, splits, damaged teeth, excessive fraying, color changes or degradation from fumes..
- ☐ Check the X- axis, Theta- axis and Z- axis FFC cables for cracks, incorrect positioning, damaged edge or damaged connectors.

NOTE: The autosampler requires no extra lubrication throughout its lifetime.
For further details refer to the SPS 4 service manual G8410-90050.

Sample preparation system (SPS 3) accessory

- ☒ Section NOT Applicable
- ☐ Check the x-axis and z-axis timing belts – Replace if there is any cracks, splits or color deterioration and belt tension.
- ☐ Check belt tensions - adjust if required
- ☐ Check the lubrication pad for single x-axis shaft. If pad is dry or customer has observed any vibration or erratic movements of the x-axis carriage, add 1 mL of Dow Corning 200 ̈ Fluid, 200 CS into the well.
- ☐ Check the auto-sampler ability to find tube positions - Calibrate if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Revision: 10.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

Vapor generation accessory VGA (hydride generator)

- ☒ Section NOT Applicable
- ☐ Inspect VGA gas supply hose.
- ☐ Inspect/replace VGA pump tubing
- ☐ Check low gas pressure interlock setting – adjust if required.
- ☐ Check precision orifice gas flow setting – adjust if required.
- ☐ Check gas regulator pressure to 46 psi (325 kPa) – adjust if required.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

UltraAA lamp accessory (external)

- ☒ Section NOT Applicable
- ☐ Check the condition of the power cable.
- ☐ Clean the exterior surfaces of the accessory with soft lint free cloth. This cloth can be dampened with warm water or a mild detergent. Do not use organic solvents or abrasive cleaning agents.

Restore System

- ☐ If you have altered the customer's instrumentation during the course of PM, restore to the original status to allow the customer to conduct their normal activities (e.g., reload the customer's method.)

Guidance

If the PM service is performed prior to a qualification service, then use the qualification procedure as a guide for final instrument set up and checkout.

Revision: 10.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

Signature Page

Service Review

- ☒ Attach available reports/printouts of all tests to this documentation.
- ☒ Record the Preventive Maintenance service activity in the customer's records/logbook.
- ☒ Update/reset instrument maintenance counters as appropriate.
- ☒ Affix the PM sticker to the system or instrument logbook based on the customer's request
- ☒ Complete the Service Engineer Comments section if there are additional comments.
- ☒ Review this service, parts replaced, and test results obtained with the customer.
- ☒ If the instrument firmware was updated, record the details of the change in the Service Engineer's Comments box or if necessary, in the customer's IQ records.

Test Results

Test Description	Expected Test Result	Actual Test Result
Flame optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	49 %
Flame performance test with 5 ppm copper sample		
Air /acetylene, mixing paddle removed	Abs value > 0.5	0.559%
Air /acetylene, mixing paddle installed, 10 replicates	%RSD < 1.0	0.2 %
Deuterium furnace optics PMT Gain test		
For copper at 324.8 nm, 4 mA, 0.5 nm slit width	< 55 %	-
Deuterium furnace performance test with 25 ppb copper sample (324.8 nm)		
Precision %RSD	≤ 4.0 %	-
Abs value	≥ 0.15	-
Zeeman furnace analytical performance: 25 ppb copper sample (327.4 nm)		
Precision %RSD	≤ 4.0 %	-
Abs value	≥ 0.10	-
MSR%	≥ 70 %	-

Revision: 10.00, Issued: November 2021

© Agilent Technologies, Inc. 2021

AA consumable and parts list table

Part Description	Part Number	Product/Model # where used	PM supplied or Consumable	Instrument Type
Test Solution - Cu 5ppm solution	6610030100	50 55 140 240 280	PM supplied	Common
Test Solution - Blank solution	5190-7001	50 55 140 240 280	PM supplied	Common
Copper, 1000 ug/ml, 100ml	5190-8279	50 55 140 240 280	*	Common
Kit, Mix 7 O-rings, aqueous, complete set	9910053400	50 55 140 240 280	PM supplied	Flame
Organic Kit	9910053500	50 55 140 240 280	PM supplied	Flame
Wire Nebulizer Cleaning	9910024700	50 55 140 240 280	consumable	Flame
Tubing Capillary Std Neb	9910024800	50 55 140 240 280	consumable	Flame
Capillary Tube Hvac Neb (3) (organics only)	9910044000	50 55 140 240 280	consumable	Flame
Glass impact beads (5/pkg)	9910025700	50 55 140 240 280	consumable	Flame
Teflon impact beads (5/pkg) (organics only)	9910053300	50 55 140 240 280	consumable	Flame
Burner cleaning strip (100/pkg)	9910053900	50 55 140 240 280	consumable	Flame
Window UV silica - round (right side)	2010026000	50 55 140 240 280	PM supplied	Common
Window UV silica - rectangular (left side)	2010026500	50 55 140 240 280	PM supplied	Common
Pad adhesive window (round)	4910012700	50 55 140 240 280	PM supplied	Common
Pad adhesive window (rectangular)	4910012800	50 55 140 240 280	PM supplied	Common
Electrode kit (1 pr) (D2)	6310003400	GTA120	PM supplied	Furnace
Shroud (D2)	6310003100	GTA120	PM supplied	Furnace
Zeeman electrode kit (1 pr)	6310003500	GTA120	PM supplied	Furnace
Zeeman shroud	6310003600	GTA120	PM supplied	Furnace
O-ring PGD rinse bottle	6910025900	PSD120	PM supplied	Furnace

* For engineers who only service AA instruments 5190-8279 can be used as a cheaper alternative for 6610030100.

Items classified as PM supplied in the above table are included in the standard PM. Those classified as consumable should be provided by the customer or charged to the customer if supplied by the Agilent service engineer.

Revision: 13.00, Issued: November 2021
© Agilent Technologies, Inc. 2021

Service Engineer Comments (optional)

If there are any comments, please fill in this box. If there are no comments, please leave it blank. This box is for the service engineer's use only.

Service Completion

Service request number: 6003549143 Date service completed: 30 Jan 2025
Agilent signature: Kanyakorn S. Customer signature: Sunda Y.
Total number of pages in this document: 13

Revision: 13.00, Issued: November 2021
© Agilent Technologies, Inc. 2021

SVD Results Report

Report ID: 0100056 Start Time: 1/30/2025 9:14:28 AM End Time: 1/30/2025 12:14:28 PM
Customer: UAE Service Engineer: Kanyakorn S.
Address: Spa Udomsuk 41, Sukhumvit Rd, Bangkok Contact Details: 026378353#1

Configuration:

Serial Number: MY13160001 Turret Type: Automatic
Instrument Model: Varian AA140/240/280 Number Of Lamps: 4
Flame Instrument: True Mono Type: Automatic
Furnace Instrument: True Gasbox Type: 'Y' Gas Box
Zeeman Present: False Auto Burner Adjuster: False
Internal Zeeman: False Mains Frequency: 50
Internal UltraAA: False Firmware Version: 2.11
Optics Type: Quartz Beam Photomultiplier Type: Normal(900nm)
D2 CO Correction Filtered: True PWB Version: 45
Boot Clock Version: 1.00

EEPROM Data:

Instrument Run Time: 67718.180 D2 Run Hours: 53398.500
Zero Wavelength Offset: 0.133 D2 Serial Number: not set
Mono Correction: 0.770 D2 Install Date: 1/1/1970
Flame Hours: 32411.834 D2 Original Intensity: 1.000
D2 Last Intensity: 475.000

Frequency:

Averaging Period: 30.0
Datapoint Count: 20

Upper Limit: 51.00 Highest Measured Frequency: 50.00
Average Frequency: 50.00
Lower Limit: 49.00 Lowest Measured Frequency: 50.00

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
Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Sensitivity Test 01
Comment
Methods Cu
Computer name DESKTOP-RSUFPS
Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL ZIRCO	0.000	38.8	0.0002
Readings			
0.0002	0.0000	0.0001	1/30/2025 10:51:45 AM

STANDARD 1	5.000	0.1	0.5571
Readings			
0.5571	0.5563	0.5575	1/30/2025 10:52:22 AM

Curve Fit = Linear Origin
Characteristic Conc = 0.030 mg/L
r = 1.0000
Calculated Conc = 0.002 mg/L
Residuals = -0.0007 0.0000

5 ppm Cu	Abs = 0.11141 x C	5.000	0.3	0.6558
Readings				
0.6552	0.5595	0.5515	1/30/2025 10:52:54 AM	

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
Analyst
Date Started 1/30/2025 10:33 AM GMT: 1/30/2025 3:33 AM
Worksheet Precision Test
Comment
Methods Cu
Computer name DESKTOP-RSUFPS
Serial Number: MY13160001

Method: Cu (Flame)

Sample ID	Conc. mg/L	%RSD	Mean Abs
CAL ZIRCO	0.000	64.1	-0.0002
Readings			
-0.0003	-0.0003	-0.0001	1/30/2025 10:48:52 AM

STANDARD 1	5.000	0.3	0.6052
Readings			
0.6030	0.6073	0.6047	1/30/2025 10:47:24 AM

Curve Fit = Linear Origin
Characteristic Conc = 0.1705 mg/L
r = 1.0000
Calculated Conc = -0.0025 0.0000
Residuals = 0.002 0.0000

5 ppm Cu	Abs = 0.12105 x C	5.000	0.7	0.6261
Readings				
0.6035	0.5952	0.6047	0.6042	0.6079
0.6035	0.5776	0.6064	0.6079	1/30/2025 10:48:32 AM

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PinAAcle 900F Preventive Maintenance (PM)

Company Name:	United Analyst and Engineering Consultant Co., LTD.		
Address (Instrument Location):	41 Sukumvit Rd., Phra Khanong, Bangkok 10260		
Serial Number:	PFBS20031902	PM Number:	1 of 2
Customer Name (if applicable):	K. Yinda	Telephone Number:	095-5580049
Customer Support Engineer Name:	K. Chayanon	Service Order Number:	WO-03126047
Date PM Performed: (DD-MMM-YYYY)	29-Apr-2025	Next PM Due Date: (DD-MMM-YYYY)	29-Oct-2025
Standard Labor Hours to Complete PM :		5 hours	

Part Number	Release	Publication Date	
09370145 Rev.9	A	January 2018	

Scope

The purpose of this PM is to ensure the continued functionality of the PinAAcle 900F by inspecting and replacing any worn or damaged parts. This service should only be performed by a trained representative of PerkinElmer.

The customer should save their method before the PM begins.

General Instructions:

The customer must provide the engineer operational data to demonstrate recent instrument performance prior to starting the PM.

Always check with the customer before making any changes that may affect the customer's analysis or calibration, including a current back-up of system software and/or data files.

The completed document should be signed by an authorized PerkinElmer and customer representative and left with the customer.

Update the PM sticker and instrument logbook as required.

Copyright Information

This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc. Copyright © 2013 PerkinElmer, Inc.

Trademarks

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law. PerkinElmer is a registered trademark of PerkinElmer, Inc. All other trademarks and registered trademarks not owned by PerkinElmer, Inc. or its subsidiaries that are depicted herein are the property of their respective owners. Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no Warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. PerkinElmer shall not be liable for incidental or consequential damages in connection with the furnishing or use of this document.

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Parts Lists

Parts Included with the PM		
Part Number (if applicable)	Description	Quantity
B0501696	Fan Filters	2
N3160156	O-Ring Kits for Sampling Introduction (Stainless Steels Nebulizer)	N/A
N3160157	O-Ring Kits for Sampling Introduction (Plastic Nebulizer)	1
N9301714	Replacement Acetylene Filter Cartridge	1
TH001022	Replacement Air Filter Cartridge	1

Additional Reagents and Standards Required for PM

Part Number (if applicable)	Description	Quantity	Batch/Lot #	Expired Date (MM/YY)
N9300183	1000 mg/L Copper Standard	AR	27-39CUY1	Nov 2025

Additional Reagents and Standards Required for PM (Customer Support Solution)

Part Number (if applicable)	Description	Quantity	Batch/Lot #	Expiration Date (MM/YY)
N/A	DI Water	250 mL	AR	AR
N/A	0.5% HNO ₃	250 mL	AR	AR

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Procedure Checklist

Use (✓) to check off those steps in the checklist that have been completed.

1. General:

- ☒ Review the instrument performance with the customer and document any recent problems.
- ☒ Inspect the customer log book and make any appropriate PM entries.
- ☒ Perform general inspection of system for cleanliness.

2. PC Instrument Software:

- ☒ Instrument Software user files/databases archived, packed, and/or deleted as needed.

3. Mechanical:

- ☒ Inspect and clean all fans and filters. Replace filters if necessary.
- ☒ Inspect all gas lines for leaks and/or wear. Replace if needed.
- ☒ Clean exterior of the instrument.
- ☒ Inspect the burner head, burner chamber, and nebulizer. Clean if needed as stated in the Hardware Guide.
- ☒ Check burner head dimensions with the feeler gauge as stated in the Hardware Guide in the Maintenance chapter section on cleaning the burner head and checking slot width. Replace if out of specification.
- ☒ Check the condition of the end cap, burner head, and nebulizer O-rings. Replace if necessary.
- ☒ Check the drain system for signs of wear. Replace worn or damaged parts.
- ☒ Visually check for proper flame conditions when lighting the Air-C2H2 and N2O-C2H2 flames (if applicable).

4. Electrical:

- ☒ Inspect PC boards. Clean if necessary.
- ☒ Carefully check all internal and external cable connections.
- ☒ Check instrument firmware revisions upgrade to current levels (if necessary).
- ☒ Run Diagnostics Test within the Advanced function of the Spectrometer page. Check the results in the service log folder in the Spectrometer BM Log Viewer.

5. Optics:

- ☒ Inspect and clean the sample compartment windows, if needed.
- ☒ Inspect optics. Clean or replace if necessary.

6. Gasses:

- ☒ Verify that the Gasses supplied to the instrument are within the pressure and purity specifications found in the PinAAcle 900 Series Pre-Installation Checklist SDB.
- ☒ Verify that the acetylene filter and air filter element is dry. Replace if necessary.

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7. Flame Interlock Check:

Description: Check to ensure that all safety interlocks are closed.

Parameter	Specification	Test Results	Pass/Fail
Flame Sensor	Air/C ₂ H ₄ Flame correctly shuts down	Active	Passed
Drain Sensor	Air/C ₂ H ₄ Flame correctly shuts down	Active	Passed
Nebulizer Sensor	Air/C ₂ H ₄ Flame correctly shuts down	Active	Passed
C ₂ H ₄ Pressure Sensor	Air/C ₂ H ₄ Flame correctly shuts down	Active	Passed
Air Pressure Sensor	Air/C ₂ H ₄ Flame correctly shuts down	Active	Passed
Burner Head Sensor	Choosing Nitrous Oxide as the oxidant should trigger an interlock shut down	Active	Passed

8. After PM Performance tests:

8.1 Detector Linearity with Barium

Description: Ensures that the detector is linear in the Visible Range.

Parameter	Specification	Certificate Value at 553.6 nm (Abs.)	Test Results	Pass/Fail
1.0 A ND Filter	± 5% from Cert.	0.9668	0.9878	Passed
0.2 A ND Filter	± 5% from Cert.	0.1953	0.1876	Passed

8.2 Baseline Noise at 1.0 Absorbance with Barium

Description: Ensures that a high absorbance will not produce excessive noise.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.010	0.0005	Passed

8.3 AA Baseline Noise with Copper

Description: Check baseline noise.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.001	0.0001	Passed

8.4 D₂ Background Compensation with Copper

Description: Verifies the instruments ability to compensate for Background absorption.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.010	0.009	Passed

8.5 AA-BG Baseline Noise with Copper

Description: Ensures that background correction does not produce excessive noise.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.005	0.0001	Passed

8.6 AA-BG Baseline Noise with Arsenic

Description: Ensures that background correction does not produce excessive noise at a low wavelength.

Parameter	Specification	Results	Pass/Fail
Standard Deviation	≤ 0.005	0.0004	Passed

8.7 Flame Sensitivity

Description: Instrument Sensitivity checked against Copper standard.

Standard Copper Sensitivity	Specification	Results (Abs.)	Pass/Fail
5 mg/L Sensitivity SS Neb (if applicable)	> 0.250 Abs.	N/A	Not Applicable
2 mg/L Sensitivity HS Neb (if applicable)	> 0.250 Abs.	N/A	Passed

10. Review:

- ☒ Review with the customer PM work performed.
- ☒ Review with the customer routine maintenance procedures.
- ☒ Discuss recommended customer supplied materials to have on hand.
- ☒ Attach PM sticker.

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Additional Comments

Additional Comments Regarding the PM

Review

The preventive maintenance checks and if applicable performance tests for PinAAcLe 900F have been completed.

This PinAAcLe 900F ☒ Passes ☐ Fails the preventive maintenance.

Review of Preventive Maintenance:

Authorized PerkinElmer Representative:	Date:
<i>Chaynam k</i>	29 Apr 2025
Authorized Customer Representative:	Date:
<i>Chaynam k</i>	29 Apr 2025

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Work Description

- PM 22 (AA-BG)
- Cleaning Cell Chamber, Filter
- Wavelength Calibration / Pass
- Wavelength Scan AA Cu, Ba, K, Ni ; Pass
- Align cell with Hg ; OK

Start Date	End Date	Work Description
29/04/2025	09/05/2025	

Quantity	Calibrated Tool	Description	Serial Number	Last Calibration Date	Next Calibration Date

Part Number	Part Description	Note	Lot/Serial Number	Quantity

Part Number	Part Description	Start Date	Quantity
SV000013	Preventive maintenance	29/04/2025	4

Work Complete	Customer Signature	Technician Signature
Yes ☒ No ☐	<i>Chaynam k</i>	<i>Chaynam k</i>
PM/OP/PPV Left with Customer		
Yes ☐ No ☒		

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Customer Acknowledgment of receipt of the above repair / replacement.

Special Terms and Conditions: This is not an invoice.

Terms will be applied to your invoice if applicable.

Atomic Absorption/FIAS 100/400 Preventive Maintenance (PM)			
Company Name:	United Analyst and Engineering Consultant Co., LTD.		
Address (Instrument Location):	41 Sukumvit Rd., Phra Khanong, Bangkok 10260		
Room Number:	Lab		
Asset Number (if applicable):	2 of 2W	Customer System ID:	K. Yenida
Service Engineer Name:	K. Chayanan	Service Order Number:	WO-03051971
Date PM Performed: (DD-MM-YYYY)	29-Apr-2025	Next PM Due Date: (DD-MM-YYYY)	29-Oct-2025

Part Number	Release	Publication Date	
09370005	C	January 2013	

Scope

The purpose of this PM is to ensure the continued functionality of the Atomic Absorption/FIAS 100/400 by inspecting and replacing any worn or damaged parts. This service should only be performed by a trained representative of PerkinElmer. The customer should save their method before the PM begins.

General Instructions:

Always check with the customer before making any changes that may affect the customer's analysis or calibration. The completed document should be signed by an authorized PerkinElmer and customer representative and left with the customer. Update the PM sticker and instrument logbook as required.

Copyright Information

This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this publication may be reproduced in any form whatsoever or translated into any language without the prior, written permission of PerkinElmer, Inc. Copyright © 2013 PerkinElmer, Inc.

Trademarks

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law. PerkinElmer is a registered trademark of PerkinElmer, Inc. All other trademarks and registered trademarks not owned by PerkinElmer, Inc. or its subsidiaries that are depicted herein are the property of their respective owners. Except as specifically set forth in its terms and conditions of sale, PerkinElmer makes no warranty of any kind with regard to this document, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. PerkinElmer shall not be liable for incidental or consequential damages in connection with the furnishing or use of this document.

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Procedure Checklist

Use (✓) to check off those steps in the checklist that have been completed.

1. General:

☒ Review the instrument performance with the customer and document any recent problems.

☒ Is the Working Environment Acceptable? If not, document.

☐ Visual Damage (if yes, describe)

☒ Check incoming AC line voltage for proper levels and grounding.
☒ Verify Voltage switch on back of instrument is correct.
☒ Perform general inspection of system for cleanliness. Clean if needed.
☒ Gas supply cylinders secured, lines leak checked and argon or nitrogen supply pressure verified (45 – 58 psi).
☒ Inspect the customer log book and make any appropriate PM entries.
☒ Fan checked and filter cleaned
☒ Heating mantle or Universal Cell Holder checked

2. Instrument components

☒ Non-return valve checked/repaired/replaced if needed (B019 B111). Clean the valve if there is any liquid in it. Replace the rubber sleeve (B013 5123) if it is worn. Check the flow meter for any signs of fluid in it. Clean the flow meter if needed.
☒ Verify condition of pump pressure adjustment levers (B050 7794 - look for cracks or problems with the springs), pump rollers (B300 0251 check for wear), and thumb screws (B050 7796).
☒ Check the Multipot valve for proper switching, flow, and insure there are no leaks. Clean valve parts and replace o-rings if needed (large o-ring: B050 1250, small o-ring: B004 5095). Use a squirt bottle & fishing line to try to dislodge clogs.
☒ Firmware Version checked. Latest is 2.20.

3. Mixing/Separation Assembly & Pump Tubing:

☒ Mixing separator assembly checked
☒ Filters/membranes checked (B050 8306)
☒ Condition of the pump tubing (replace if necessary), correct pump tubing for the solutions being run. Make sure the correct magazines are being used. B050 7791 for 0.13 – 1.80 mm tubing; B050 7792 for 1.60 – 3.18 mm tubing.

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4. Cell, Cell Windows, Transfer Line:

- ☒ Cell checked
- ☒ Cell windows checked
- ☒ Transfer line checked for moisture (if moisture is a problem, the Nafion dryer might be needed)

5. Operational Tests:

- ☒ Run DI water through the carrier/reductant/sample system. Verify smooth flow of liquid throughout without leaks. Replace tubing & fittings if needed.

6. Review:

- ☒ Review with the customer PM work performed.
- ☒ Review with the customer routine maintenance procedures.
- ☒ Discuss recommended customer-supplied materials to have on hand.
- ☒ Attach PM sticker.
- ☒ Update Logbook.

Additional Comments

Additional Comments Regarding the PM

Review

The preventive maintenance checks and if applicable performance tests for FIAS 100/400 have been completed.

This FIAS 100/400 Passes ☒ Fails ☐ the preventive maintenance.

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Document History

Revision	Description of Change	Page(s)	Date
A	First release		May 2008
B	Addition of Batch/Lot Number, Expiration Date, and Report Fields.	2,7	February 2009
C	Update to new format	All	January 2013

PerkinElmer Scientific (Thailand) Co., Ltd.
299 Sol Bounplai 4
Kwang Bangkapi, Khut Huay Kwang
Bangkok 10310
Thailand
Tel: 66 2719 8420 ; Fax: +66 2 319 7900
http://www.perkinelmer.com

Service Report

Work Order Number	Activity Code	Billing Type	Requested Start Date	Model	Serial Number
W0203051271	Planned Maintenance	Contract	10/03/2563 23.08 h.	80505370	10023400521
Service Representative Name	Contract Number	Expiry Date	Equipment ID	System ID	
Kanan, Chayman	SC-0033659090	24/03/2025	N/A	N/A	
Customer Location					
N/A					
Equipment Location					
N/A					
Customer Contact					
K. Nithakorn surit	Phone Number	085-6580048	Fax Number	N/A	Email
					perkin@k12@gmail.com
					Purchase Order
					HPD-240400211

Work Description

- PM 2/2 Warranty
- Cleaning Port Vals, Mandrol, Taring
- Run (up) test, Pass

Start Date	End Date
29/04/2025	29/04/2025
29/04/2025	29/04/2025

Work Description

--	--

Tools Used

Quantity	Calibrated Tool	Description	Serial Number	Last Calibration Date	Next Calibration Date
*** No Calibrated Tools Used ***					

Material Used

Part Number	Part Description	None	Lot/Serial Number	Quantity
*** No Parts Used ***				

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Remark : - UUC = Unit Under Calibration

- N/A = Not Available

- The result expanded uncertainty of measurement U is stated as the standard uncertainty of measurement multiplied by the coverage factor k, which for a normal distribution corresponds to a coverage probability of approximately 95%

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TM-708-02 R01 1/11/2021

TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
134-4 PATTANAKARN ROAD, SUKHUMVIT 21, SUKHUMVIT, BANGKOK 10250
TEL: 0-2717-9430/4-2414 FAX: 0-2716-9444

Certificate of Calibration

Certificate No.: 25H806
Page: 1 of 2

Equipment: Dial Thermo-Hygrometer

Manufacturer: Bangjo

Model: -

Serial No.: -

ID No.: UAE-ANV-094/2548

Condition As-Received: Used Item

Received Date: 04 April 2025

Calibration Date: 10 April 2025

Reference: 2504-0193WSC

Ambient Temperature: $(25 \pm 3) ^\circ\text{C}$

Relative Humidity: $(50 \pm 20) \%$

This certificate may not be reproduced other than in full,
except with the prior written approval of the head of
Corporate Services 3: Equipment Calibration and Testing Services.

Submitted by: United Analyst and Engineering Consultant Co., Ltd.

81 Soi Udomsak 41, Sukhumvit Road, Bangkok,
PhraNakhon, Bangkok 10260

Procedure used: Calibration were conducted using in-house calibration procedure CP-H02 according to comparison
with standard chilled mirror sensor for humidity measurement function and comparison with standard
temperature probe for temperature measurement function into humidity / temperature chamber.

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1. Reference standards Instruments

Instrument	Model	Serial No.	Certificate No.	Due Date
1) Chilled Mirror Hygrometer	Dew Master	44730	22685	10 Sep 2025
2) Handheld Thermometer With Sensor	1521	ASA339	2411176	25 Oct 2025

2. The certificate is valid only to the item calibrated on date and place of calibration.

3. This Certification is traceable to the International System of Unit maintained through:

- Thunder Scientific Corporation, NVLAP Accreditation No. Calibration 200582-9
- Technology Promotion Association (Thailand-Japan), NSC-ONSC Accredited No. Calibration 0008

Calibrated by: Kriengkong Chrit
Issue Date: 10 April 2025

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Cert. No.: 25H806
Page: 2 of 2

Result of Calibration: Without Adjustment
Function: Humidity Measurement

Reference Temperature ($^{\circ}\text{C}$)	Standard Humidity (%R.H.)	UUC* Reading (%R.H.)	Correction (%R.H.)	Uncertainty of Measurement ($\pm\%$ R.H.)
25.0	40.1	40	0.1	1.7
25.0	60.0	60	0.0	1.8
25.0	80.0	76	4.0	1.9

Result of Calibration: Without Adjustment
Function: Temperature Measurement

Standard Temperature ($^{\circ}\text{C}$)	UUC* Reading ($^{\circ}\text{C}$)	Correction ($^{\circ}\text{C}$)	Uncertainty of Measurement ($\pm^{\circ}\text{C}$)
20.012	20.5	-0.488	0.72
25.034	25.0	0.034	0.72
30.032	30.0	0.032	0.72
35.022	34.5	0.522	0.72
40.040	39.0	1.040	0.72

UUC*: Unit Under Calibration

The reported uncertainty of measurement was based on standard uncertainty multiplied
by coverage factor $k = 2.00$, providing confidence level approximately 95%

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
134-4 PATTANAKARN ROAD, SUKHUMVIT 21, SUKHUMVIT, BANGKOK 10250
TEL: 0-2717-9430/4-2414 FAX: 0-2716-9444

Certificate of Calibration

Certificate No.: 25P1384
Page: 1 of 2

Equipment: Aneroid Barometer

Manufacturer: Bangjo

Model: -

Serial No.: -

ID No.: UAE-EMA2-109/2055

Condition As-Received: Used Item

Received Date: 04 April 2025

Calibration Date: 17 April 2025

Reference: 2504-0196WBC

Ambient Temperature: $(29 \pm 2) ^\circ\text{C}$

Relative Humidity: $(50 \pm 15) \%$

Atmospheric Pressure: 1000 mbar

This certificate may not be reproduced other than in full,
except with the prior written approval of the head of
Corporate Services 3: Equipment Calibration and Testing Services.

Submitted by: United Analyst and Engineering Consultant Co., Ltd.

81 Soi Udomsak 41, Sukhumvit Road, Bangkok,
PhraNakhon, Bangkok 10260

Procedure used: The calibration was conducted by direct comparison method against Pressure Measuring Instruments
Standard according to calibration procedure CP-P10, using "DKD-R G-1", Calibration of Pressure Gauge " as
a guidelines.

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1. Reference standards Instruments

Instrument	Model	Serial No.	Certificate No.	Due Date
1) Standard Barometer	DP142	1422605046	MP-0153-24	15 May 2025

2. This instrument was installed in vertical orientation and center of the dial was used as the reference level

3. This result of calibration was made on requested at the point specified by customer

4. Scale and conversion factor is 1 kPa = 7.50062 mmHg

5. This result of calibration instrument was in absolute pressure

6. This instrument was used clean air as pressure media

7. The certificate is valid only to the item calibrated on date and place of calibration.

8. This Certification is traceable to the International System of Unit maintained through:

- National Institute of Metrology Thailand (NIMT)

Calibrated by: Kriengkong Saivichai
Issue Date: 21 April 2025

Approved Signatory:

Atthaporn
[] Phatinee Prapapal
[] Surp Sunomman
[x] Atthaporn Paturach

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Result of calibration: Without adjustment
Function: Absolute Pressure Measurement

Range: 770 mmHg to 770 mmHg
Scale Interval: 1 mmHg (The Fifth Estimate)

Increasing Pressure	Applied Pressure (mmHg)	716.76	720.83	742.09	750.76	761.21	772.15
UUC* Indication (mmHg)	720.0	739.0	740.6	760.0	760.0	770.0	
Error (mmHg)	0.22	0.17	0.00	-0.23	1.21	-2.15	

Decreasing Pressure	Applied Pressure (mmHg)	772.15	761.23	750.91	740.12	720.82	710.63
UUC* Indication (mmHg)	770.0	760.0	750.0	740.0	730.0	720.0	
Error (mmHg)	-2.15	-1.23	-0.91	-0.12	0.18	0.17	

The uncertainty of measurement was ± 0.24 mmHg

* UUC = Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied
by a coverage factor $k = 2$ providing a level of confidence of approximately 95 %

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4353 Sukhumvit, Bangna, Bangkok 10260 Tel. 081-454-2804, 0-2399-0469

Calibration Certificate

Issued by Calibration & Test Section Meteorological Instruments Bureau

Date of Issue 19 February, 2025

Certification No. 106/25

Page 1 of 5

Object Wind Speed & Wind Direction Data Logger

Manufacturer SCARLETT/LECH

Type WL-21

Mfg Code Wireless Receiver 2205DR0114

Wind Sensor 2205DT0114

Customer United Analyst and Engineering Consultant Co., Ltd.

81 Soi Udomsak 41, Sukhumvit Road,

Bangchak, Prakanong Bangkok 10260

Calibration Condition Temperature 25.1 °C Barometric Pressure 1011.3 hPa

NATIONAL STANDARD WIND TUNNEL Wind Aloft Plotting Board

Micromanometer Theodor Friedrich, F6014 Serial No. 9310119 HOOK GAGE NO. 1425

N.I.S.T. Test Reference Number 731/241460 Standard Velocity at 20 - 30 m/sec

Ultrasonic Anemometer Model DA 659 3TV (Sensor TR 90AH)

Serial Number 110730029 (Sensor 120629586)

JAPAN QUALITY ASSURANCE ORGANIZATION

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Vane Angel Bench Stand Model 16112	
Young Meteorological Instruments	
WIND DIRECTION	TESTED WIND DIRECTION
0	0
90	90
180	180
270	

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
19 February, 2025

Page 5 of 5

Standard Temp °C	Temperature Sensor Reading	
	Reading °C	Correction °C
40.4	40	0.4
39.1	39	0.4
14.6	16	-0.4

Calibrated by:
Mr. Watcharapol Subwat
Mechanical Engineer

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4353 Sukhumvit, Bangna, Bangkok 10260 Tel. 081-454-2804,0-2399-0469

Calibration Certificate

Issued by Calibration & Test Section Meteorological Instruments Bureau

Date of Issue 17 March, 2025

Certification No. 159/25

Page 1 of 2


Object Wind Speed & Wind Direction
Manufacturer Met One Instruments
Type Sensor 054B ID No UAE LFM 1992551
Data Logger 539
Mig Code Sensor Y11374
Data Logger X20032
Customer United Analyst and Engineering Consultant Co., Ltd.
61 Soi Udomsak 41 Sukhumvit Road,
Bangchak, Prakanong, Bangkok 10260

Calibration Condition Temperature 25.1 °C Barometric Pressure 1009.6 hPa

NATIONAL STANDARD WIND TUNNEL Vane Angel Bench Stand Model 18112
Micromanometer Tesdor Friedrich FC914 Serial No. 9310119
HOOK GAGE NO 1425 Hand Hook Tesdor Friedrich Type 9320 0960 Serial 0925
NIST Test Reference Number 731/241460 Standard Velocity 0.20 - 30 m/sec
Ultrasonic Anemometer Model DA-650-3TV (Sensor TR-00AH)
Serial Number 110732525 (Sensor 120629589)

JAPAN QUALITY ASSURANCE ORGANIZATION Standard Velocity 0.20 m/sec

Calibrated by:
Mr. Watcharapol Subwat
Mechanical Engineer

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
Vane Angel Bench Stand Model 18112 Young Meteorological Instruments	
WIND DIRECTION	TESTED WIND DIRECTION
0	0
90	91
180	180
270	271

Calibrated by:
Mr. Watcharapol Subwat
Mechanical Engineer

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4353 Sukhumvit, Bangna, Bangkok 10260 Tel. 081-454-2804,0-2399-0469

Calibration Certificate

Issued by Calibration & Test Section Meteorological Instruments Bureau

Date of Issue 18 March, 2025

Certification No. 162/25

Page 1 of 5

Object Wind Speed & Wind Direction Data Logger
Manufacturer SCARLETT LICH
Type WL 21
Mig Code Wireless Receiver 2111DR0041
Wind Sensor 2111DT0041
Customer United Analyst and Engineering Consultant Co., Ltd.
61 Soi Udomsak 41 Sukhumvit Road,
Bangchak, Prakanong, Bangkok 10260

Calibration Condition Temperature 25.1 °C Barometric Pressure 1011.4 hPa


NATIONAL STANDARD WIND TUNNEL Vane Angel Bench Stand Model 18112
Micromanometer Tesdor Friedrich FC914 Serial No. 9310119 HOOK GAGE NO 1425
NIST Test Reference Number 731/241460 Standard Velocity 0.20 - 30 m/sec
Ultrasonic Anemometer Model DA-650-3TV (Sensor TR-00AH)
Serial Number 110732525 (Sensor 120629589)

JAPAN QUALITY ASSURANCE ORGANIZATION Standard Velocity 0.20 m/sec

STANDARD THERMOMETER Tesdor Friedrich Day No 8590/94 Wet No 0383/94
Serial No 8590/94 Thermochronix No 915802

STANDARD BAROMETER Digital Barometer Vacuola Type 81000 No. V1220015
Digital Barometer Vacuola Type 81000 No. V1220015

Calibrated by:
Mr. Watcharapol Subwat
Mechanical Engineer

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Vane Angel Bench Stand Model 18112	
Young Meteorological Instruments	
WIND DIRECTION	TESTED WIND DIRECTION
0	0
90	90
180	180
270	180

Calibrated by

Mr. Watcharapol Subwat
Mechanical Engineer

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18 March, 2025

Page 5 of 5

Standard	Temperature Sensor Reading	
	Reading	Correction
Temp	°C	°C
45.3	45	0.3
30.4	30	0.4
15.4	15	0.4

Calibrated by

Mr. Watcharapol Subwat
Mechanical Engineer

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4353 Sukhumvit, Bangna, Bangkok 10260 Tel. 081-454-2804, 0-2399-0469

Calibration Certificate

Issued by : Calibration & Test Section : Meteorological Instruments Bureau

Date of Issue : 28 April, 2025

Certification No. 302/25

Page : 1 of 5

Object : Wind Speed & Wind Direction Data Logger

Manufacturer : SCARLET/TECH

Type : WL-21

Mfg Code : Wireless Receiver 2301DR0024

Wind Sensor 2301DT0024

Customer : United Analyst and Engineering Consultant Co., Ltd.

81 Soi Udomsuk 41 Sukhumvit Road,

Bangchak, Prakanong, Bangkok 10260

Calibration Condition : Temperature 25.1 °C Barometric Pressure 1006.7 hPa

NATIONAL STANDARD WIND TUNNEL : Vane Angel Bench Stand Model 18112

: Micromanometer Theodor Friedrich FC014 Serial No. 9319119 : HOOK GAGE NO. 1425

N.I.S.T. Test Reference Number 731241460 Standard Velocity of 20 - 30 m/sec

: Ultrasonic Anemometer Model DA-650-3TV (sensor TR-90AH)

Serial Number 110730029 (sensor 120620686)

JAPAN QUALITY ASSURANCE ORGANIZATION



: Standard Velocity of 0 - 20 m/sec

STANDARD THERMOMETER : Theodor Friedrich Dry No. 839094 Wet No. 838954

Iusto, Iusto L45 Serial No. 02846567 : Thermochocid No. 918802

STANDARD BAROMETER : Digital Barometer Vaisala Type PTB220 No. V1220015

Digital Barometer Vaisala Type PTB220 No. V1220015

Calibrated by :  Signed : 

Mr. Watchapol Subwat

Mr. Watchapol Subwat

Mechanical Engineer

THAI METEOROLOGICAL DEPARTMENT

4353 Sukhumvit, Bangna, Bangkok 10260 Tel. 081-454-2804, 0-2399-0469

The Result of Calibration


Certification No. 302/25

28 April, 2025

Page : 2 of 5

Standard Ultrasonic Anemometer	HOOK GAGE NO. 1425			TESTED ANEMOMETER	
	Pressure mbar	Vacuum mbar	Velocity m/sec	Velocity m/sec	Correction m/sec
1.00	-	-	-	1.0	0.00
3.02	-	-	-	3.0	0.02
5.00	-	-	-	5.0	0.00
7.04	-	-	-	7.0	0.04
9.02	-	-	-	9.0	0.02
11.02	-	-	-	11.0	0.02
13.01	-	-	-	13.0	0.01
15.01	-	-	-	15.0	0.01
17.02	-	-	-	17.0	0.02
20.02	-	-	-	20.0	0.02

Vane Angel Bench Stand Model 18112 Young Meteorological Instruments	
WIND DIRECTION	TESTED WIND DIRECTION
0	0
90	90
180	180
270	

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

Calibrated by: Hobkrapel
Mr. Watcharapol Subwai
Mechanical Engineer

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3 : EQUIPMENT CALIBRATION AND TESTING SERVICES
53/44 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
TEL. 0-2717-3000 FAX. 0-2719-9484

Certificate of Testing

Cert.No.: 25TW108
Page.: 1 of 2

Equipment : DO Meter
Manufacturer : Horiba
Model : LAQUA-DO210
Serial No. : HE2L0036
ID No. : UAE.EFM.018/2566(EFM.DO.03/66)
Received Date : 20 May 2025
Test Date : 21 May 2025
Reference : 2505-0605WSC-7
Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
3 Soi Udomsuk 41, Sukhumvit Road, Bangkok,
Phrakhanong, Bangkok 10260
Laboratory Condition : Temperature (25 ± 5) °C
Humidity (50 ± 20) %
Test Procedure : In - house method : CP-CH9
by Comparison Technique with Azide Modification Method
Tested by : Walalak Sirihean
Approved by : Sailthip
Approved Signatory
() Chakrit Waewwanjua
() Ponpan Palpim
(✓) Sailthip Meangmal
Issue Date : 23 May 2025

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1. Reference Standard Instruments :

This measurement result is traceable to the International System of Unit through the reference standards laboratory of Industrial Calibration Center, Technology Promotion Association (Thailand-Japan).

Instruments	Serial No.	ID No.	Certificate No.	Due Date
1. Burette	-	130BU10	25CG1126	18 Mar 2027
2. Balance	14233821	110RC001	24MM131	04 July 2025

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Result : Dissolved Oxygen Meter Adjustment With Air 100 %
Dissolved Oxygen Probe No.: 9K0LD009

Titration Method (Azide Modification Method) (mg/L)	DO Meter Reading (mg/L)	Standard Deviation (mg/L)
8.22	8.22	0.0071

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3 : EQUIPMENT CALIBRATION AND TESTING SERVICES
53/44 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
TEL. 0-2717-3000-29 FAX. 0-2719-9484

Certificate of Calibration

Cert. No.: 25LM80
Page.: 1 of 2

Equipment : DO Meter with Sensor
Manufacturer : Horiba
Model : LAQUA-DO210
Serial No. : HE2L0036
ID No. : UAE.EFM.018/2566(EFM.DO.03/66)
Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
3 Soi Udomsuk 41, Sukhumvit Road,
Bangkok, Phrakhanong,
Bangkok 10260
Location : TPA On Site Calibration Laboratory
Received Order : 20 May 2025
Calibrated Date : 22 May 2025
Ambient Temperature : (26 ± 10) °C
Relative Humidity : (50 ± 30) %
AC Line Voltage : (220 ± 22) V
Calibrated by : Preecha Hahib
Approved by : Kunchit
Approved Signatory
() Chakrit Waewwanjua
() Suwil Imjai
(✓) Kunchit Promprat
Issue Date : 28 May 2025

The Uncertainties are for a confidence probability of approximately 95%
This certificate may not be reproduced other than in full, except with the prior written
Approval of the head of Corporate Services 3 : Equipment Calibration and Testing Services.

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Equipment : DO Meter with Sensor
Condition As-Received : Used Item
Reference : 2505-0505WSC-8
Procedure Used :-

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This instrument was connected with thermistor sensor, S/N : 9K0L0009

Calibration Point (°C)	Immersion Depth (mm)	Standard Temperature (°C)	UUC* Reading (°C)	Error (°C)	Uncertainty (± °C)	Coverage Factor k
15.0	100	15.004	15.2	0.196	0.16	2.00
30.0	100	30.005	30.1	0.095	0.16	2.00
45.0	100	45.008	45.0	-0.008	0.16	2.00

UUC* : Unit Under Calibration

The reported uncertainty of measurement was based on a standard uncertainty multiplied by a coverage factor *k*, providing a level of confidence of approximately 95 %.

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List of Instruments Certification for Air

No.	Instrument/Equipment	Parameter	Manufacturer	Model/Serial No.	Calibrator	Certification No.	Date of Calibration	Due date of Calibration	Remark
1	DO Meter	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	20-Aug-25	19-Aug-26	
2	Speed Sensor	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	17-Apr-25	16-Apr-26	
3	DO Meter	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	10-Apr-25	9-Apr-26	
4	Speed Sensor	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	10-Apr-25	9-Apr-26	
5	Speed Sensor	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	10-Apr-25	9-Apr-26	
6	Speed Sensor	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	10-Apr-25	9-Apr-26	
7	Speed Sensor	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	10-Apr-25	9-Apr-26	

List of Instruments Certification for Water Quality Analysis

No.	Instrument/Equipment	Parameter	Manufacturer	Model/Serial No.	Calibrator	Certification No.	Date of Calibration	Due date of Calibration	Remark
1	DO Meter	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	20-Aug-25	19-Aug-26	
2	DO Meter	DO	Hioki	MB-9220-01M	Hioki International Technology Co., Ltd. (Thailand Branch)	250508116-0003	20-Aug-25	19-Aug-26	

MIRACLE INTERNATIONAL TECHNOLOGY CO., LTD.
214 Bangrak Rd. Bangrak Bangkok 10160
Tel: 0-2665-4647-8 Fax: 0-2665-4649 http://www.mit.co.th

CALIBRATION CERTIFICATE

Page 1 of 3

Certificate No. : I.202508116-0003

Date Issued : 21-Aug-25

Customer : United Analyst and Engineering Consultant Co., Ltd.
81 Soi Udomsuk 41, Sukhumvit Road, Bangkok, Phrakong, Bangkok 10260

Equipment : Mass Flow Meter

Manufacturer : Alicat Scientific

Model : MB-50SCCM-D/5M

Serial No. : 202984

ID No./Tag No. : UAE.BFM.1972562

Date Received : 20-Aug-25

Date Calibrated : 20-Aug-25

Calibrated by : Saruth Srirachikul

Calibration Method or Calibration Procedure Used

In-house method : CP-34 by comparison against mass flow calibrator.

This certificate is traceable to national standards, which realize the units of measurement according to the International System of Units (SI).

Result of Calibration

The reported uncertainty of measurement was based on standard uncertainty multiplied by a coverage factor *k* = 2, providing a level confidence approximately 95 percent

This certificate may not be reproduced other than in full except with the prior written approval of the Miracle International Technology Company Limited.

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Calibration Media : Air


Type : Mass Flowmeter

Unit Under Calibration Reference Condition : At atmospheric pressure and room temperature condition

Temperature (°C)	Pressure (kPa)	UUC Reading (ml/min)	STD Reading (ml/min)	Error (ml/min)	Uncertainty (± ml/min)
22.87	106.42	0.00	0.000 *	0.000	0.063
22.81	100.66	3.00	3.054	-0.054	0.13
22.81	100.97	3.33	3.372	-0.042	0.14
22.83	101.09	4.00	4.038	-0.038	0.16
22.81	101.91	9.00	9.071	-0.071	0.17
22.75	102.07	10.00	10.088	-0.088	0.19

Error - Unit Under Calibration - Standard

Marked * are not included in the NSC-ONS accreditation schedule for our laboratory

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Note : The actual flow rate is determined by the equation :

$$Q_{Meas} = Q_{Ref} \times \frac{P_{Ref}}{P_{Meas}} \times \frac{T_{Meas}}{T_{Ref}}$$

; Q = Flow rate
 ; P = Absolute pressure
 ; T = Absolute temperature
 ; Subscript "Meas" = Measurement condition
 ; Subscript "Ref" = Reference condition

Condition As-Received : Used Item

The measurement results and statements of conformity with specification only relate to the item calibrated.

Traceability of Certificate :

The International System of Units (SI) through

NIMT Certificate No. MW-0648-24 for Gas Flow meter Serial No. MS209179C/MS209179A, Due 02-Jul-26

NIMT Certificate No. MW-0747-24 for Gas Flow meter Serial No. MS209179B/MS209179A, Due 28-Jun-26

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TECHNOLOGY PROMOTION ASSOCIATION (THAILAND-JAPAN)
CORPORATE SERVICES 3: EQUIPMENT CALIBRATION AND TESTING SERVICES
53/4 PATTANAKARN ROAD SOI 18, SUANLUANG, SUANLUANG BANGKOK 10250
TEL.0-2717-3000-29 FAX.0-2719-9484

Certificate of Calibration

Cert.No.: 25CH736
Page.: 1 of 3

Equipment : pH Meter
Manufacturer : Horiba
Model : LAQUA-PH210
Serial No. : HA0D0081
ID No. : UAEEFM.074/2564(EFM.pH.07/64)
Condition As-Received: Used Item
Received Date : 18 June 2025
Calibration Date : 23 June 2025
Reference : 2506-0591WSC-3
Submitted by : United Analyst and Engineering Consultant Co.,Ltd.
 3 Soi Udomsuk 41, Sukhumvit Road, Bangchak,
 Phrakhanong, Bangkok 10260

Ambient Temperature : (25 ± 2.5) °C
Relative Humidity : (50 ± 15) %
Calibration Procedure : In-house method :
 - CP-CH5 by direct measurement with DC voltage
 standard and direct measurement with
 certified reference material (CRM)
 - CP-CH6 by comparison with temperature standard

Calibrated by : Waleak Sinithean
Approved by :
 Approved Signatory
 () Chakrit Waewwanjua
 () Ponpan Palpim
 (✓) Sathip Meangmai
Issue Date : 24 June 2025

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2. Certified Reference Materials :The measurement results are traceable to SI through Hach Lange GmbH Ltd.,
 Deutsche Akkreditierungsstelle, Accredited No.D-RM-15184-01-00
 :The measurement results are traceable to SI through CPA chem Ltd.,
 ANSI-ASQ National Accreditation Board, Accredited No. AR-1635

Buffer Solution	Manufacturer	Lot No.	Exp. date
pH 4.007	CPA chem	1066665	18 Jan 2027
pH 7.000	Hach Lange GmbH	C03232	02 Dec 2026
pH 10.010	CPA chem	1066669	18 Jan 2026

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Calibration Results

Function : pH Measurement

Performing three buffers standard curve by using buffer nominal pH (4,7)(7,10)

Unit Under Calibration	Standard pH Buffer Solution	Actual pH Reading	Actual mV Reading (mV)	Uncertainty of pH Measurement (±)	Coverage factor k
pH Electrode S/N.: -	4.007	4.01	176.5	0.0079	2.00
	7.000	6.98	0.8	0.0085	2.00
	7.000	7.00	-0.7	0.0095	2.00
	10.010	10.01	-172.7	0.0092	2.00

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
1. APHA, AWWA, WEF, Standard Methods for the Examination of Water and Wastewater, 24th ed. Washington, DC: APHA, 2023.
2. United States Environmental Protection Agency. Standards of Performance for New Stationary Sources. 40 CFR 60. Appendix A, 2023.

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3. United States...

3. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Acid Digestion of Sediments, Sludges, and Soils. SW-846 Method 3050B, 1996.
4. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Volatile Organic Compounds in Various Sample Matrices Using Equilibrium Headspace Analysis. SW-846 Method 5021A, 2014.
5. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Inductively Coupled Plasma-Optical Emission Spectrometry. SW-846 Method 6010D, 2018.
6. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Flame Atomic Absorption Spectrophotometry. SW-846 Method 7000B, 2007.
7. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry. SW-846 Method 8260D, 2018.
8. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Soil and Waste pH. SW-846 Method 9045D, 2004.

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ลำดับ	สารมลพิษ	วิธีวิเคราะห์
1	Aldrin	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
2	Arsenic	1) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
3	Barium	Digestion, Inductively Coupled Plasma Method ^(a)
4	α-BHC	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
5	β-BHC	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
6	δ-BHC	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
7	γ-BHC	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
8	Biochemical Oxygen Demand	1) 5-Day BOD Test, Azide Modification Method ^(a) 2) 5-Day BOD Test, Membrane Electrode Method ^(a)
9	Cadmium	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a) 1) Closed Reflux, Titrimetric Method ^(a) 2) Closed Reflux, Colorimetric Method ^(a) 3) Open Reflux, Titrimetric Method ^(a)
10	Chemical Oxygen Demand	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
11	Chlordane	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
12	Chromium	ADMI Weighted-Ordinate Spectrophotometric Method ^(a)
13	Color	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
14	Copper	1) Distillation, Colorimetric Method ^(a) 2) Total Cyanide after Distillation, by Flow Injection Analysis Method ^(a)
15	Cyanide	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
16	o,p'-DDT	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
17	4,4'-DDD	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
18	4,4'-DDE	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
19	4,4'-DDT	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
20	Dieldrin	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
21	Endosulfan I	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
22	Endosulfan II	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
23	Endosulfan sulfate	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
24	Endrin	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)

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UAE
UNION ANALYST AND ENGINEERING
CONSULTANT COMPANY LIMITED

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ลำดับ	สารเคมี	วิธีวิเคราะห์
25	Endrin aldehyde	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
26	Formaldehyde	Distillation, Colorimetric Method ⁽²⁾
27	Free Chlorine	1) Iodometric Method ^(a) 2) DPD Ferrous Titrimetric Method ^(a)
28	Heptachlor	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
29	Heptachlor Epoxide	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
30	Hexavalent Chromium	Colorimetric Method ^(a)
31	Lead	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
32	Manganese	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
33	Mercury	1) Digestion, Inductively Coupled Plasma Method ^(a) Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ^(a)
34	Methoxychlor	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
35	Nickel	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
36	Oil & Grease	1) Liquid-Liquid, Partition-Gravimetric Method ^(a) 2) Soxhlet Extraction Method ^(a)
37	pH	Electrometric Method ^(a)
38	Phenols	1) Distillation, Chloroform Extraction Method ^(a) 2) Distillation, Direct Photometric Method ^(a)
39	Selenium	1) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
40	Sulfide	1) Iodometric Method ^(a) 2) Methylene Blue Method ^(a)
41	Temperature	Laboratory and Field Methods ^(a)
42	Total Dissolved Solids	Dried at 180 °C ^(a)
43	Total Kjeldahl Nitrogen	Semi-Micro-Kjeldahl Method ^(a)
44	Total Suspended Solids	Dried from 103 to 105 °C ^(a)
45	Trivalent Chromium	1) Digestion, Direct Air-Acetylene Flame Method; Colorimetric Method, Calculation ^(a) 2) Digestion, Inductively Coupled Plasma Method; Colorimetric Method, Calculation ^(a)
46	Zinc	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)

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ลำดับ	สารเคมี	วิธีวิเคราะห์
1	Acenaphthene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
2	Acetone	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
3	Aldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
4	Anthracene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
5	Antimony	Digestion, Inductively Coupled Plasma Method ^(a)
6	Arsenic	1) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
7	Atrazine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
8	Barium	Digestion, Inductively Coupled Plasma Method ^(a)
9	Benz(a)anthracene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
10	Benzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
11	Benzo(b)fluoranthene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
12	Benzo(k)fluoranthene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
13	Benzolic acid	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)

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14 Benzo(a)pyrene...

ลำดับ	สารเคมี	วิธีวิเคราะห์
14	Benz(a)pyrene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
15	Benzo(g,h,i)perylene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
16	Beryllium	Digestion, Inductively Coupled Plasma Method ^(a)
17	Bis(2-chloroethyl)ether	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
18	Bis(2-ethylhexyl)phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
19	Bromodichloromethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
20	Bromoform	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
21	Butanol	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
22	Butyl benzyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
23	Cadmium	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Electrothermal Atomic Absorption Spectrometric Method ^(a)
24	Carbazole	3) Digestion, Inductively Coupled Plasma Method ^(a) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
25	Carbon disulfide	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
26	Carbon tetrachloride	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
27	Chlordane	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
28	p-Chloroaniline	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)

29 Chlorobenzene...

ลำดับ	สารเคมี	วิธีวิเคราะห์
29	Chlorobenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
30	Chlorodibromomethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
31	Chloroform	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
32	2-Chlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
33	Chromium	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
34	Chromium (III)	1) Digestion, Direct Air-Acetylene Flame Method; Colorimetric Method; Calculation ^(a) 2) Digestion, Inductively Coupled Plasma Method; Colorimetric Method; Calculation ^(a)
35	Chromium (VI)	Colorimetric Method ^(a)
36	Chrysene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
37	Cyanide	Distillation, Colorimetric Method ^(a)
38	2,4-D	Liquid-Liquid Extraction, Gas Chromatographic Method ^(a)
39	DDD	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
40	DDE	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
41	DDT	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
42	Dibenz(a,h)anthracene	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)

43 Di-n-butyl phthalate...

ลำดับ	สารเคมี	วิธีวิเคราะห์
43	Di-n-butyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
44	1,2-Dichlorobenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
45	1,3-Dichlorobenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
46	1,4-Dichlorobenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
47	3,3'-Dichlorobenzidine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
48	1,1-Dichloroethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
49	1,2-Dichloroethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
50	1,1-Dichloroethylene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
51	cis-1,2-Dichloroethylene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
52	trans-1,2-Dichloroethylene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
53	2,4-Dichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
54	1,2-Dichloropropane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
55	1,3-Dichloropropane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
56	1,3-Dichloropropene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
57	Dieldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
58	Diethyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
59	2,4-Dimethylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
60	2,4-Dinitrophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)

Signature

Stamp

61 2,4-Dinitrotoluene...

61 2,4-Dinitrotoluene...

ลำดับ	สารเคมี	วิธีวิเคราะห์
61	2,4-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
62	2,6-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
63	Di-n-Octyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
64	Endosulfan	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
65	Endrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
66	Ethylbenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
67	Fluoranthene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
68	Fluorene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
69	Heptachlor	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
70	Heptachlor epoxide	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
71	Hexachlorobenzene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
72	Hexachloro-1,3-butadiene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
73	n-Hexane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)

Signature

Stamp

61 2,4-Dinitrotoluene...

61 2,4-Dinitrotoluene...

ลำดับ	สารเคมี	วิธีวิเคราะห์
74	α-HCH	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
75	β-HCH	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
76	γ-HCH	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
77	Hexachlorocyclopentadiene	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
78	Hexachloroethane	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
79	Indeno(1,2,3-cd)pyrene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
80	Isophorone	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
81	Lead	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Electrothermal Atomic Absorption Spectrometric Method ^(a)
82	Manganese	3) Digestion, Inductively Coupled Plasma Method ^(a) 1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
83	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ^(a)
84	Methanol	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
85	Methoxychlor	1) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
86	Methyl bromide	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)

87 Methylene chloride...

ลำดับ	สารเคมี	วิธีวิเคราะห์
87	Methylene chloride	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
88	2-Methylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
89	2-Methylnaphthalene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
90	Methyl tert-butyl ether	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
91	Naphthalene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
92	Nickel	1) Digestion, Direct Air-Acetylene Flame Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
93	Nitrobenzene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
94	N-Nitrosodiphenylamine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
95	N-Nitrosodi-n-propylamine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
96	Polychlorinated Biphenyls - PCB 1016 - PCB 1221 - PCB 1232 - PCB-1242 - PCB-1248 - PCB-1254 - PCB-1260	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
97	Pentachlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
98	pH	Electrometric Method ^(a)
99	Phenanthrene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)

100 Phenol...

ลำดับ	สารเคมี	วิธีวิเคราะห์
100	Phenol	1) Distillation, Chloroform Extraction Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
101	Pyrene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
102	Selenium	1) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(a) 2) Digestion, Inductively Coupled Plasma Method ^(a)
103	Silver	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
104	Styrene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
105	1,1,2,2-Tetrachloroethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
106	Tetrachloroethylene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
107	Toluene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
108	Toxaphene	1) Liquid-Liquid Extraction, Gas Chromatographic Method ^(a) 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(a)
109	TPH (C ₅ - C ₆)	1) Purge and Trap, Gas Chromatographic Method ^(12,21) 2) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(12,21)
110	TPH (C ₅ - C ₁₆)	Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(9,22)
111	TPH (C ₅ - C ₃₅)	Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(9,22)
112	1,2,4-Trichlorobenzene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
113	1,1,1-Trichloroethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
114	1,1,2-Trichloroethane	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)
115	Trichloroethylene	Purge and Trap Gas Chromatographic/Mass Spectrometric Method ^(a)

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Chromium (ต่อ)...

ลำดับ	สารเคมี	วิธีวิเคราะห์
6	Chromium (ค่า)	2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
7	Cobalt	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
8	Copper	1) Isokinetic Sampling, Digestion, Direct Air-Acetylene Flame Method ^[5] 2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
9	Cresol	Absorption Sampling, Gas Chromatographic Method ^[5]
10	Dioxins/Furans	Isokinetic Sampling ^[5]
11	Hydrogen Chloride	Isokinetic Sampling, Ion Chromatographic Method ^[5]
12	Hydrogen Fluoride	Isokinetic Sampling, Ion Chromatographic Method ^[5]
13	Hydrogen Sulfide	Absorption Sampling, Iodometric Method ^[5]
14	Lead	1) Isokinetic Sampling, Digestion, Direct Air-Acetylene Flame Method ^[5] 2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
15	Manganese	1) Isokinetic Sampling, Digestion, Direct Air-Acetylene Flame Method ^[5] 2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
16	Mercury	Isokinetic Sampling, Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ^[5]
17	Nickel	1) Isokinetic Sampling, Digestion, Direct Air-Acetylene Flame Method ^[5] 2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
18	Opacity	Ringelmann's Method ^[5]
19	Oxides of Nitrogen	1) Absorption Sampling, Phenoldisulfonic acid Method ^[5] 2) Instrumental Analyzer Method ^[5]
20	Selenium	1) Isokinetic Sampling, Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^[5] 2) Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method ^[5]
21	Sulfur Dioxide	1) Absorption Sampling, Barium-Thorin Titrimetric Method ^[5] 2) Instrumental Analyzer Method ^[5]
22	Sulfuric Acid	Isokinetic Sampling, Barium-Thorin Titrimetric Method ^[5]

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ลำดับ	สารเคมี	วิธีวิเคราะห์
1	Aldrin	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^[3,9,23] 2) Ultrasonic Extraction, Gas Chromatographic Method ^[10,23]
2	Antimony	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^[3,6,14] 2) Digestion, Inductively Coupled Plasma Method ^[7,14]
3	Arsenic	1) Waste Extraction, Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^[3,6,16] 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^[3,6,14] 3) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^[7,16]
4	Barium	4) Digestion, Inductively Coupled Plasma Method ^[7,14] 1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^[3,6,14]
5	Beryllium	2) Digestion, Inductively Coupled Plasma Method ^[7,14] 1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^[3,6,14]
6	Cadmium	2) Digestion, Inductively Coupled Plasma Method ^[7,14] 1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method ^[3,6,15] 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^[3,6,14]
7	Chlordane	3) Digestion, Flame Atomic Absorption Spectrometric Method ^[7,15] 4) Digestion, Inductively Coupled Plasma Method ^[7,14] 1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^[3,9,23] 2) Ultrasonic Extraction, Gas Chromatographic Method ^[10,23]

8 Chromium...

ลำดับ	สารเคมี	วิธีวิเคราะห์
8	Chromium	1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method ^(3.6.15) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6.14) 3) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.1.5) 4) Digestion, Inductively Coupled Plasma Method ^(7.14)
9	Chromium (III)	1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method; Waste Extraction, Colorimetric Method; Calculation ^(3.6.15.17) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method; Waste Extraction, Colorimetric Method; Calculation ^(3.6.15.17) 3) Digestion, Flame Atomic Absorption Spectrometric Method; Alkaline Digestion, Colorimetric Method; Calculation ^(7.8.15.17) 4) Digestion, Inductively Coupled Plasma Method; Alkaline Digestion, Colorimetric Method; Calculation ^(7.8.15.17)
10	Chromium (VI)	1) Waste Extraction, Colorimetric Method ^(3.17) 2) Alkaline Digestion, Colorimetric Method ^(6.17)
11	Cobalt	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6.14)
12	Copper	2) Digestion, Inductively Coupled Plasma Method ^(7.14) 1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method ^(3.6.15) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6.14)
13	2,4-D	3) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.15) 4) Digestion, Inductively Coupled Plasma Method ^(7.14) 1) Waste Extraction, Gas Chromatographic Method ^(3.26) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(3.9.23)
14	DDD	Method ^(10.23) 1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)

R/M

๑๖๖

๑๖๖

๑๖๖

๑๖๖

๑๖๖

๑๖๖

ลำดับ	สารเคมี	วิธีวิเคราะห์
15	DDE	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
16	DDT	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
17	Dieldrin	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
18	Endrin	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
19	Heptachlor	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
20	Lead	1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method ^(3.6.15) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6.14) 3) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.15)
21	Lindane	4) Digestion, Inductively Coupled Plasma Method ^(7.14) 1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9.23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
22	Mercury	1) Waste Extraction, Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ^(3.19) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6.14) 3) Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ⁽¹⁹⁾ 4) Digestion, Inductively Coupled Plasma Method ^(7.14)

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Polychlorinated Biphenyls (ต่อไป)...

ลำดับ	สารเคมี	วิธีวิเคราะห์
	Polychlorinated Biphenyls (ต่อไป) - 2,2',3,4,4',5'-Hexachlorobiphenyl - 2,2',3,4,5,5'-Hexachlorobiphenyl - 2,2',3,5,5',6'-Hexachlorobiphenyl - 2,2',4,4',5,5'-Hexachlorobiphenyl - 2,2',3,3',4,4',5'-Heptachlorobiphenyl - 2,2',3,4,4',5,5'-Heptachlorobiphenyl - 2,2',3,4,4',5,6'-Heptachlorobiphenyl - 2,2',3,3',4,4',5,5',6'-Nonachlorobiphenyl Pentachlorophenol	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method ^(3,9,28) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28) Electrometric Method ^(3,1,32) 1) Waste Extraction, Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(3,6,21) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3,6,14) 3) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(7,21) 4) Digestion, Inductively Coupled Plasma Method ^(7,14)
27	pH Selenium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3,6,14) 2) Digestion, Inductively Coupled Plasma Method ^(3,6,14)
28	Silver	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3,6,14) 2) Digestion, Inductively Coupled Plasma Method ^(3,6,14)
29	Thallium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3,6,14) 2) Digestion, Inductively Coupled Plasma Method ^(7,14)

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32 Toxaphene...

ลำดับ	สารเคมี	วิธีวิเคราะห์
32	Toxaphene	1) Waste Extraction, Separatory Funnel Liquid-Liquid Extraction, Gas Chromatographic Method ^(3.9,23) 2) Ultrasonic Extraction, Gas Chromatographic Method ^(10,23) 1) Waste Extraction, Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(3.12,27) 2) Waste Extraction, Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(3.11,27) 3) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27) 4) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11,27)
33	Trichloroethylene	1) Waste Extraction, Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(3.12,27) 2) Waste Extraction, Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(3.11,27) 3) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27) 4) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11,27)
34	Vanadium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6,14) 2) Digestion, Inductively Coupled Plasma Method ^(7.14)
35	Zinc	1) Waste Extraction, Digestion, Flame Atomic Absorption Spectrometric Method ^(3.6,15) 2) Waste Extraction, Digestion, Inductively Coupled Plasma Method ^(3.6,14) 3) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.15) 4) Digestion, Inductively Coupled Plasma Method ^(7.14)

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Heptachlor epoxide (๓๖)...

ลำดับ	สารเคมี	วิธีวิเคราะห์
70	Heptachlor epoxide (๗๖)	2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
71	Hexachlorobenzene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
72	Hexachloro-1,3-butadiene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
73	n-Hexane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
74	α-HCH	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23)
75	β-HCH	2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
76	γ-HCH	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
77	Hexachlorocyclopentadiene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
78	Hexachloroethane	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
79	Indeno(1,2,3-cd)pyrene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
80	Isophorone	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
81	Lead	1) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.14) 2) Digestion, Inductively Coupled Plasma Method ^(7.14)
82	Manganese	1) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.13) 2) Digestion, Inductively Coupled Plasma Method ^(7.14)

83 Mercury...

ลำดับ	สารเคมี	วิธีวิเคราะห์
83	Mercury	1) Digestion, Cold-Vapor Atomic Absorption Spectrometric Method ⁽¹⁹⁾ 2) Thermal Decomposition Amalgamation and Atomic Absorption Spectrometric Method ⁽¹⁹⁾
84	Methanol	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
85	Methoxychlor	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
86	Methyl bromide	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
87	Methylene chloride	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)
88	2-Methylphenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
89	2-Methylnaphthalene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
90	Methyl tert-butyl ether	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
91	Naphthalene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.23) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
92	Nickel	1) Digestion, Flame Atomic Absorption Spectrometric Method ^(7.13) 2) Digestion, Inductively Coupled Plasma Method ^(7.14)
93	Nitrobenzene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
94	N-Nitrosodiphenylamine	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
95	N-Nitrosodi-n-propylamine	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
96	Polychlorinated Biphenyls - Aroclor 1016	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10.24)

Polychlorinated Biphenyls(๓๖)...

ลำดับ	สารเคมี	วิธีวิเคราะห์
96	Polychlorinated Biphenyls(พีซี) - Aroclor 1221 - Aroclor 1232 - Aroclor 1242 - Aroclor 1248 - Aroclor 1254 - Aroclor 1260 Polychlorinated Biphenyls - 2-Chlorobiphenyl - 2,3-Dichlorobiphenyl - 2,2',5-Trichlorobiphenyl - 2,4',5-Trichlorobiphenyl - 2,2',3,5'-Tetrachlorobiphenyl - 2,2',5,5'-Tetrachlorobiphenyl - 2,3',4,4'-Tetrachlorobiphenyl - 2,2',3,4,5'- Pentachlorobiphenyl - 2,2',4,5,5'- Pentachlorobiphenyl - 2,3',3',4',6'- Pentachlorobiphenyl - 2,2',3,4,4',5'- Hexachlorobiphenyl - 2,2',3,4,5,5'- Hexachlorobiphenyl - 2,2',3,5,5',6'- Hexachlorobiphenyl - 2,2',4,4',5,5'- Hexachlorobiphenyl - 2,2',3,3',4,4',5,5'- Heptachlorobiphenyl - 2,2',3,4,4',5,5'- Heptachlorobiphenyl - 2,2',3,4,4',5',6- Heptachlorobiphenyl - 2,2',3,4',5,5',6- Heptachlorobiphenyl - 2,2',3,3',4,4',5,5',6- Nonachlorobiphenyl	2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28) Ultrasonic Extraction, Gas Chromatographic Method ^(10,24)

WAE
UNITED ANALYST AND ENGINEERING
CONSULTANT COMPANY LIMITED

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97 Pentachlorophenol...

ลำดับ	สารเคมี	วิธีวิเคราะห์
97	Pentachlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28)
98	Phenanthrene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10,25) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28)
99	Phenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28)
100	Pyrene	1) Ultrasonic Extraction, Gas Chromatographic Method ^(10,25) 2) Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10,28)
101	Selenium	1) Digestion, Hydride Generation/Atomic Absorption Spectrometric Method ^(7,21) 2) Digestion, Inductively Coupled Plasma Method ^(7,14)
102	Silver	Digestion, Inductively Coupled Plasma Method ^(7,14)
103	Styrene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11,27)
104	1,1,1,2,2-Tetrachloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27)
105	Tetrachloroethylene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11,27)
106	Toluene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11,27)
107	Toxaphene	Ultrasonic Extraction, Gas Chromatographic Method ^(10,25)
108	TPH (C ₈ -C ₁₆)	1) Purge and Trap, Gas Chromatographic Method ^(13,27) 2) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13,27)
109	TPH (C ₈ -C ₁₆)	Ultrasonic Extraction, Gas Chromatographic Method ^(10,22)
110	TPH (C ₁₅ -C ₃₃)	Ultrasonic Extraction, Gas Chromatographic Method ^(10,22)

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ลำดับ	สารเคมี	วิธีวิเคราะห์
111	1,2,4-Trichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
112	1,1,1-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
113	1,1,2-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
114	Trichloroethylene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)
115	2,4,5-Trichlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
116	2,4,6-Trichlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method ^(10.28)
117	1,3,5-Trimethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
118	Vanadium	Digestion, Inductively Coupled Plasma Method ^(7.14)
119	Vinyl acetate	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
120	Vinyl chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27)
121	m-Xylene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)
122	o-Xylene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)
123	p-Xylene	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)
124	Xylene (Total)	1) Purge and Trap, Gas Chromatographic/Mass Spectrometric Method ^(13.27) 2) Equilibrium Headspace, Gas Chromatographic/Mass Spectrometric Method ^(11.27)

Signature: *[Signature]*
Name: **สมชาย งามวิจิตร**
Title: **Senior Analyst**
Company: **CONSULTANT COMPANY LIMITED**

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- กระทรวงอุตสาหกรรม. ประกาศกระทรวงอุตสาหกรรม, พ.ศ. 2549. เรื่อง กำหนดค่าปริมาณขั้นต่ำที่เจือปนในอากาศที่ระเหยจากปล่องของหม้อไอน้ำที่ใช้ถ่านเป็นเชื้อเพลิง. ราชกิจจานุเบกษา. 4 ธันวาคม 2549 เล่มที่ 123 ตอนพิเศษ 125จ.
- สมาคมวิศวกรรมสิ่งแวดล้อมแห่งประเทศไทย. คู่มือวิเคราะห์น้ำเสีย. พิมพ์ครั้งที่ 4. กรุงเทพฯ: เรือนแก้วการพิมพ์. 2547.
- กระทรวงอุตสาหกรรม. ประกาศกระทรวงอุตสาหกรรม, พ.ศ. 2566. เรื่อง การจัดการสิ่งปฏิกูลหรือวัสดุที่ไม่ใช้แล้ว. ราชกิจจานุเบกษา. 31 พฤษภาคม 2566. เล่มที่ 140 ตอนพิเศษ 126 จ.
- APHA, AWWA, WEF. Standard Methods for the Examination of Water and Wastewater. 24th ed. Washington, DC: APHA, 2023.
- United States Environmental Protection Agency. Standards of Performance for New Stationary Sources. 40 CFR 60. Appendix A, 2020.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. SW-846, 2014.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods Acid Digestion of Sediments, Sludges, and Soils. SW-846 Method 3050B, 1996.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Alkaline Digestion for Hexavalent Chromium. SW-846 Method 3060A, 1996.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Separatory Funnel Liquid-Liquid Extraction. SW-846 Method 3510C, 1996.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Ultrasonic Extraction. SW-846 Method 3550C, 2007.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Volatile Organic Compounds in Various Sample Matrices Using Equilibrium Headspace Analysis. SW-846 Method 5021A, 2014.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Purge-and-Trap for Aqueous Samples. SW-846 Method 5030C, 2003.
- United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Closed System Purge and Trap and Extraction for Volatile Organics in Soil and Waste Sample. SW-846 Method 5035A, 2000.
- United States...

Signature: *[Signature]*
Name: **สมชาย งามวิจิตร**
Title: **Senior Analyst**
Company: **CONSULTANT COMPANY LIMITED**

14. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Inductively Coupled Plasma-Optical Emission Spectrometry. SW-846 Method 6010D, 2014.
15. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Flame Atomic Absorption Spectrophotometry. SW-846 Method 7000B, 2007.
16. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Arsenic (Atomic Absorption, Gaseous Hydride). SW-846 Method 7061A, 1992.
17. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Chromium, Hexavalent (Colorimetric). SW-846 Method 7196A, 1992.
18. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Mercury in Liquid Waste (Manual Cold Vapor Technique). SW-846 Method 7470A, 1994.
19. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Mercury in Solid or Semisolid Waste (Manual Cold-Vapor Technique). SW-846 Method 7471B, 1998.
20. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Mercury in Solids and Solutions by Thermal Decomposition, Amalgamation, and Atomic Absorption Spectrophotometry. SW-846 Method 7473, 2007.
21. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Selenium (Atomic Absorption, Borohydride Reduction). SW-846 Method 7742, 1994.
22. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Nonhalogenated Organics Using GC/FID. SW-846 Method 8015D, 2003.
23. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Organochlorine Pesticides by Gas Chromatography. SW-846 Method 8081B, 2007.
24. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Polychlorinated Biphenyls (PCBs) by Gas Chromatography. SW-846 Method 8082A, 2007.
25. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Organochlorine Pesticides by Gas Chromatography. SW-846 Method 8081B, 2007.
26. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. Chlorinated Herbicides by GC Using Methylation or Pentafluorobenzoylation Derivatization. SW-846 Method 8151A, 1996.
27. United States...