



ที่ อก ๐๓๒๐/๑๑๓๕๒

กรมโรงงานอุตสาหกรรม  
ถนนพระรามที่ ๖ แขวงทุ่งพญาไท  
เขตราชเทวี กรุงเทพฯ ๑๐๕๐๐

๒๗ กรกฎาคม ๒๕๖๖

เรื่อง ค่ออายุหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน

เรียน กรรมการผู้จัดการ บริษัท อีสเทิร์น ไทย คอนซัลตติ้ง ๑๙๙๒ จำกัด

อ้างถึง คำขอค่ออายุของห้องปฏิบัติการวิเคราะห์เอกชน ลงวันที่ ๗ มิถุนายน ๒๕๖๖

สิ่งที่ส่งมาด้วย เอกสารแนบท้ายหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน

๑. รายชื่อผู้ควบคุมดูแลห้องปฏิบัติการวิเคราะห์ จำนวน ๔๐ ราย
๒. รายชื่อเจ้าหน้าที่ประจำห้องปฏิบัติการวิเคราะห์ จำนวน ๒๕ ราย
๓. ขอบข่ายสารมลพิษที่ได้รับขึ้นทะเบียนจากกรมโรงงานอุตสาหกรรม จำนวน ๒๙๒ รายการ จำนวน ๑๙ แผ่น

ตามหนังสือที่อ้างถึง บริษัท อีสเทิร์น ไทย คอนซัลตติ้ง ๑๙๙๒ จำกัด ขอค่ออายุหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน เลขทะเบียน ๖-๐๐๓ สถานที่ตั้งเลขที่ ๖๘๓ หมู่ที่ ๑๑ ถนนสุขุมวิท ๘ ตำบลหนองแขม อำเภอศรีราชา จังหวัดชลบุรี ต่อกรมโรงงานอุตสาหกรรม นั้น

กรมโรงงานอุตสาหกรรมพิจารณาแล้ว ให้บริษัท อีสเทิร์น ไทย คอนซัลตติ้ง ๑๙๙๒ จำกัด ค่ออายุหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน โดยมีองค์ประกอบดังนี้

- ก. ผู้ควบคุมดูแลห้องปฏิบัติการวิเคราะห์ จำนวน ๔๐ ราย ตามสิ่งที่ส่งมาด้วย ๑
- ข. เจ้าหน้าที่ประจำห้องปฏิบัติการวิเคราะห์ จำนวน ๒๕ ราย ตามสิ่งที่ส่งมาด้วย ๒
- ค. ขอบข่ายสารมลพิษที่ได้รับขึ้นทะเบียนให้วิเคราะห์ในน้ำเสีย จำนวน ๔๗ รายการ อากาศเสีย (ปล่องระบาย) จำนวน ๒๒ รายการ น้ำใต้ดิน จำนวน ๑๑๑ รายการ สิ่งปฏิกูลหรือวัสดุที่ไม่ใช้แล้ว จำนวน ๑๘ รายการ และดิน จำนวน ๙๕ รายการ รวมทั้งสิ้นจำนวน ๒๙๒ รายการ ตามสิ่งที่ส่งมาด้วย ๓

หนังสือฉบับนี้จะหมดอายุในวันที่ ๕ กรกฎาคม ๒๕๖๙ หากประสงค์จะค่ออายุหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน ให้ยื่นคำขอค่ออายุพร้อมเอกสารประกอบคำขอต่อกรมโรงงานอุตสาหกรรม ภายใน ๓๐ วัน ก่อนวันสิ้นอายุของหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน ทั้งนี้ สามารถยื่นคำขอผ่านระบบอิเล็กทรอนิกส์ได้ที่หน้าเว็บไซต์กรมโรงงานอุตสาหกรรม

จึงเรียนมาเพื่อทราบ

ขอแสดงความนับถือ

(นายทวี อำพาพันธ์)

ผู้อำนวยการศูนย์วิจัยและเตือนภัยมลพิษโรงงานภาคตะวันออก

ปฏิบัติราชการแทนอธิบดีกรมโรงงานอุตสาหกรรม

ศูนย์วิจัยและเตือนภัยมลพิษโรงงานภาคตะวันออก

โทร. ๐ ๓๓๑๓ ๖๐๕๙ ต่อ ๕๐๐๓-๒

ไปรษณีย์อิเล็กทรอนิกส์ ewr@diw.mail.go.th



"อุตสาหกรรมก้าวไกล ประเทศไทยก้าวหน้า ร่วมกันพัฒนา อุตสาหกรรมสีเขียว"

COPY



เอกสารแนบท้ายหนังสือรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน

บริษัท อีสเทิร์น ไทย คอนซัลตติ้ง ๑๙๙๒ จำกัด เลขทะเบียน ๖-๐๐๓

ที่ อก ๐๓๒๐/๑๑๓๕๒

ลงวันที่ ๒๗ กรกฎาคม ๒๕๖๖

ก. ผู้ควบคุมดูแลห้องปฏิบัติการวิเคราะห์ จำนวน ๔๐ ราย

๑) นางสาวมาลีเกษ เลพะวัจกุล	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๑
๒) นายวัฒนา โคตรกล้า	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๒
๓) นางวรรณเพ็ญ เหลาจินดาวัฒน์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๓
๔) นายกะวีร์ สุธาทิตย์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๔
๕) นางสาวนันท์นภัส แบนุนทด	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๕
๖) นางสาวพรนภา หลงคำหังษ์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๖
๗) นางสาวอภิลิ ขันอารมย์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๗
๘) นางสาวอัจฉรี จิตตะยโสธร	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๘
๙) นางสาวจิรพร ปานคง	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๐๙
๑๐) นายสุทธา สองธนิย์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๐
๑๑) นางสาวนันทประภา อุยสูงเนิน	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๑
๑๒) นายธงไชย บุญศักดิ์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๒
๑๓) นางสาวนันทพร กลิ่นโสภณ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๓
๑๔) นายธีระพงษ์ นวลอินทร์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๔
๑๕) นางสาวแพรว พลเสน	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๕
๑๖) นายทรงพล ผิวอ้วน	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๖
๑๗) นายภาคภูมิ บัวสวัสดิ์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๗
๑๘) นางสาวจันทน์ สายพันธ์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๘
๑๙) นายภาณุพงศ์ บำรุงรส	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๑๙
๒๐) นางสาวปภาณิน จันต๊ะสอน	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๐
๒๑) นายวรากร ไหวทะเสวี	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๑
๒๒) นางสาววรรณภา ไชยศิริ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๒
๒๓) นางสาวพรพิมล ภูมิคอนสาร	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๓
๒๔) นางสาวธมลวรรณ ผลอ้อ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๔
๒๕) นางสาวบุญเรือง บุญถม	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๕
๒๖) นางสาวกสิณันท์ บ่อมน้อย	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๖
๒๗) นายชานูวัฒน์ โชติวงค์	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๗
๒๘) นางสาวพจณีย์ งามวิสัย	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๘
๒๙) นายวิญญ์สวัสดิ์ สิงโต	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๒๙
๓๐) นางสาวนุกุล อารมศรี	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๐
๓๑) นายศุภฤกษ์ พาดกลาง	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๑
๓๒) นายณิชาพล ทองหล่อ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๒
๓๓) นายธรรมรัตน์ โพธิ์ตันคำ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๓
๓๔) นายโอชา ขวัญศิริมงคล	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๔
๓๕) นายเมธี สุขประเสริฐ	ทะเบียนเลขที่	๖-๐๐๓-ค-๐๐๓๕

COPY

๓๖) นางสาวพรทิพย์ วัลย์กุลกุล	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๓๖
๓๗) นางสาวอภิญญา เสริมสนธิ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๓๗
๓๘) นางสาวนันทพรธมภ์ ประดิษฐ์บุษ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๓๘
๓๙) นางสาวสุวิษา เอ็งเสียง	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๓๙
๔๐) นางสาวระพีณ อินัน	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๔๐

ข. เจ้าหน้าที่ประจำห้องปฏิบัติการวิเคราะห์ จำนวน ๒๕ ราย

๑) นางสาวดวงกมล เนื้อทอง	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๑
๒) นางสาววิษรารักษ์ อินทสุข	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๒
๓) นางสาวกัญจน์ธวิภา จันทร์ชอดแก้ว	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๓
๔) นางสาวณัฏฐา มงคลโกชน	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๔
๕) นางสาวณัฐวิภา อัมมัทธน์	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๕
๖) นางสาวนิอรธมา ประ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๖
๗) นางสาวธัญลักษณ์ ชื่นไธ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๗
๘) นางสาวสุทธิดา สร้างแก้ว	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๘
๙) นายอุดมทรัพย์ เจริญจริง	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๐๙
๑๐) นายนาธิป สงวนศิลป์	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๐
๑๑) นายวิรัชชัย พอใจ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๑
๑๒) นายอัญชลี ทะพงษ์	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๒
๑๓) นางสาวสุมิตรรา มีแก่น	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๓
๑๔) นางสาวสรวรรยา เพชรประไพ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๔
๑๕) นางสาวจุฬามาศ เจริญพรหม	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๕
๑๖) นางสาวนิภาพร คำชมภู	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๖
๑๗) นางสาวอรุษา พันธเมือง	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๗
๑๘) นายกิตติ ไทโรจน์	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๘
๑๙) นายชาญณรงค์ ตั้งธรรมรักษ์	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๑๙
๒๐) นางสาวปาริศา เอื้อนทิยะ	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๐
๒๑) นางสาวจุฑาทิพย์ กิจดี	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๑
๒๒) นางสาวสุภาวดี ศรีละออง	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๒
๒๓) นางสาวณัฐชยา บรรพบุตร	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๓
๒๔) นางสาวณัฏฐนิช นนตานอก	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๔
๒๕) นางสาวดวงสุดา แสนวันดี	ทะเบียนเลขที่	๗-๐๐๓-๗-๐๐๒๕

COPY

เอกสารแนบท้ายหนังสือรับต่ออายุขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชน  
บริษัท อีสเทิร์น ไทย คอนซัลติง ๑๙๙๒ จำกัด เลขทะเบียน ๗-๐๐๓  
ที่ ๐๓๒๐/๑๑๓๔๒ ลงวันที่ ๒๗ กรกฎาคม ๒๕๖๖

ขอขยายสารมลพิษที่ได้รับขึ้นทะเบียนจากกรมโรงงานอุตสาหกรรม จำนวน ๒๕๒ รายการ

น้ำเสีย จำนวน ๔๗ รายการ

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
1	Aldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
2	Arsenic	1) Continuous Hydride Generation/Atomic Absorption Spectrometric Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
3	Barium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
4	α-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
5	β-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
6	δ-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
7	γ-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
8	Biochemical Oxygen Demand	1) 5-Day BOD Test, Membrane Electrode Method <sup>[4]</sup> 2) 5-Day BOD Test, Azide Modification Method <sup>[4]</sup>
9	Cadmium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
10	Chemical Oxygen Demand	Closed Reflux, Titrimetric Method <sup>[4]</sup>
11	cis-Chlordane	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
12	trans-Chlordane	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
13	Chromium	1) Digestion, Direct Air-Acetylene Flame Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>

COPY

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
14	Color	ADMI Weighted-Ordinate Spectrophotometric Method <sup>[4]</sup>
15	Copper	1) Digestion, Direct Air-Acetylene Flame Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
16	Cyanide	Distillation, Colorimetric Method <sup>[4]</sup>
17	4,4'-DDD	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
18	4,4'-DDE	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
19	DDT	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
20	Dieldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
21	Endosulfan I	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
22	Endosulfan II	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
23	Endosulfan sulfate	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
24	Endrin	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
25	Endrin aldehyde	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
26	Endrin ketone	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
27	Formaldehyde	Distillation, Colorimetric Method <sup>[3]</sup>
28	Free Chlorine	1) Iodometric Method <sup>[4]</sup> 2) Colorimetric Method <sup>[4]</sup>

COPY

29 Heptachlor...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
29	Heptachlor	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
30	Heptachlor Epoxide	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
31	Hexavalent Chromium	Filtration, Colorimetric Method <sup>[4]</sup>
32	Lead	1) Digestion, Direct Air-Acetylene Flame Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
33	Manganese	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
34	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>[4]</sup>
35	Methoxychlor	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[4]</sup>
36	Nickel	1) Digestion, Direct Air-Acetylene Flame Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
37	Oil and Grease	Liquid-Liquid, Partition-Gravimetric Method <sup>[4]</sup>
38	pH	Electrometric Method <sup>[4]</sup>
39	Phenols	Distillation, Direct Photometric Method <sup>[4]</sup>
40	Selenium	Digestion, Hydride Generation/Atomic Absorption Spectrometric Method <sup>[4]</sup>
41	Sulfide	ZnS Precipitation, Iodometric Method <sup>[4]</sup>
42	Temperature	Field Method <sup>[4]</sup>
43	Trivalent Chromium	1) Digestion, Direct Air-Acetylene Flame Method; Filtration, Colorimetric Method; Calculation <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method; Colorimetric Method; Calculation <sup>[4]</sup>
44	Total Dissolved Solids	Dried at 180 °C <sup>[4]</sup>
45	Total Kjeldahl Nitrogen	Macro Kjeldahl Method <sup>[4]</sup>
46	Total Suspended Solids	Dried at 103-105 °C <sup>[4]</sup>
47	Zinc	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>

COPY

อากาศเสีย...

อากาศเสีย (ปล่อยระบาย) จำนวน 21 รายการ

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
1	Antimony	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
2	Arsenic	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
3	Cadmium	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
4	Carbon Monoxide	1) Bag, Non-Dispersive Infrared Method <sup>[5]</sup> 2) Instrumental Analyzer Method <sup>[5]</sup>
5	Chromium	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
6	Cobalt	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
7	Copper	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
8	Hydrogen Sulfide	Absorption Sampling, Iodometric Method <sup>[5]</sup>
9	Lead	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
10	Manganese	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
11	Mercury	Isokinetic Sampling, Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>[5]</sup>
12	Nickel	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
13	Opacity	Ringelmann's Method <sup>[1,5]</sup>
14	Oxides of Nitrogen	1) Absorption Sampling, Phenoldisulfonic Acid Method <sup>[8]</sup> 2) Instrumental Analyzer Method <sup>[7]</sup>
15	Selenium	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
16	Sulfur Dioxide	1) Absorption Sampling, Barium-Thorin Titrimetric Method <sup>[5]</sup> 2) Instrumental Analyzer Method <sup>[5]</sup>
17	Sulfuric Acid	Isokinetic Sampling, Barium-Thorin Titrimetric Method <sup>[5]</sup>
18	Tin	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>

COPY

19 Total Suspended Particulate...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
19	Total Suspended Particulate	Isokinetic Sampling, Gravimetric Method <sup>[6]</sup>
20	Vanadium	Isokinetic Sampling, Digestion, Inductively Coupled Plasma Method <sup>[5]</sup>
21	Xylene	Adsorption Sampling, Gas Chromatographic Method <sup>[6]</sup>

น้ำใต้ดิน จำนวน 111 รายการ

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
1	Acenaphthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
2	Acetone	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
3	Aldrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
4	Anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
5	Antimony	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
6	Arsenic	1) Continuous Hydride Generation/Atomic Absorption Spectrometric Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
7	Barium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
8	Benz(a)anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
9	Benzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
10	Benzo(b)fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
11	Benzo(k)fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
12	Benzo(a)pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
13	Benzo[g,h,i]perylene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
14	Beryllium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>

COPY

15 Bis(2-chloroethyl)ether...



ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
15	Bis(2-chloroethyl)ether	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
16	Bis(2-ethylhexyl)phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
17	Bromodichloromethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
18	Bromoform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
19	Butanol	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
20	Butyl benzyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
21	Cadmium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
22	Carbazole	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
23	Carbon disulfide	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
24	Carbon tetrachloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
25	Chlordane	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
26	p-Chloroaniline	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
27	Chlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
28	Chlorodibromomethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
29	Chloroform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
30	2-Chlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
31	Chromium	1) Digestion, Direct Air-Acetylene Flame Method <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
32	Chromium (III)	1) Digestion, Direct Air-Acetylene Flame Method; Filtration, Colorimetric Method; Calculation <sup>[4]</sup> 2) Digestion, Inductively Coupled Plasma Method; Filtration, Colorimetric Method; Calculation <sup>[4]</sup>

COPY

33 Chromium (VI)

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
33	Chromium (VI)	Filtration, Colorimetric Method <sup>[4]</sup>
34	Chrysene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
35	Cyanide	Distillation, Colorimetric Method <sup>[4]</sup>
36	DDD	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
37	DDE	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
38	DDT	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
39	Dibenz(a,h)anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
40	Di-n-butyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
41	1,2-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
42	1,3-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
43	1,4-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
44	1,1-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
45	1,2-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
46	1,1-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
47	cis-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
48	trans-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
49	2,4-Dichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
50	1,2-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
51	1,3-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>

COPY

52 Dieldrin...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
52	Dieldrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
53	Diethyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
54	2,4-Dimethylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
55	2,4-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
56	2,6-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
57	Di-n-octyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
58	Endosulfan	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
59	Endrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
60	Ethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
61	Fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
62	Fluorene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
63	Heptachlor	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
64	Heptachlor epoxide	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
65	Hexachlorobenzene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
66	Hexachloro-1,3-butadiene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
67	n-Hexane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
68	$\alpha$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
69	$\beta$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>

COPY

70  $\gamma$ -HCH...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
70	$\gamma$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
71	Hexachlorocyclopentadiene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
72	Hexachloroethane	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
73	Indeno(1,2,3-cd)pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
74	Isophorone	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
75	Lead	1) Digestion, Direct Air-Acetylene Flame Method <sup>(4)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>
76	Manganese	Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>
77	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>(4)</sup>
78	Methoxychlor	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
79	Methylene chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
80	2-Methylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
81	2-Methylnaphthalene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
82	Methyl tert-butyl ether	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
83	Naphthalene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
84	Nickel	Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>
85	Nitrobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
86	N-Nitrosodi-n-propylamine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
87	pH	Electrometric Method <sup>(4)</sup>
88	Phenanthrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>

COPY

89 Phenol...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
89	Phenol	1) Distillation, Direct Photometric Method <sup>(4)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
90	Pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
91	Selenium	Digestion, Hydride Generation/Atomic Absorption Spectrometric Method <sup>(4)</sup>
92	Silver	Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>
93	Styrene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
94	1,1,2,2-Tetrachloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
95	Trichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
96	1,2,4-Trichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
97	1,1,1-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
98	1,1,2-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
99	Tetrachloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
100	Toluene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
101	2,4,5-Trichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
102	2,4,6-Trichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
103	1,3,5-Trimethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
104	Vanadium	Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>
105	Vinyl acetate	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
106	Vinyl chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>

COPY

107 m-Xylene

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
107	m-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
108	o-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
109	p-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
110	Xylene (Total)	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(4)</sup>
111	Zinc	Digestion, Inductively Coupled Plasma Method <sup>(4)</sup>

สิ่งปฏิกูลหรือวัสดุที่ไม่ใช้แล้ว จำนวน 18 รายการ

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
1	Antimony	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
2	Arsenic	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
3	Barium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
4	Beryllium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
5	Cadmium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
6	Chromium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
7	Chromium (VI)	1) Waste Extraction, Digestion, Colorimetric Method <sup>(2,13)</sup> 2) Alkaline Digestion, Colorimetric Method <sup>(9,13)</sup>
8	Cobalt	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
9	Copper	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>(2,9,10)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>

COPY

10 Lead

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
10	Lead	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
11	Mercury	1) Waste Extraction, Digestion, Cold Vapor Atomic Absorption Spectrometric Method <sup>[2,11]</sup> 2) Digestion, Cold vapor Atomic Absorption Spectrometric Method <sup>[9,11]</sup>
12	Nickel	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
13	Molybdenum	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
14	Selenium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
15	Silver	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
16	Thallium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
17	Vanadium	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
18	Zinc	1) Waste Extraction, Digestion, Inductively Coupled Plasma Method <sup>[2,9,10]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>

COPY

ดิน...

## ดิน จำนวน 95 รายการ

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
1	Acenaphthene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
2	Acetone	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>
3	Anthracene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
4	Antimony	Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
5	Arsenic	Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
6	Barium	Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
7	Benz(a)anthracene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
8	Benzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>
9	Benzo(b)fluoranthene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
10	Benzo(k)fluoranthene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
11	Benzo(a)pyrene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
12	Benzo(g,h,i)perylene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
13	Beryllium	Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>
14	Bis(2-chloroethyl)ether	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
15	Bis(2-ethylhexyl)phthalate	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[15,17]</sup>
16	Bromodichloromethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>
17	Bromoform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>
18	Butanol	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>

COPY

19 Butyl benzyl phthalate..



ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
19	Butyl benzyl phthalate	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
20	Cadmium	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
21	Carbazole	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
22	Carbon disulfide	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
23	Carbon tetrachloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
24	p-Chloroaniline	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
25	Chlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
26	Chlorodibromomethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
27	Chloroform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
28	2-Chlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
29	Chromium	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
30	Chromium (III)	Digestion, Inductively Coupled Plasma Method; Filtration, Colorimetric Method; Calculation <sup>(9,10)</sup>
31	Chromium (VI)	Alkaline Digestion, Colorimetric Method <sup>(12,13)</sup>
32	Chrysene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
33	Dibenz(a,h)anthracene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
34	Di-n-butyl phthalate	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
35	1,2-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
36	1,3-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
37	1,4-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>

COPY

38 1,1-Dichloroethane

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
38	1,1-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
39	1,2-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
40	1,1-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
41	cis-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
42	trans-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
43	2,4-Dichlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
44	1,2-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
45	1,3-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
46	Diethyl phthalate	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
47	2,4-Dimethylphenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
48	2,4-Dinitrotoluene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
49	2,6-Dinitrotoluene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
50	Di-n-octyl phthalate	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
51	Ethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
52	Fluoranthene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
53	Fluorene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
54	Hexachlorobenzene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
55	Hexachloro-1,3-butadiene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>

COPY

56 n-Hexane

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
56	n-Hexane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
57	Hexachlorocyclopentadiene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
58	Hexachloroethane	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
59	Indeno(1,2,3-cd)pyrene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
60	Isophorone	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
61	Lead	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
62	Manganese	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
63	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>(9,11)</sup>
64	Methylene chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
65	2-Methylphenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
66	2-Methylnaphthalene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
67	Methyl tert-butyl ether	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
68	Naphthalene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
69	Nickel	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
70	Nitrobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
71	N-Nitrosodi-n-propylamine	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
72	Phenanthrene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
73	Phenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
74	Pyrene	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>

COPY

75 Selenium...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
75	Selenium	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
76	Silver	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
77	Styrene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
78	1,1,2,2-Tetrachloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
79	Tetrachloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
80	Toluene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
81	1,2,4-Trichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
82	1,1,1-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
83	1,1,2-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
84	Trichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
85	2,4,5-Trichlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
86	2,4,6-Trichlorophenol	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(15,17)</sup>
87	1,3,5-Trimethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
88	Vanadium	Digestion, Inductively Coupled Plasma Method <sup>(9,10)</sup>
89	Vinyl acetate	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
90	Vinyl chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
91	m-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
92	o-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>
93	p-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(14,16)</sup>

COPY

94 Xylene (Total)...

ลำดับที่	สารมลพิษ	วิธีวิเคราะห์
94	Xylene (Total)	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[14,16]</sup>
95	Zinc	Digestion, Inductively Coupled Plasma Method <sup>[9,10]</sup>

## เอกสารอ้างอิง

1. กระทรวงอุตสาหกรรม. ประกาศกระทรวงอุตสาหกรรม พ.ศ.2549 เรื่องกำหนดค่าปริมาณเขม่าควันที่เจือปนในอากาศที่ระบายออกจากปล่องของหม้อน้ำโรงสีข้าวที่ใช้แก๊สเป็นเชื้อเพลิง. ราชกิจจานุเบกษา. 4 ธันวาคม 2549. เล่มที่ 123 ตอนพิเศษ 125 ง.
2. กระทรวงอุตสาหกรรม. ประกาศกระทรวงอุตสาหกรรม พ.ศ.2548 เรื่อง การกำจัดสิ่งปฏิกูลหรือวัสดุที่ไม่ใช้แล้ว. ราชกิจจานุเบกษา. 25 มกราคม 2549. เล่มที่ 123 ตอนพิเศษ 113.
3. สมาคมวิศวกรรมสิ่งแวดล้อมแห่งประเทศไทย. คู่มือวิเคราะห์น้ำเสีย. พิมพ์ครั้งที่ 4. กรุงเทพฯ: เรือนแก้วการพิมพ์, 2547.
4. APHA, AWWA, WEF. *Standard Methods for the Examination of Water and Wastewater*. 23<sup>rd</sup> ed. Washington, DC : APHA, 2017
5. United States Environmental Protection Agency. *Standard of Performance for New Stationary Sources*. 40 CFR Part 60. Appendix A, 2017.
6. United States Environmental Protection Agency. *Standard of Performance for New Stationary Sources*. 40 CFR Part 60. Appendix A, 2019.
7. United States Environmental Protection Agency. *Standard of Performance for New Stationary Sources*. 40 CFR Part 60. Appendix A, 2020.
8. United States Environmental Protection Agency. *Standard of Performance for New Stationary Sources*. 40 CFR Part 60. Appendix A, 2023.
9. United States Environmental Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Acid Digestion of Sediments Sludge and Soils. SW-846 Method 3050B*, 1996.
10. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Inductively Coupled Plasma-Atomic Emission spectrometry. SW-846 Method 6010C*, 2007.
11. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Mercury in Solid or Semisolid Waste (Manual Cold-Vapor Technique). SW-846 Method 7471B*, 2007.
12. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Alkaline digestion for Hexavalent Chromium. SW-846 Method 3060A*, 1996.

COPY

13 United..

13. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Chromium. Hexavalent (Colorimetric). SW-846 Method 7196A*, 1992

14. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Closed-System Purge-and-Trap and Extraction for Volatile Organics in Soil and Waste Samples. SW-846 Method 5035A*, 2002

15. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Ultrasonic Extraction. SW-846 Method 3550C*, 2007

16. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS). SW-846 Method 8260D*, 2018

17. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. *Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry. SW-846 Method 8270E*, 2018



COPY

ศูนย์วิจัยและเตือนภัยมลพิษโรงงานภาคตะวันออก กรมโรงงานอุตสาหกรรม โทร. ๐ ๓๓๓๓ ๖๐๕๔ ต่อ ๕๐๐๓-๖

ที่ อก ๐๓๒๐/ ๕๖๐๕ 1



กรมโรงงานอุตสาหกรรม  
ถนนพระรามที่ ๖ แขวงทุ่งพญาไท  
เขตราชเทวี กรุงเทพฯ ๑๐๕๐๐

๑๕ พฤษภาคม ๒๕๖๗

เรื่อง เปลี่ยนแปลงเอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษ เปลี่ยนแปลงสารมลพิษในดิน และเปลี่ยนแปลงบุคลากร  
ของห้องปฏิบัติการวิเคราะห์

เรียน กรรมการผู้จัดการ บริษัท อีสเทิร์น ไทย คอนซัลติ้ง ๑๙๙๒ จำกัด

อ้างถึง คำขอขึ้นทะเบียน/ต่ออายุ/เปลี่ยนแปลงบุคลากร และชนิดสารมลพิษของห้องปฏิบัติการวิเคราะห์เอกชน  
ลงวันที่ ๕ มีนาคม ๒๕๖๗

สิ่งที่ส่งมาด้วย เอกสารแนบท้ายหนังสือเปลี่ยนแปลงเอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษ และเปลี่ยนแปลง  
สารมลพิษบริษัท อีสเทิร์น ไทย คอนซัลติ้ง ๑๙๙๒ จำกัด จำนวน ๑๒ แผ่น

ตามคำขอฯ ที่อ้างถึง บริษัท อีสเทิร์น ไทย คอนซัลติ้ง ๑๙๙๒ จำกัด ห้องปฏิบัติการวิเคราะห์  
เอกชน เลขทะเบียน ว-๐๐๓ สถานที่ตั้งเลขที่ ๖๘๓ หมู่ที่ ๑๑ ถนนสุขาภิบาล ๘ ตำบลหนองขาม  
อำเภอศรีราชา จังหวัดชลบุรี แจ้งขอเปลี่ยนแปลงเอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษในน้ำเสีย น้ำใต้ดิน  
เปลี่ยนแปลงสารมลพิษในดิน และเปลี่ยนแปลงบุคลากร นั้น

กรมโรงงานอุตสาหกรรมพิจารณาแล้ว มีความเห็นดังนี้

๑. ให้ยกเลิกผู้ควบคุมดูแลห้องปฏิบัติการวิเคราะห์ จำนวน ๑ ราย

นายวัฒนา โคตรหล้า ทะเบียนเลขที่ ว-๐๐๓-ค-๐๐๐๒

๒. ให้ยกเลิกเจ้าหน้าที่ประจำห้องปฏิบัติการวิเคราะห์ จำนวน ๓ ราย

๑) นางสาวอัญชลี ทะพงษ์ ทะเบียนเลขที่ ว-๐๐๓-จ-๐๐๑๒

๒) นางสาวจุฑามาศ เจริญพรหม ทะเบียนเลขที่ ว-๐๐๓-จ-๐๐๑๕

๓) นางสาวณัฐธิชา นนตานอก ทะเบียนเลขที่ ว-๐๐๓-จ-๐๐๒๔

๓. ให้ยกเลิกขอขยายรายการสารมลพิษในน้ำเสีย และน้ำใต้ดินตามรายการเอกสารแนบท้าย  
หนังสือต่ออายุรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์เอกชนที่ อก ๐๓๒๐/๑๓๓๔๖ ลงวันที่ ๒๗ กรกฎาคม ๒๕๖๖

๔. ให้วิเคราะห์สารมลพิษตามขอขยายที่ได้รับขึ้นทะเบียนให้วิเคราะห์ในน้ำเสีย จำนวน ๔๗ รายการ  
และน้ำใต้ดิน จำนวน ๑๑๑ รายการ รวมทั้งสิ้นจำนวน ๑๕๘ รายการ ตามเอกสารแนบท้ายหนังสือเปลี่ยนแปลง  
เอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษ เปลี่ยนแปลงสารมลพิษในดิน และเปลี่ยนแปลงบุคลากร ดังสิ่งที่ส่งมาด้วย

๕. ให้วิเคราะห์สารมลพิษตามขอขยายที่ได้รับขึ้นทะเบียนให้วิเคราะห์เพิ่มเติมในดิน จำนวน  
๑๒ รายการ ตามเอกสารแนบท้ายหนังสือเปลี่ยนแปลงเอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษเปลี่ยนแปลงสารมลพิษ  
ในดิน และเปลี่ยนแปลงบุคลากร ดังสิ่งที่ส่งมาด้วย

อนึ่ง หนังสือ ....

COPY



"อุตสาหกรรมก้าวไกล ประเทศไทยก้าวหน้า ร่วมกันพัฒนา อุตสาหกรรมสีเขียว"



-๒-

อนึ่ง หนังสือฉบับนี้จะหมดอายุพร้อมหนังสือต่ออายุรับขึ้นทะเบียนห้องปฏิบัติการวิเคราะห์  
เอกชนในวันที่ ๕ กรกฎาคม ๒๕๖๙  
จึงเรียนมาเพื่อทราบ

ขอแสดงความนับถือ

10

(นายพรยศ กลั่นกรอง)  
รองอธิบดี ปฏิบัติราชการแทน  
อธิบดีกรมโรงงานอุตสาหกรรม

ศูนย์วิจัยและเตือนภัยมลพิษโรงงานภาคตะวันออก

โทร. ๐ ๓๓๓๓ ๖๐๕๕ ต่อ ๕๐๐๑-๒

ไปรษณีย์อิเล็กทรอนิกส์ sirwad@w.mae.go.th

COPY



"อุตสาหกรรมก้าวไกล ประเทศไทยก้าวหน้า ร่วมกันพัฒนา อุตสาหกรรมสีเขียว"





เอกสารแนบท้ายหนังสือเปลี่ยนแปลงเอกสารอ้างอิงวิธีวิเคราะห์สารมลพิษ

บริษัท อีสเทิร์น ไทย คอนซัลติง จำกัด

เลขทะเบียน 7-003

ที่ กก ๐๓๒๐/

ลงวันที่

ขอขยายสารมลพิษที่ได้รับขึ้นทะเบียนจากกรมโรงงานอุตสาหกรรม จำนวน ๑๗๐ รายการ

น้ำเสีย จำนวน 47 รายการ

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
1	Aldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
2	Arsenic	1) Continuous Hydride Generation/Atomic Absorption Spectrometric Method <sup>(1)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
3	Barium	Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
4	α-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
5	β-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
6	δ-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
7	γ-BHC	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
8	Biochemical Oxygen Demand	1) 5-Day BOD Test, Membrane Electrode Method <sup>(1)</sup> 2) 5-Day BOD Test, Azide Modification Method <sup>(1)</sup>
9	Cadmium	Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
10	Chemical Oxygen Demand	Closed Reflux, Titrimetric Method <sup>(1)</sup>
11	cis-Chlordane	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>

COPY 12 trans-Chlordane ...

-๒-

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
12	trans-Chlordane	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
13	Chromium	1) Digestion, Direct Air-Acetylene Flame Method <sup>(1)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
14	Color	ADMI Weighted-Ordinate Spectrophotometric Method <sup>(1)</sup>
15	Copper	1) Digestion, Direct Air-Acetylene Flame Method <sup>(1)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
16	Cyanide	Distillation, Colorimetric Method <sup>(1)</sup>
17	4,4'-DDD	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
18	4,4'-DDE	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
19	DDT	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
20	Dieldrin	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
21	Endosulfan I	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
22	Endosulfan II	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
23	Endosulfan sulfate	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>(1)</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>
24	Endrin	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>(1)</sup>

COPY 25 Endrin aldehyde ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
25	Endrin aldehyde	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[1]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[1]</sup>
26	Endrin ketone	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[1]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[1]</sup>
27	Formaldehyde	Distillation, Colorimetric Method <sup>[4]</sup>
28	Free Chlorine	1) Iodometric Method <sup>[1]</sup> 2) Colorimetric Method <sup>[1]</sup>
29	Heptachlor	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[1]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[1]</sup>
30	Heptachlor Epoxide	1) Liquid-Liquid Extraction, Gas Chromatographic Method <sup>[1]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[1]</sup>
31	Hexavalent Chromium	Filtration, Colorimetric Method <sup>[1]</sup>
32	Lead	1) Digestion, Direct Air-Acetylene Flame Method <sup>[1]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
33	Manganese	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
34	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>[1]</sup>
35	Methoxychlor	Liquid-Liquid Extraction, Gas Chromatographic / Mass Spectrometric Method <sup>[1]</sup>
36	Nickel	1) Digestion, Direct Air-Acetylene Flame Method <sup>[1]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
37	Oil and Grease	Liquid-Liquid, Partition-Gravimetric Method <sup>[1]</sup>
38	pH	Electrometric Method <sup>[1]</sup>
39	Phenols	Distillation, Direct Photometric Method <sup>[1]</sup>
40	Selenium	Digestion, Hydride Generation/Atomic Absorption Spectrometric Method <sup>[1]</sup>

**COPY**

41 Sulfide ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
41	Sulfide	ZnS Precipitation, Iodometric Method <sup>[1]</sup>
42	Temperature	Field Method <sup>[1]</sup>
43	Trivalent Chromium	1) Digestion, Direct Air-Acetylene Flame Method; Filtration, Colorimetric Method; Calculation <sup>[1]</sup> 2) Digestion, Inductively Coupled Plasma Method; Colorimetric Method; Calculation <sup>[1]</sup>
44	Total Dissolved Solids	Dried at 180 °C <sup>[1]</sup>
45	Total Kjeldahl Nitrogen	Macro Kjeldahl Method <sup>[1]</sup>
46	Total Suspended Solids	Dried at 103-105 °C <sup>[1]</sup>
47	Zinc	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>

**น้ำใต้ดิน จำนวน 111 รายการ**

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
1	Acenaphthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
2	Acetone	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
3	Aldrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
4	Anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
5	Antimony	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
6	Arsenic	1) Continuous Hydride Generation/Atomic Absorption Spectrometric Method <sup>[1]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
7	Barium	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
8	Benz(a)anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>

**COPY**

9 Benzene ...

-๕-

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
9	Benzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
10	Benzo(b)fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
11	Benzo(k)fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
12	Benzo(a)pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
13	Benzo(g,h,i)perylene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
14	Beryllium	Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
15	Bis(2-chloroethyl)ether	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
16	Bis(2-ethylhexyl)phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
17	Bromodichloromethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
18	Bromoform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
19	Butanol	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
20	Butyl benzyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
21	Cadmium	Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
22	Carbazole	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
23	Carbon disulfide	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
24	Carbon tetrachloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>

**COPY**

25 Chlordane ...

-๖-

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
25	Chlordane	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
26	p-Chloroaniline	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
27	Chlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
28	Chlorodibromomethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
29	Chloroform	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
30	2-Chlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
31	Chromium	1) Digestion, Direct Air-Acetylene Flame Method <sup>(1)</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>(1)</sup>
32	Chromium (III)	1) Digestion, Direct Air-Acetylene Flame Method; Filtration, Colorimetric Method; Calculation <sup>(1)</sup> 2) Digestion, Inductively Coupled Plasma Method; Filtration, Colorimetric Method; Calculation <sup>(1)</sup>
33	Chromium (VI)	Filtration, Colorimetric Method <sup>(1)</sup>
34	Chrysene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
35	Cyanide	Distillation, Colorimetric Method <sup>(1)</sup>
36	DDD	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
37	DDE	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
38	DDT	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
39	Dibenz(a,h)anthracene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>

**COPY**

40 Di-n-butyl phthalate ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
40	Di-n-butyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
41	1,2-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
42	1,3-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
43	1,4-Dichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
44	1,1-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
45	1,2-Dichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
46	1,1-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
47	cis-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
48	trans-1,2-Dichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
49	2,4-Dichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
50	1,2-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
51	1,3-Dichloropropane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
52	Dieldrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
53	Diethyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
54	2,4-Dimethylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>

COPY

55 2,4-Dinitrotoluene ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
55	2,4-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
56	2,6-Dinitrotoluene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
57	Di-n-octyl phthalate	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
58	Endosulfan	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
59	Endrin	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
60	Ethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
61	Fluoranthene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
62	Fluorene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
63	Heptachlor	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
64	Heptachlor epoxide	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
65	Hexachlorobenzene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
66	Hexachloro-1,3-butadiene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
67	n-Hexane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
68	$\alpha$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>
69	$\beta$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>(1)</sup>

COPY

70  $\gamma$ -HCH ...



ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
70	$\gamma$ -HCH	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
71	Hexachlorocyclopentadiene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
72	Hexachloroethane	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
73	Indeno(1,2,3-cd)pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
74	Isophorone	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
75	Lead	1) Digestion, Direct Air-Acetylene Flame Method <sup>[1]</sup> 2) Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
76	Manganese	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
77	Mercury	Digestion, Cold-Vapor Atomic Absorption Spectrometric Method <sup>[1]</sup>
78	Methoxychlor	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
79	Methylene chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
80	2-Methylphenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
81	2-Methylnaphthalene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
82	Methyl tert-butyl ether	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
83	Naphthalene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
84	Nickel	Digestion, Inductively Coupled Plasma Method <sup>[1]</sup>
85	Nitrobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>

86 N-Nitrosodi-n-propylamine ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
86	N-Nitrosodi-n-propylamine	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[1]</sup>
87	pH	Electrometric Method <sup>[4]</sup>
88	Phenanthrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
89	Phenol	1) Distillation, Direct Photometric Method <sup>[4]</sup> 2) Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
90	Pyrene	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
91	Selenium	Digestion, Hydride Generation/Atomic Absorption Spectrometric Method <sup>[4]</sup>
92	Silver	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
93	Styrene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
94	1,1,2,2-Tetrachloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
95	Trichloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
96	1,2,4-Trichlorobenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
97	1,1,1-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
98	1,1,2-Trichloroethane	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
99	Tetrachloroethylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
100	Toluene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
101	2,4,5-Trichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>

102 2,4,6-Trichlorophenol ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
102	2,4,6-Trichlorophenol	Liquid-Liquid Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
103	1,3,5-Trimethylbenzene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
104	Vanadium	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>
105	Vinyl acetate	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
106	Vinyl chloride	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
107	m-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
108	o-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
109	p-Xylene	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
110	Xylene (Total)	Purge and Trap, Gas Chromatographic/Mass Spectrometric Method <sup>[4]</sup>
111	Zinc	Digestion, Inductively Coupled Plasma Method <sup>[4]</sup>

**ดิน จำนวน 12 รายการ**

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
1	$\alpha$ -HCH	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
2	$\beta$ -HCH	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
3	$\gamma$ -HCH	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
4	Heptachlor	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>

**COPY** 5 Aldrin ...

ลำดับ ที่	สารมลพิษ	วิธีวิเคราะห์
5	Aldrin	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
6	Heptachlor epoxide	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
7	Chlordane	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
8	Dieldrin	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
9	Endrin	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
10	DDD	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
11	DDT	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>
12	Methoxychlor	Ultrasonic Extraction, Gas Chromatographic/Mass Spectrometric Method <sup>[2,3]</sup>

**เอกสารอ้างอิง**

1. APHA, AWWA, WEF. **Standard Methods for the Examination of Water and Wastewater**. 24<sup>th</sup> ed. Washington, DC: APHA, 2023.
2. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. **Ultrasonic Extraction**. SW-846 Method 3550C, 2007
3. United States Environment Protection Agency. Test Methods for Evaluation Solid Waste Physical/Chemical Methods. **Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry**. SW-846 Method 8270E, 2018
4. สมาคมวิศวกรรมสิ่งแวดล้อมแห่งประเทศไทย. คู่มือวิเคราะห์น้ำเสีย. พิมพ์ครั้งที่ 4. กรุงเทพฯ: เรือนแก้วการพิมพ์, 2547.

**COPY**

**ANALYTICAL BALANCE (DU)**

**Model : XS205DU**

**Serial No. : 1126323724**

Certificate No. : 23-148799  
Sample Code : 23-56200-001

## CERTIFICATE OF CALIBRATION

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapibarn 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Analytical Balance Room)

Equipment : ELECTRONIC BALANCE

Manufacturer : METTLER TOLEDO

Model : XS205DU

Serial No. : 1126323724

ID No. : LABE 05/1

Date of Receipt : 22 December 2023

Date of Calibration : 22 December 2023

Calibrated by Mr. Somwang Sangdee  
Scientist

Approved by (Mr. Somchai Neampunt)  
Signed for Director

Issue date 25 December 2023

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).

Certificate No. : 23-148799  
Sample Code : 23-56200-001

## REPORT OF CALIBRATION

Equipment : ELECTRONIC BALANCE  
Manufacturer : METTLER TOLEDO  
Model : XS205DU  
Capacity : Max 81 g / 220 g  
Resolution : 0.01 mg / 0.1 mg  
Serial No. : 1126323724  
ID No. : LABE 05/1

### Result of Calibration

#### 1. Test weight and repeatability of reading

Repeatability is a measure of the ability of a balance to supply the same result in repetitive weighings with one and the same load under the same measurement condition. The measurement of the repeatability must include both the balance specifications and the ambient (vibration, fluctuating air current/temperature/humidity, etc.) Operator handling of the balance is also included in the standard deviation.

Unit : g	Range : 80	<input checked="" type="checkbox"/> Before adjustment	<input checked="" type="checkbox"/> After adjustment
<input type="checkbox"/> No adjustment	Nominal value	40	80
<input checked="" type="checkbox"/> Adjustment	Standard weight	40.000054	80.000048
	Average reading of indicator	40.000026	80.000037
	Standard deviation	0.000015	0.000016
		0.000008	0.000009

Unit : g	Range : 200	<input checked="" type="checkbox"/> Before adjustment	<input checked="" type="checkbox"/> After adjustment
<input type="checkbox"/> No adjustment	Nominal value	100	200
<input checked="" type="checkbox"/> Adjustment	Standard weight	100.000042	200.000041
	Average reading of indicator	100.00003	200.00004
	Standard deviation	0.00005	0.00005
		0.00003	0.00005

**COPY**



Certificate No. : 23-148799  
Sample Code : 23-56200-001

Page 3 of 4

## REPORT OF CALIBRATION

## Result of Calibration

## 2. Sensitivity or value of a scale division

Change in the output variable of a measuring instrument divided by the associated change in the input variable.

Unit : g

Range : 80

Range : 200

Test Point	Sensitivity, S	Test Point	Sensitivity, S
0	1.00748	0	1.0274
40	0.98753	100	0.9975
80	0.99751	200	0.9975

## 3. Departure of indication from nominal value, Linearity

Unit : g

Nominal Value	Standard Value	Average Reading of Indicator	Correction Value	Expanded Uncertainty	Coverage Factor (k)
Unload	0.0000000	0.00000	0.00000	0.000012	2.05
0.01	0.0100025	0.01000	0.00000	0.000012	2.05
0.1	0.1000019	0.10001	-0.00001	0.000013	2.03
1	1.0000125	1.00001	0.00000	0.000015	2.02
5	5.0000208	5.00004	-0.00002	0.000021	2.00
10	10.000004	10.00008	-0.00008	0.000026	2.00
20	20.000030	20.00011	-0.00008	0.000036	2.00
50	50.000014	50.00014	-0.00013	0.000068	2.00
100	100.000042	100.0001	-0.0001	0.00016	2.00
150	150.000056	150.0001	0.0000	0.00022	2.00
200	200.000041	200.0002	-0.0002	0.00027	2.00

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.

Certificate No. : 23-148799  
Sample Code : 23-56200-001

Page 4 of 4


## REPORT OF CALIBRATION

## Result of Calibration :

## 4. Eccentric or off-centre loading

Deviation of the measurement value through off-center (eccentric) loading. The corner load increases with the weight of the load and its removal from the center of the pan support.

Weighting pan	<input type="radio"/> Circle <input type="radio"/> Triangular <input checked="" type="radio"/> Rectangular	Test weight : 50 and 100 Unit : g
Range	80	200
Position	Reading of indicator	Reading of indicator
1	50.00015	100.0001
2	50.00022	100.0001
3	50.00008	100.0001
4	50.00002	100.0000
5	50.00016	100.0002
6	50.00014	100.0001
Maximum difference	0.00013	0.0001



## Condition of Calibration

1. Calibration Method : WI-DL-004 base on UKAS LAB 14: 2019

2. This result of calibration was found accurate as shown on date and place of calibration only.

3. Condition of Calibration tem: Normal

4. This certification is traceable to the International System of Unit maintained at :

Through the reference standard laboratory of Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (Instrument number 1).

## 5. Reference standard instrument :

Instrument	Class	ID No.	Certificate No.	Due Date
1) STANDARD WEIGHT 1 mg to 1 kg	E2	LB-WE-79	23-105642	10 September 2024

End of Report

6. Ambient conditions	Min	Max
Temperature (°C)	22.8	23.0
Relative Humidity (%Rh)	43.5	51.1
Air pressure (hPa)	1012.5	1014.5

**ANALYTICAL BALANCE**

**Model : SECURA224-1S**

**Serial No. : 0036707137**



Certificate No. : 23-148800  
Sample Code : 23-56200-002

## CERTIFICATE OF CALIBRATION

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapibarn 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Analytical Balance Room)

Equipment : ELECTRONIC BALANCE

Manufacturer : SARTORIUS

Model : SECURA224-1S

Serial No. : 0036707137

ID No. : LABE 05/2

Date of Receipt : 22 December 2023

Date of Calibration : 22 December 2023

Calibrated by Mr. Somwang Sangdee  
Scientist

Approved by (Mr. Somchai Neampunt)  
Signed for Director

Issue date : 25 December 2023

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



Certificate No. : 23-148800  
Sample Code : 23-56200-002

## REPORT OF CALIBRATION

Equipment : ELECTRONIC BALANCE  
Manufacturer : SARTORIUS  
Model : SECURA224-1S  
Capacity : Max 220 g  
Resolution : 0.0001 g  
Serial No. : 0036707137  
ID No. : LABE 05/2


## Result of Calibration

## 1. Test weight and repeatability of reading

Repeatability is a measure of the ability of a balance to supply the same result in repetitive weighings with one and the same load under the same measurement condition. The measurement of the repeatability must include both the balance specifications and the ambient (vibration, fluctuating air current/temperature/humidity, etc.) Operator handling of the balance is also included in the standard deviation.

Unit : g	Range : 220	<input checked="" type="checkbox"/> Before adjustment	<input checked="" type="checkbox"/> After adjustment
	Nominal value	100	200
<input type="checkbox"/> No adjustment	Standard weight	100.000042	200.000041
<input checked="" type="checkbox"/> Adjustment	Average reading of indicator	99.9998	199.9998
	Standard deviation	0.00006	0.00007

Unit :	Range :	<input type="checkbox"/> Before adjustment	<input type="checkbox"/> After adjustment
	Nominal value		
<input type="checkbox"/> No adjustment	Standard weight		
<input type="checkbox"/> Adjustment	Average reading of indicator		
	Standard deviation		



Certificate No. : 23-148800  
Sample Code : 23-56200-002

## REPORT OF CALIBRATION

## Result of Calibration

## 2. Sensitivity or value of a scale division

Change in the output variable of a measuring instrument divided by the associated change in the input variable.

Unit : g

Range : 220

Range : -

Test Point	Sensitivity, S	Test Point	Sensitivity, S
0	0.7980	-	-
100	0.8978	-	-
200	0.8978	-	-

## 3. Departure of indication from nominal value, Linearity

Unit : g

Nominal Value	Standard Value	Average Reading of Indicator	Correction Value	Expanded Uncertainty	Coverage Factor (k)
Unload	0.0000000	0.0000	0.0000	0.000086	2.00
0.01	0.0100025	0.0100	0.0000	0.000086	2.00
0.1	0.1000019	0.1000	0.0000	0.000087	2.00
1	1.0000125	1.0000	0.0000	0.000087	2.00
2	2.0000089	2.0000	0.0000	0.000087	2.00
5	5.0000208	5.0001	-0.0001	0.000088	2.00
10	10.000004	10.0000	0.0000	0.000090	2.00
20	20.000030	20.0000	0.0000	0.000093	2.00
50	50.000014	50.0000	0.0000	0.00011	2.00
100	100.000042	100.0000	0.0000	0.00016	2.00
200	200.000041	200.0000	0.0000	0.00028	2.00

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.



Certificate No. : 23-148800  
Sample Code : 23-56203-002


## REPORT OF CALIBRATION

## Result of Calibration :

## 4. Eccentric or off-centre loading

Deviation of the measurement value through off-center (eccentric) loading. The corner load increases with the weight of the load and its removal from the center of the pan support.

Weighing pan	<input checked="" type="radio"/> Circle <input type="radio"/> Triangular <input type="radio"/> Rectangular	Test weight : 100 Unit : g
Range	220	
Position	Reading of indicator	Reading of indicator
1	100.0000	-
2	100.0000	-
3	100.0000	-
4	99.9999	-
5	100.0000	-
6	100.0000	-
Maximum difference	0.0001	-



## Condition of Calibration

1. Calibration Method : WI-CL-004 base on UKAS LAB 14: 2019
2. This result of calibration was found accurate as shown on date and place of calibration only.

3. Condition of Calibration item: Normal

4. This certification is traceable to the International System of Unit maintained at : \*

\* Through the reference standard laboratory of Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (Instrument number 1)

5. Reference standard instrument :

Instrument	Class	ID No.	Certificate No.	Due Date
1) STANDARD WEIGHT 1 mg to 1 kg	E2	LB-WE-79	23-105642	10 September 2024

End of Report

6. Ambient conditions	Min	Max
Temperature (°C)	24.4	24.8
Relative Humidity (%Rh)	39.9	41.1
Air pressure (hPa)	1012.2	1012.8



**AUTOCLAVE**

**Model : FLA-1000**

**Serial No. : 55169083**

## CERTIFICATE OF CALIBRATION

Certificate No. : 24-089294

Sample Code : 24-35676-004

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
663 Moo 11, Sukhapibarn 8 Rd., Nongkham, Sriracha,  
Chonburi 20230

Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Autoclave Room)

Equipment : Autoclave  
Manufacturer : TOMY Model : FLS-1000  
Serial No. : 55169083 ID No. : LABE 43/2  
Date of Receipt : 16 July 2024 Date of Calibration : 16 July 2024

## Condition of Calibration

1. Environment  
1.1 Ambient temperature : Maximum 31.3 °C ; Minimum 30.4 °C  
1.2 Relative humidity : Maximum 67.7 % ; Minimum 65.1 %  
1.3 Line voltage supplied : Maximum 219.9 VAC ; Minimum 216.8 VAC

## 2. Calibration method

The calibration use in-house method: WI-CL-025 based on BS 2646-1: 2021

## 3. Reference standard instrument

Instrument	Model	ID No.	Certificate No.	Due Date
3.1 Temperature Data Logger	HiTemp 140	LB-TEM-22	24-030263	18 March 2025
3.2 Temperature Data Logger	HiTemp 140	LB-TEM-23	24-030264	18 March 2025
3.3 Temperature Data Logger	HiTemp 140	LB-TEM-24	24-030265	18 March 2025

## 4. This certificate is traceable to the international system of unit (SI Unit).

The measurement is traceable to Asia Medical and Agricultural Laboratory and Research Center Public Company Limited.

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

## 6. Condition of calibration item : Normal

Calibrated by Mr. Pattanapong Pulngern

Scientist

Approved by

(Mr. Somchai Neampunt)

Signed for Director

Issue date 17 July 2024

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC)

## REPORT OF CALIBRATION

Certificate No. : 24-089294

Sample Code : 24-35676-004

## Results of Calibration

Resolution : 1 °C

## 1. Reporting of Temperature

Calibration point (°C)	UUC* setting (°C)	UUC* reading		Measured Temperature at each positions (°C)			Uncertainty ± (°C)	Coverage factor k
		Temperature (°C)	Pressure (kPa)	# 1	# 2 Ref	# 3		
121	121	121	120	121.41	121.30	121.32	1.2	2.00

## 2. Characterization results

Calibration Point (°C)	Stability ± (°C)	Uniformity (°C)	Overall Variation (°C)
121	0.88	0.51	1.92

## Notes

1. UUC\* = Unit Under Calibration
2. The quoted uncertainty includes "Stability of chamber and loading effect in chamber at 20% of uniformity".
3. Uniformity - the maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time.
4. Stability - one-half of the greatest maximum difference of measured temperatures at any one sensor.
5. Overall variation - the difference of the maximum and the minimum measured temperatures throughout observation time.
6. UUC\* reading - the average reading of indicating device that forms the integral part of the autoclave.
7. Calibration results without adjustment.




Figure: Example of sensor  
installation Positions

- Standard 1 : In the chamber drain, within 100 mm  
Standard 2 : In the upper half of the chamber.  
Standard 3 : Attached to the load temperature probe, within 15 mm

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003

- End of Report -

## **BOD INCUBATOR**

**Model : LABE 19/3**



Page 1 of 3

## CERTIFICATE OF CALIBRATION

Certificate No. : 24-089291

Sample Code : 24-35676-001

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapibarn 8 Rd., Nongkham, Sriracha,  
Chonburi 20230

Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Laboratory)

Equipment : Temperature controlled enclosures (Incubator)

Manufacturer : พัด เครื่องเย็น Model : N/A

Serial No. : S43020027 ID No. : LABE 19/3

Date of Receipt : 16 July 2024 Date of Calibration : 16 July 2024

## Condition of Calibration

1. Environment
- |                           |   |         |           |   |         |           |
|---------------------------|---|---------|-----------|---|---------|-----------|
| 1.1 Ambient temperature   | : | Maximum | 30.6 °C   | : | Minimum | 28.9 °C   |
| 1.2 Relative humidity     | : | Maximum | 76.9 %    | : | Minimum | 69.4 %    |
| 1.3 Line voltage supplied | : | Maximum | 219.8 VAC | : | Minimum | 217.1 VAC |

## 2. Calibration method

TLAS-G-20: Guidelines for calibration and checks of temperature controlled enclosures.

## 3. Reference standard instrument

Instrument	ID No.	Certificate No.	Due Date
Data acquisition with sensor (RTD-P1100)	LB-DA-12 (RTD-168 to RTD-176)	24-045389	28 April 2025

## 4. This certificate is traceable to the international system of unit (SI Unit).

The measurement is traceable to Asia Medical and Agricultural Laboratory and Research Center Public Company Limited.

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

## 6. Condition of calibration item : Normal

Calibrated by : Mr. Pattanapong Pulngern  
Scientist

Approved by : (Mr. Somchai Neampunt)  
Signed for Director

Issue date : 17 July 2024

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is Issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



Page 2 of 3

## REPORT OF CALIBRATION

Certificate No. : 24-089291

Sample Code : 24-35676-001

## Results of Calibration

Resolution : 0.1 °C

## 1. Reporting of Temperature

Calibration point (°C)	UUC* setting (°C)	UUC* reading (°C)	Measured temperature at each positions (°C)									Uncertainty ± (°C)	Coverage factor k
			# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9 <sup>Rev</sup>		
20	20.0	20.0	20.56	20.45	20.01	19.85	20.21	20.25	20.17	20.05	20.11	0.24	2.00

## 2. Characterization results

Calibration point (°C)	Stability ± (°C)	Uniformity (°C)	Overall variation (°C)
20	0.08	0.50	0.87

## Notes

- UUC\* = Unit Under Calibration



## REPORT OF CALIBRATION

Certificate No. : 24-089291

Sample Code : 24-35676-001

## Results of Calibration

## Notes

1. Sensor installation locations
  - 1.1 All sensors at any corners or walls should be positioned 5 cm (a x b x c) from the wall.
  - 1.2 The reference sensor is preferably located of the geometric center of the chamber.
2. Interior dimensions approx of chamber :  
W = 70 cm ; D = 55 cm ; H = 140 cm
3. Air valve or fresh air level : Off
4. Fan level : Open
5. The quoted uncertainty includes "Stability of chamber and loading effect in chamber at 20% of uniformity".
6. Uniformity - the maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time.
7. Stability - one-half of the greatest maximum difference of measured temperatures at any one sensor.
8. Overall variation - the difference of the maximum and the minimum measured temperatures throughout observation time.
9. UUC\* reading - the average reading of indicating device that forms the integral part of the enclosure.
10. Calibration results without adjustment.




Figure: Example of sensor  
installation Positions

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.

- End of Report -

**COPY**



# **BOD INCUBATOR**

**Model : LABE 19/5**

## CERTIFICATE OF CALIBRATION

Certificate No. : 24-046203  
Sample Code : 24-18906-002

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapibarn 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Laboratory)

Equipment : Temperature controlled enclosures (Incubator)

Manufacturer : Lovibond Model : Tc445S

Serial No. : 0520/005227 ID No. : LABE 19/5

Date of Receipt : 18 April 2024 Date of Calibration : 18 April 2024

## Condition of Calibration

1. Environment
- |                           |           |           |           |           |
|---------------------------|-----------|-----------|-----------|-----------|
| 1.1 Ambient temperature   | : Maximum | 35.0 °C   | : Minimum | 33.7 °C   |
| 1.2 Relative humidity     | : Maximum | 69.1 %    | : Minimum | 50.0 %    |
| 1.3 Line voltage supplied | : Maximum | 222.5 VAC | : Minimum | 218.8 VAC |

## 2. Calibration method

TLAS-G-20: Guidelines for calibration and checks of temperature controlled enclosures.

## 3. Reference standard instrument

Instrument	ID No.	Certificate No.	Due Date
Data Acquisition With Sensor (RTD-Pt100)	LB-DA-08 (RTD-248 to RTD-256)	23-084070	06 August 2024

## 4. This certificate is traceable to the international system of unit (SI Unit).


The measurement is traceable to Asia Medical and Agricultural Laboratory and Research Center Public Company Limited.

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

## 6. Condition of calibration item : Normal

Calibrated by Mr. Sarawoot Thammo  
Scientist  
Issue date 19 April 2024

Approved by

  
(Mr. Nuttaput Timula)  
Signed for Director

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).

## REPORT OF CALIBRATION

Certificate No. : 24-046203  
Sample Code : 24-18906-002

## Results of Calibration

Resolution : 0.1 °C

## 1. Reporting of Temperature

Calibration point (°C)	UUC* setting (°C)	UUC* reading (°C)	Measured temperature at each positions (°C)									Uncertainty ± (°C)	Coverage factor k
			# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9 <sup>Ref</sup>		
20	20.5	20.0	20.28	19.86	19.90	19.91	19.82	20.10	20.01	19.89	19.75	0.59	2.00

## 2. Characterization results


Calibration point (°C)	Stability ± (°C)	Uniformity (°C)	Overall variation (°C)
20	0.45	0.85	1.31

## Notes

- UUC\* = Unit Under Calibration

Calibrated by Mr. Sarawoot Thammo  
Scientist  
Issue date 19 April 2024

Approved by

  
(Mr. Nuttaput Timula)  
Signed for Director

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).

## REPORT OF CALIBRATION

Page 3 of 3

Certificate No. : 24-046203

Sample Code : 24-18906-002

### Results of Calibration

#### Notes

1. Sensor installation locations
  - 1.1 All sensors at any corners or walls should be positioned 5 cm (a x b x c) from the wall.
  - 1.2 The reference sensor is preferably located of the geometric center of the chamber.
2. Interior dimensions approx of chamber :  
W = 60 cm ; D = 56 cm ; H = 146 cm
3. Air valve or fresh air level : Off
4. Fan level : Open
5. The quoted uncertainty includes "Stability of chamber and loading effect in chamber at 20% of uniformity".
6. Uniformity - the maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time.
7. Stability - one-half of the greatest maximum difference of measured temperatures at any one sensor.
8. Overall variation - the difference of the maximum and the minimum measured temperatures throughout observation time.
9. UUC\* reading - the average reading of indicating device that forms the integral part of the enclosure.
10. Calibration results without adjustment.




Figure: Example of sensor  
installation Positions

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.

End of Report

**COPY**

**DO**

**Model : YSI 5000**

**Serial No. : 18E101961**



Harikul Science Co.,Ltd.  
694 Soi Ratchadanivet 24, Pracharatbamphen,  
Samsaennok, Huaikhwang, Bangkok 10310  
Tel: 0-2274-2456 Fax: 0-2274-2443  
Email: info@harikul.com www.harikul.com  
Certificate of Calibration

CERT.No.: HS-V053H

Calibration Date : 13 Aug 24	Model : YSI 5000
Submitted by : Eastern Thai Consulting 1992 Company Limited	S/N : 18E101961
683 Moo.11 Sukaphibal8 Rd., Nongkham, Sriracha,	Probe : YSI 5010
Chonburi 20230	S/N : 18A100724
	ID NO. : -
Avg Room Temp : 20 °C	Air Temp ref : S/N. F8065C26
Avg Water Temp : 20 °C	Barometric ref : S/N. F8065C26
Air Pressure : 760.00 mmHg	Water Temp ref : -
Salinity : 0 ppt	ID NO. HS001
	Technician : Kittipong M.

#### Calibration Details

Calibration Point	100% air sat. (@20 °C, DO = 9.09 mg/l)	(status)	(status)
Measurement 1 (mg/l)	9.08	(PASS)	-
Measurement 2 (mg/l)	9.08	(PASS)	-
Measurement 3 (mg/l)	9.09	(PASS)	-
Measurement 4 (mg/l)	9.09	(PASS)	-
Measurement 5 (mg/l)	9.09	(PASS)	-
Measurement 6 (mg/l)	9.09	(PASS)	-
Measurement 7 (mg/l)	9.07	(PASS)	-
Measurement 8 (mg/l)	9.07	(PASS)	-
Measurement 9 (mg/l)	9.07	(PASS)	-
Measurement 10 (mg/l)	9.07	(PASS)	-

Mean Measurement	9.08	mg/l	-	-
Inaccuracy	0.01	mg/l	-	-

Overall Status (PASS)

#### Manufacturer Specification

Accuracy = +/- 0.02 mg/l

- 1) This certificate is issued based on the result that are found as shown on date and place of test only.
- 2) The calibration procedure followed in accordance with Harikul Science Co., Ltd.
- 3) This result shall not be used for advertising purpose.

Technician Signature  
(Kittipong Maekwong)

Laboratory Manager  
(Natenapha Pisatkunchon)

COPY



**Hot Air Oven**

**Model : UM 400**

**Serial No. : 900982**

## CERTIFICATE OF CALIBRATION

Certificate No. : 24-001944  
Sample Code : 24-00963-001Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapibarn 8 Rd., Nongkham,  
Sriracha, Chonburi 20230Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Hot Lab)

Equipment : Temperature controlled enclosures (Hot air oven)

Manufacturer : Memmert Model : UM 400

Serial No. : 900982 ID No. : LABE 17/1

Date of Receipt : 09 January 2024 Date of Calibration : 09 January 2024

## Condition of Calibration

1. Environment
- |                           |           |           |           |           |
|---------------------------|-----------|-----------|-----------|-----------|
| 1.1 Ambient temperature   | : Maximum | 30.6 °C   | : Minimum | 29.2 °C   |
| 1.2 Relative humidity     | : Maximum | 57.5 %    | : Minimum | 46.4 %    |
| 1.3 Line voltage supplied | : Maximum | 229.5 VAC | : Minimum | 222.5 VAC |

## 2. Calibration method

TLAS-G-20: Guidelines for calibration and checks of temperature controlled enclosures.

## 3. Reference standard instrument

Instrument	ID No.	Certificate No.	Due Date
Data Acquisition With Sensor (RTD-Pt100)	LB-DA-10 (RTD-257 to RTD-265)	23-066256	29 June 2024

## 4. This certificate is traceable to the international system of unit (SI Unit).

The measurement is traceable to Asia Medical and Agricultural Laboratory and Research Center Public Company Limited.

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

## 6. Condition of calibration item : Normal

Calibrated by Mr. Sarawoot Thammo  
Scientist

Approved by

(Mr. Somchai Nearnpunt)

Signed for Director

Issue date 09 January 2024

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).

## REPORT OF CALIBRATION

Certificate No. : 24-001944  
Sample Code : 24-00963-001

## Results of Calibration

Resolution : 0.1 °C

## 1. Reporting of Temperature

Calibration point (°C)	UUC* setting (°C)	UUC* reading (°C)	Measured temperature at each positions (°C)									Uncertainty ± (°C)	Coverage factor k
			# 1	# 2	# 3	# 4	# 5	# 6	# 7	# 8	# 9 <sup>Ref</sup>		
60	60.0	60.0	60.04	59.90	59.81	59.84	59.47	59.91	60.08	59.98	59.87	0.25	2.00
85	85.0	85.0	86.07	85.75	85.58	85.62	84.69	85.83	86.28	85.94	85.77	0.34	2.00

## 2. Characterization results

Calibration point (°C)	Stability ± (°C)	Uniformity (°C)	Overall variation (°C)
60	0.11	0.49	0.80
85	0.09	1.13	1.72

## Notes

- UUC\* = Unit Under Calibration



## REPORT OF CALIBRATION

Certificate No. : 24-001944

Sample Code : 24-00963-001

## Results of Calibration

## Notes

1. Sensor installation locations
  - 1.1 All sensors at any corners or walls should be positioned 5 cm (a x b x c) from the wall.
  - 1.2 The reference sensor is preferably located of the geometric center of the chamber.
2. Interior dimensions approx of chamber :  
W = 40 cm ; D = 28 cm ; H = 39 cm
3. Air valve or fresh air level : Off
4. Fan level : Open
5. The quoted uncertainty includes "Stability of chamber and loading effect in chamber at 20% of uniformity".
6. Uniformity - the maximum difference of measured temperatures at any sensors and the measured temperature at the reference location which are observed at the same time.
7. Stability - one-half of the greatest maximum difference of measured temperatures at any one sensor.
8. Overall variation - the difference of the maximum and the minimum measured temperatures throughout observation time.
9. UUC\* reading - the average reading of indicating device that forms the integral part of the enclosure.
10. Calibration results without adjustment.




Figure: Example of sensor  
installation Positions

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.

- End of Report -



COPY

**LIQUID IN GLASS THERMOMETER**

**Model : Total Immersion**

**Serial No. : 43560**



QUALITY CALIBRATION CO.,LTD.  
235 Petchkasem 63/2 Road, Laksong, Bangkae, Bangkok 10160  
Tel (662) 421-5402, (662) 444-0152-3, Fax (662) 809-4584  
www.qcalibration.com



CERTIFICATE No : 23T10864  
REFERENCE No : 71117-1

PAGE : 1 OF 2

## Certificate of Calibration

EQUIPMENT : LIQUID IN GLASS THERMOMETER  
MANUFACTURER : PRECISION  
MODEL : 0 °C TO 100 °C  
SERIAL No : 43560  
ID No : LABE 16/1  
RESOLUTION : 0.1 °C  
TYPE : TOTAL IMMERSION  
CONDITION AS RECEIVED : USED ITEM  
SUBMITTED BY : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 MOO 11, SUKHAPIBAN 8 ROAD, NONGKHAM,  
SRIRACHA, CHONBURI 20230

CALIBRATED BY : CHARUKIT L.  
CALIBRATION DATE : 09-Nov-23  
APPROVED BY : PONGSAK J.  
ISSUED DATE : 09-Nov-23  
RECEIVED DATE : 02-Nov-23

**COPY**

THIS CERTIFICATE MAY NOT BE REPRODUCED OTHER THAN IN FULL EXCEPT WITH THE PRIOR WRITTEN APPROVAL OF  
QUALITY CALIBRATION CO., LTD.



QUALITY CALIBRATION CO.,LTD.  
235 Petchkasem 63/2 Road, Laksong, Bangkae, Bangkok 10160  
Tel (662) 421-5402, (662) 444-0152-3, Fax (662) 809-4584  
www.qcalibration.com

CERTIFICATE No : 23T10864

PAGE : 2 OF 2

## Calibration Report

EQUIPMENT : LIQUID IN GLASS THERMOMETER  
MANUFACTURER : PRECISION  
MODEL : 0 °C TO 100 °C  
ID No : LABE 16/1  
RESOLUTION : 0.1 °C  
RECEIVED DATE : 02-Nov-23  
AMBIENT TEMPERATURE : 23 °C ± 3 °C  
SERIAL NUMBER : 43560  
TYPE : TOTAL IMMERSION  
CALIBRATION DATE : 09-Nov-23  
RELATIVE HUMIDITY : 50 %RH ± 20 %RH

### CONDITION OF THIS RESULTS OF CALIBRATION

- THIS INSTRUMENT WAS CALIBRATED BASED ON ASTM E77:1992 BY COMPARISON WITH STANDARD PLATINUM RESISTANCE THERMOMETER (SPRT) INTO LIQUID BATH TEMPERATURE CONTROLLER. THE TEMPERATURE SCALE USED WAS BASED ON ITS-90.
- REFERENCE STANDARD INSTRUMENTS :-

INSTRUMENT	MODEL	SERIAL No	CERTIFICATE No	DUE DATE
1) STANDARD THERMOMETER	1502	77964	23T3927	08-Mar-24
2) SPRT PROBE	5614	636636	23T3927	08-Mar-24
3) PRECISION BATH	7320	A21105	22T13199	14-Dec-23
4) PRECISION BATH	CTR-40	A68155	22T13198	09-Dec-23
5) PRECISION BATH	6045	3C023	22T13200	19-Dec-23
- THE CERTIFICATE IS VALID FOR THE ITEM CALIBRATED AS SHOWN ON THE DATE AND PLACE OF CALIBRATION ONLY.
- THIS RESULT EXCLUDE LONG TERM STABILITY OF THE UNIT UNDER CALIBRATION.
- THIS CERTIFICATE IS TRACEABLE TO THE INTERNATIONAL SYSTEM OF UNIT MAINTAINED AT:-  
- NATIONAL INSTITUTE OF METROLOGY (THAILAND).

### RESULT OF CALIBRATION : WITHOUT ADJUSTMENT

STANDARD READING (°C)	UUC* READING (°C)	IMMERSION DEPTH (mm)	CORRECTION (°C)	EMERGENT STEM TEMPERATURE (°C)	UNCERTAINTY OF MEASUREMENT (±°C)
0.009	0.0	60	0.0090	N/A	0.26
25.01	25.0	165	0.0050	N/A	0.26
50.00	50.0	275	0.0040	N/A	0.26
99.991	100.0	360	-0.009	29.3	0.26

### UUC\* : UNIT UNDER CALIBRATION

THE REPORTED UNCERTAINTY OF MEASUREMENT WAS BASED ON A STANDARD UNCERTAINTY MULTIPLIED BY A COVERAGE FACTOR k=2, PROVIDING A LEVEL OF CONFIDENCE APPROXIMATELY 95%.  
END OF CALIBRATION REPORT

**COPY**



**pH Meter**

**Model : SevenCompact S220**

**Serial No. : B448305208**



## CERTIFICATE OF CALIBRATION

Page 1 of 3

Supersede to Calibration Certificate No. 24-001949

Certificate No. : 24-001949/1

Sample Code : 24-00963-006

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.

683 Moo 11, Sukhapibarn 8 Rd., Nongkham,  
Sriracha, Chonburi 20230Location of Calibration : EASTERN THAI CONSULTING 1992 CO., LTD.  
(Laboratory)

Equipment : pH Meter

Manufacturer : METTLER TOLEDO

Model : SevenCompact S220

Serial No. : B448305208

ID No. : LABE 11/4

Date of Receipt : 09 January 2024

Date of Calibration : 09 January 2024

## Condition of Calibration

## 1. Environment

1.1 Ambient temperature : 22.4 ± 0.2 °C 1.2 Relative humidity : 56.4 % ± 2.1 %

## 2. Calibration method

In house method WI-CL-019: based on direct measurement by using standard voltage calibrator and using certified reference material (CRM).

## 3. Reference standard / Certified reference material

Instrument	ID No.	Certificate No.	Due Date
3.1 Voltage Calibrator	LB-AMC-01	23E3244	03 October 2024
3.2 Digital Thermometer	LB-TH-33	23-098974	25 August 2024
Certified Reference Material	Lot. No.	Ref No.	Expire Date
3.3 Buffer Solution pH 4.008	919273	PH216.L5	24 September 2025
3.4 Buffer Solution pH 6.986	941727	PH107.L5	06 November 2024
3.5 Buffer Solution pH 9.997	919278	PH220.L5	24 September 2024

## 4. This certificate is traceable to the international system of unit (SI Unit).

4.1 Instrument No. 3.1 through Technology Promotion Association (Thailand-Japan).

4.2 Instrument No. 3.2 through Asia Medical and Agricultural Laboratory and Research Center Public Company Limited.

4.3 Buffer Solution No. 3.3 and No. 3.5 traceable to CPA chem (through primary measurement method-Harned cell using calibrated thermometer, barometer, and nanovoltmeter Accredited laboratory ISO/IEC 17025 and ISO/IEC 17034).

4.4 Buffer Solution No. 3.4 traceable to CPA chem (CPA RefN HARNED CELL LotN 61275737; CPA RefN HARNED CELL LotN 61273986 Accredited laboratory ISO/IEC 17025 and ISO/IEC 17034).

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

## 6. Condition of calibration item : Normal

Calibrated by Mr. Nuttaput Timula

Approved by

(Mr. Sornchai Neampunt)

Scientist

Signed for Director

Issue date 31 January 2024

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



## REPORT OF CALIBRATION

Page 2 of 3

Supersede to Calibration Certificate No. 24-001949

Certificate No. : 24-001949/1

Sample Code : 24-00963-006

Equipment : pH Meter Resolution : 0.01 pH ; 0.1 mV ; 0.1 °C  
Manufacturer : METTLER TOLEDO Model : SevenCompact S220  
Serial No. : B448305208 ID No. : LABE 11/4  
Range : -2.000 pH to 20.000 pH ; ± 2000.0 mV ; -5.0 °C to 130.0 °C

## Results of Calibration

## Part 1. DC Voltage measurement

pH Meter Serial No. : B448305208

Nominal Value	Applied DC Voltage	Average indicator reading		Uncertainty	Coverage factor
		mV	pH		
0	414.113	413.9	0.00	± 0.083	2.00
4	177.477	177.4	4.00	± 0.083	2.00
7	0.000	0.1	7.00	± 0.083	2.00
10	-177.477	-177.3	10.00	± 0.083	2.00
14	-414.113	-413.8	14.00	± 0.083	2.00

## Part 2. Performance of Electrode system

Electrode Manufacturer : METTLER TOLEDO Model : InLab Expert Pro-ISM

Electrode Serial No. : 2453982

Three-Point Calibration at pH4, pH7 and pH10 Percent Slope : 98.3

Standard Buffer Solution	Average indicator reading		Error Value	Uncertainty	Coverage factor
	pH	mV			
pH (@ 25 °C)			pH	pH	k
4.008	4.01	182.1	0.002	± 0.010	2.00
6.986	7.00	7.8	0.014	± 0.011	2.00
9.997	10.01	-167.2	0.013	± 0.011	2.00

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003.

NSC-TISI-TIS17025  
CALIBRATION 0152

Page 3 of 3

## REPORT OF CALIBRATION

Supersede to Calibration Certificate No. 24-001949

Certificate No. : 24-001949/1

Sample Code : 24-00963-006

Equipment : pH Meter (Digital Thermometer with sensor)

## Thermometer readout

Manufacturer : METTLER TOLEDO Model : SevenCompact S220  
Serial No. : B448305208 ID No. : LABE 11/4  
Resolution : 0.1 °C Range : -5.0 °C to 130.0 °C

## Thermometer sensor

Manufacturer : METTLER TOLEDO Model : InLab Expert Pro-ISM  
Serial No. : 2453982 ID No. : N/A

## Condition of Calibration

1. Environment 1.1 Ambient temperature : 22.6 °C ± 0.1 °C  
1.2 Relative humidity : 55.1 % ± 3.3 %

## 2. Calibration method

- 2.1 The calibration use in house method WI-CL-021 : by comparison with standard thermometer  
2.2 The calibration by comparison unit under calibration (UUC) to the standard thermometer in a calibration bath at the controlled temperature.  
2.3 The temperature scale in use of this laboratory is the international temperature scale of 1990 (ITS-90).

## 3. Reference standard instrument

Instrument	Model	ID. No.	Certificate No.	Due date
3.1 Resistance Thermometer	PT-100	RTD-90	23-098974	25 August 2024
3.2 Thermometer Readout	GT-11	LB-TH-33	23-098974	25 August 2024

## 4. This certificate is traceable to the international system of unit (SI Unit).

Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (Accreditation Under TLAS Laboratory Calibration No.0152)

## 5. This result of calibration was found accurate as shown on date and place of calibration only.

6. Condition of Calibration item : Normal

## Results of Calibration

Calibration point °C	Average of standard reading °C	Unit under calibration			Expanded uncertainty °C	Coverage factor k
		Immersion depth mm	Average reading °C	Correction value °C		
25	25.000	120	25.0	0.000	± 0.14	2.00

## Notes

- Calibration results without adjustment

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003

- End of report -

**STANDARD WEIGHT 50 g**



Certificate No. : 22-052238  
Sample Code : 22-19150-003

## CERTIFICATE OF CALIBRATION

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapiban 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : Asia Medical and Agricultural Laboratory and Research Center Public Company Limited  
(Calibration Laboratory)

Equipment : Standard Weight 50 g

Manufacturer : METTLER TOLEDO

Class : F1

Serial No. : N/A

ID No. : LABE 10/1

Date of Receipt : 18 May 2022

Date of Calibration : 30 May 2022

Calibrated by : Mr. Somwang Sangdee  
Scientist

Approved by : ( Mr. Somchai Neampunt )  
Signed for Director

Issue date : 31 May 2022

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



Certificate No. : 22-052238  
Sample Code : 22-19150-003

## REPORT OF CALIBRATION

Equipment : Standard Weight 50 g  
Manufacturer : METTLER TOLEDO  
Class : F1  
Serial No. : N/A  
ID No. : LABE 10/1

Result of Calibration : ☒ Without adjustment ☐ Adjustment

Conventional value of the result of weighing in air. For a weight taken at a reference temperature ( $t_{ref}$ ) of 20°C, the conventional mass is the mass of a reference weight of a density ( $\rho_{ref}$ ) of 8000 kg.m<sup>-3</sup> which it balances in air of a reference density ( $\rho_0$ ) of 1.2 kg.m<sup>-3</sup>

Description	Deviation	Conventional	Expanded	Maximum	ID No.
		Mass	Uncertainty	Permissible Error	
	(mg)		(mg)	± (mg)	
50 g	-0.324	49.999676 g	0.10	0.30	LABE 10/1

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k = 2.0$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003

COPY



Certificate No. : 22-052238

Sample Code : 22-19150-003

Page 3 of 3

## REPORT OF CALIBRATION

## Condition of Calibration:

1. Ambient Conditions : Temperature  $20^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$ , Relative humidity  $50\% \pm 10\%$  and air density  $1.20 \text{ kg/m}^3$ 

2. Calibration Method : Direct comparison weighing according to OIML R111-1 : 2004(E)

3. Reference standard instrument

Instrument	Class	ID No.	Certificate No.	Due Date
1) Standard Weight 1 mg to 1 kg	E2	LB-WE-79	21-079366	22 September 2022

4. This certification is traceable to the International System of Unit maintained at : -

Asia Medical and Agricultural Laboratory and Research Center Public Company Limited

(Instrument number 1).

5. Condition of Calibration item: Normal

## 6. Description of Calibrated Item :

Type and Nominal Value :	Standard Weight 50 g
Shape :	Cylindrical weight with knob
Material :	Stainless steel
Case :	Wooden Box
Comments :	Recalibration

- End of Report -

COPY

**STANDARD WEIGHT 100 g**



Certificate No. : 22-052239  
Sample Code : 22-19150-004

## CERTIFICATE OF CALIBRATION

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapiban 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : Asia Medical and Agricultural Laboratory and Research Center Public Company Limited  
(Calibration Laboratory)

Equipment : Standard Weight 100 g

Manufacturer : N/A

Class : N/A

Serial No. : N/A

ID No. : LABE 10/2

Date of Receipt : 18 May 2022

Date of Calibration : 30 May 2022

Calibrated by : Mr. Somwang Sangdee  
Scientist

Approved by

( Mr. Somchai Neampunt )

Signed for Director

Issue date : 31 May 2022

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



Certificate No. : 22-052239  
Sample Code : 22-19150-004

## REPORT OF CALIBRATION

Equipment : Standard Weight 100 g  
Manufacturer : N/A  
Class : N/A  
Serial No. : N/A  
ID No. : LABE 10/2

### Result of Calibration :

☒ Without adjustment

☐ Adjustment

Conventional value of the result of weighing in air. For a weight taken at a reference temperature ( $t_{ref}$ ) of 20°C, the conventional mass is the mass of a reference weight of a density ( $\rho_{ref}$ ) of 8000 kg.m<sup>-3</sup> which it balances in air of a reference density ( $\rho_0$ ) of 1.2 kg.m<sup>-3</sup>

Description	Deviation	Conventional	Expanded	Maximum	ID No.
		Mass	Uncertainty	Permissible Error	
	(mg)		(mg)	± (mg)	
100 g	-0.171	99.999829 g	0.16	0.50	LABE 10/2

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k = 2.0$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003

COPY



Certificate No. : 22-052239

Sample Code : 22-19150-004

Page 3 of 3

## REPORT OF CALIBRATION

## Condition of Calibration

1. Ambient Conditions : Temperature  $20^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$ , Relative humidity  $50\% \pm 10\%$  and air density  $1.18 \text{ kg/m}^3$ 

2. Calibration Method : WI-CL-007 base on OIML R 111-1 : 2004(E)

3. Reference standard instrument

Instrument	Class	ID No.	Certificate No.	Due Date
1) Standard Weight 1 mg to 1 kg	E2	LB-WE-79	21-079366	22 September 2022

4. This certification is traceable to the International System of Unit maintained at : -

Asia Medical and Agricultural Laboratory and Research Center Public Company Limited

( Instrument number 1).

5. Condition of Calibration item: Normal

## 6. Description of Calibrated Item :

Type and Nominal Value :	Standard Weight 100 g
Shape :	Cylindrical weight with knob
Material :	Stainless steel
Case :	Wooden Box
Comments :	Recalibration

- End of Report -

COPY

**STANDARD WEIGHT 50 g**





Certificate No. : 22-052237  
Sample Code : 22-19150-002

## CERTIFICATE OF CALIBRATION

Customer : EASTERN THAI CONSULTING 1992 CO., LTD.  
683 Moo 11, Sukhapiban 8 Rd., Nongkham,  
Sriracha, Chonburi 20230

Location of Calibration : Asia Medical and Agricultural Laboratory and Research Center Public Company Limited  
(Calibration Laboratory)

Equipment : Standard Weight 50 g

Manufacturer : N/A

Class : N/A

Serial No. : N/A

ID No. : LABE 10/4

Date of Receipt : 18 May 2022

Date of Calibration : 30 May 2022

Calibrated by : Mr. Somwang Sangdee  
Scientist

Approved by : ( Mr. Somchai Neampunt )  
Signed for Director

Issue date : 31 May 2022

The uncertainties are for a confidence probability of approximately 95%.

The calibration result is applied only to the above calibrated item and was found accurate as shown on date and place of calibration only.

This Certificate is issued in accordance with the conditions of accreditation granted by the Thai Laboratory Accreditation scheme which has assessed the measurement capability of the laboratory and its traceability to recognized national standards and to the unit of measurement realized at the corresponding national standards laboratory. This certificate may not be reproduced other than in full except with the prior written approval of the Asia Medical and Agricultural Laboratory and Research Center Public Company Limited (AMARC).



Certificate No. : 22-052237  
Sample Code : 22-19150-002

## REPORT OF CALIBRATION

Equipment : Standard Weight 50 g  
Manufacturer : N/A  
Class : N/A  
Serial No. : N/A  
ID No. : LABE 10/4

### Result of Calibration :

☒ Without adjustment

☐ Adjustment

Conventional value of the result of weighing in air. For a weight taken at a reference temperature ( $t_{ref}$ ) of 20°C, the conventional mass is the mass of a reference weight of a density ( $\rho_{ref}$ ) of 8000 kg.m<sup>-3</sup> which it balances in air of a reference density ( $\rho_a$ ) of 1.2 kg.m<sup>-3</sup>

Description	Deviation	Conventional Mass	Expanded Uncertainty	Maximum Permissible Error	ID No.
	(mg)		(mg)	± (mg)	
50 g	-0.111	49.999889 g	0.10	0.30	LABE 10/4

The result expanded uncertainty of measurement  $U$  is stated as the standard uncertainty of measurement multiplied by the coverage factor  $k = 2.0$ , which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with UKAS M3003

COPY

Certificate No. : 22-052237  
Sample Code : 22-19150-002

Page 3 of 3

## REPORT OF CALIBRATION

### Condition of Calibration

1. Ambient Conditions : Temperature  $20^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$ , Relative humidity  $50\% \pm 10\%$  and air density  $1.18 \text{ kg/m}^3$
2. Calibration Method : WI-CL-007 base on OIML R 111-1 : 2004(E)
3. Reference standard instrument

Instrument	Class	ID No.	Certificate No.	Due Date
1) Standard Weight 1 mg to 1 kg	E2	LB-WE-79	21-079366	22 September 2022

4. This certification is traceable to the International System of Unit maintained at : -

Asia Medical and Agricultural Laboratory and Research Center Public Company Limited

( Instrument number 1).

5. Condition of Calibration item: Normal

### 6. Description of Calibrated Item :

Type and Nominal Value :	Standard Weight 50 g
Shape :	Cylindrical weight with knob
Material :	Stainless steel
Case :	Wooden Box
Comments :	Recalibration

- End of Report -



COPY